

Application News

SSI-LCMS-128

Liquid Chromatography Mass Spectrometry

Separation of 4 Tetrahydrocannabinol Isomers: $\Delta 6a/10a$, $\Delta 8$, $\Delta 9$, and $\Delta 10$ on a Single Quad LCMS

■ Introduction

The cannabis/hemp market continues to grow each year with more states legalizing recreational marijuana, as well as the 2018 Farm Bill removing hemp from the controlled substance list. This bill also defines that any cannabis sativa L. strain with a total tetrahydrocannabinol (THC) concentration of 0.3% or less can be considered hemp and not cannabis. Due to this definition, there is an even higher demand to differentiate hemp from cannabis by determining the correct concentration of THC.

■ LCMS Instrumentation

A Shimadzu LCMS-2020 single quadrupole mass spectrometer coupled with a modular Nexera 40 series UHPLC system was employed for this evaluation. The Nexera 40 series was equipped with a PDA detector for simultaneous analysis of analyte absorbance and compound ionization.

■ LCMS Method Development

Utilizing a Restek Raptor C18 column (150 mm x 2.1 mm, 2.7 μ m) and gradient elution with an overall run time of 8 minutes, four tetrahydrocannabinol isomers, tetrahydrocannabinolic acid, and cannabidiol were chromatographically separated. Mobile phase A consisted of water with 5mM ammonium formate and 0.1% formic acid and Mobile phase B was 50:50 methanol and acetonitrile with 0.1% formic acid. A final flow rate of 0.7mL/min and sample injection volume of 5 μ L was used. Exact gradient conditions are shown in Figure 1.

The MS interface and temperature settings used can be found in Table 1. The final MS method included multiple selected ion monitoring (SIM) events to monitor for each ion of interest as well as scan events for each polarity with a Qarray voltage of 55V. The scan events monitor for in source collision induced disassociation (CID).

Neat calibration curves were run from 0.01ng/ μ L to 10ng/ μ L. Two hemp samples were extracted and analyzed against the calibration curve. Sample extraction was completed using 100mg dry flower geno-grinded for 5 minutes at 1000rpm. Ten mL methanol was added and the sample was vortexed for 1 minute before centrifuging. Then, 1 mL of the supernatant was aliquoted into a sample vial for injection. All calibrators and hemp samples were spiked with delta-9 D3-THC for internal standard quantitation.




Figure 1: LC gradient parameters

Table 1: LCMS-2020 method parameters

Drying Gas	15.0 L/min
Interface Temperature	350 °C
DL Temperature	250 °C
Heat Block Temperature	400 °C
Flow rate	0.7 mL/min
Column Oven Temperature	40 °C
Sample Tray Temperature	15 °C
SIM Channels (m/z)	315.15 (+), 357.15 (-)

■ Results and Discussion

Separation of $\Delta 6a/10a$ THC, $\Delta 8$ THC, $\Delta 9$ THC, and $\Delta 10$ THC was achieved in a single chromatographic run. A comparison of the wavelength 220nm PDA chromatogram and the MS TIC can be seen in Figure 2. Representative SIM chromatograms for all four isomers plus CBD and THCA can be seen in Figure 3.

In addition to a SIM method demonstrating chromatographic separation, a scan event was used to monitor in-source fragmentation. A single Qarray voltage of 55V was applied to fragment all analytes of interest. All four THC isomers plus CBD have a nominal mass of 314.5 g/mol and were differentiated using retention time.

Initial retention time references were confirmed using an individual neat standard. The in-source CID allowed for secondary confirmation of different THC isomers by comparing the fragmentation patterns produced. All four THC isomers produced the m/z fragments of 193 and 259, but m/z 247 is only present with $\Delta 8$ and $\Delta 9$ THC. Further, the fragment m/z 299 is only present for $\Delta 6a/10a$ THC, and fragment m/z 217 is only present for $\Delta 10$ THC (Figure 4). Individual confirmation between $\Delta 8$ THC and $\Delta 9$ THC was completed by the presence of a m/z fragment of 239 (Figure 5).




Figure 2: Representative Chromatograms. Top: PDA wavelength 220nm at 10µg/mL Bottom: SIM 315.15 TIC and SIM 357.15 TIC at 0.5µg/mL




Figure 3: Representative SIM Chromatogram for four THC isomers, CBD, and THCA at 0.5µg/mL




Figure 4: Scan and in-source CID spectra for $\Delta 10$ THC and $\Delta 6a/10a$ THC at $10\mu\text{g/mL}$. Top spectrum is a positive scan event with 0V, bottom spectrum is a positive scan event with 55V in-source CID




Figure 5: 55V In-source CID spectra for $\Delta 8$ THC and $\Delta 9$ THC at $10\mu\text{g/mL}$

All calibration curves were run in triplicate and had at least 5 points with a correlation coefficient of $R^2=0.996$ or better. Internal standard quantitation was used with a $1/C$ weighting for all analytes. $\Delta 9$ -D3-THC was used as the internal standard and was spiked using $25\mu\text{L}$ at a $500\text{ ng}/\mu\text{L}$ concentration for all samples.

The accuracy for each calibrator was between 83.7% and 117.3% with a %RSD of 7.85 or better. The average accuracy and %RSD for each calibrator can be found in Table 2.




Figure 5: Calibration curves for each compound

Table 2: Average accuracy and %RSD for all calibration levels.

(Level 1 = $0.01\text{ ng}/\mu\text{L}$, Level 2 = $0.05\text{ ng}/\mu\text{L}$, Level 3 = $0.1\text{ ng}/\mu\text{L}$, Level 4 = $0.25\text{ ng}/\mu\text{L}$, Level 5 = $0.5\text{ ng}/\mu\text{L}$, Level 6 = $1.0\text{ ng}/\mu\text{L}$, Level 7 = $5.0\text{ ng}/\mu\text{L}$)

Analyte	Level 1 Average Accuracy (%RSD)	Level 2 Average Accuracy (%RSD)	Level 3 Average Accuracy (%RSD)	Level 4 Average Accuracy (%RSD)	Level 5 Average Accuracy (%RSD)	Level 6 Average Accuracy (%RSD)	Level 7 Average Accuracy (%RSD)	R^2
CBD	NA	91 (4.67)	101 (4.28)	98 (2.71)	110 (2.30)	103 (2.55)	99 (1.51)	0.999
$\Delta 9$ -THC	NA	90 (5.3)	101 (1.85)	98 (2.05)	110 (2.79)	102 (2.33)	99 (0.44)	0.999
$\Delta 8$ -THC	NA	90 (6.95)	108 (5.49)	97 (2.32)	107 (2.68)	100 (2.87)	99 (0.93)	0.999
$\Delta 6a/10a$ -THC	98 (7.85)	90 (2.00)	108 (1.40)	100 (2.76)	108 (2.41)	96 (3.48)	NA	0.996
$\Delta 10$ -THC	100 (2.11)	95 (1.87)	108 (2.83)	97 (0.81)	100 (2.80)	NA	NA	0.998
THCA	88 (3.52)	105 (2.58)	108 (0.65)	97 (1.43)	106 (2.15)	97 (0.04)	NA	0.997

Two hemp samples were analyzed against the internal standard calibration curves. Four injections were completed for each hemp sample. One was spiked with $\Delta 8$ -THC, $\Delta 6a/10a$ -THC, and $\Delta 10$ -THC to confirm detection and baseline separation of each isomer in hemp matrix. Representative chromatograms for both hemp samples are shown in Figure 6. The final THC content was calculated two ways. The first is the current protocol using equation 1. The second uses equation 2 and factors in all THC isomers when calculating the final potency of THC. The %THC for each isomer or THCA was calculated using Equation 3.

Final %THC and %RSD for both hemp samples can be found in Table 3 and final THC content can be found in Table 4.

Equation 1: $(\%THCA \times 0.877) + \% \Delta 9\text{-THC} = \text{potency of THC}$
Equation 2: $(\%THCA \times 0.877) + \% \Delta \text{THC sum of isomers} = \text{potency of THC}$

Equation 3: $[\text{THC}] \times (\text{VOL/MG}) \times 100 = \% \text{THC}$




Figure 5: Representative chromatogram for two hemp samples. Hemp sample 2 was spiked with 0.125 ng/ μ L THC isomers $\Delta 8$ THC, $\Delta 6a/10a$ THC, and $\Delta 10$ THC.

Table 3: THC Concentration and Weight % for Hemp Samples (All values based on an average of 4 replicates)

Analyte	Sample 1 Conc. (ng/ μ L)	Sample 1 Wt %	Sample 1 %RSD	Sample 2 Conc. (ng/ μ L)	Sample 2 Wt %	Sample 2 %RSD
CBD	5.243	0.0524	1.5	2.811	0.0281	2.7
Δ 9-THC	0.550	0.0055	3.3	0.440	0.0044	3.6
Δ 8-THC	NA	NA	NA	0.132	0.0013	3.7
Δ 61/10a-THC	NA	NA	NA	0.136	0.0014	3.8
Δ 10-THC	NA	NA	NA	0.142	0.0014	3.7
THCA	4.22	0.0422	2.8	2.37	0.0237	2.3

Table 4: Hemp Samples

	%THC content	%THC content with isomers
Hemp Sample 1	0.0424	0.0424
Hemp Sample 2	0.0252	0.0293

■ Conclusion

A Shimadzu single quadrupole mass spectrometer, LCMS-2020, with in-line PDA demonstrated its capability for simultaneous detection of four THC isomers with baseline separation. Additional confirmation was completed using in-source fragmentation of four THC isomers. Linear calibration curves were acquired for each analyte with a %RSD of 0.996 or better and a % accuracy between 80 and 120 for all calibrators.

Two hemp samples were analyzed to determine the total THC content to be below 0.3% and confirm that baseline separation of all four isomers can be seen in matrix.

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8040



LCMS-8045



LCMS-8050



LCMS-8060NX



LCMS-2020



Q-TOF LCMS-9030

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU SCIENTIFIC INSTRUMENTS, INC.
Applications Laboratory
7102 Riverwood Drive, Columbia, MD 21045
Phone: 800-477-1227 Fax: 410-381-1222
URL <https://www.ssi.shimadzu.com>

For Research Use Only. Not for use in diagnostic procedures. The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

©Shimadzu Scientific Instruments, 2021
First Edition: September 2021