

LC-MS Solutions for

PFAS in Food

Accepted AOAC First Action Method to SMPR 2023.003





Introduction

Per- and polyfluoroalkyl substances (PFAS) are a broad group of thousands of manufactured chemicals characterized by fluorinated organic chains. This includes perfluorooctanoic acid (PFOA) and perfluorooctanesulfonic acid (PFOS), the two most commonly found and historically most widely used PFAS compounds. Widely known as “forever chemicals”, PFAS are persistent pollutants taken up by plants, they can bioaccumulate in animals, impact our soils, waterways and farms, and have records of adverse health effects.

PFAS have been used in manufacturing of consumer goods like carpet, clothing and shoes, cookware, and packaging for decades. Their effectiveness as both an oil and water repellent, a lubricant, as well as a fire retardant, made these substances popular for use in furniture production, linings for take-out containers, firefighting, medicine and countless other applications. Many industries have moved away from using PFAS, but due to that historic use, and since they are extremely difficult to destroy, PFAS are prolific in our environment as residue, leachate, or in legacy products.

With years of data on exposure, scientists have determined adverse health effects stemming from contact with some of the most prolific of these chemicals, and further study is underway. The historical and widespread use has resulted in finding PFAS in the environment and in our food supply. Even after phasing these out from new products, old materials still end up in landfill, and in some industries, they’re still used where we don’t have a suitable replacement. Widespread contamination has led to a necessity of testing, and many groups globally have been working on solving this problem.

Many LC-MS/MS methods have been developed for the targeted analysis of PFAS in environmental samples. Most notably [EPA Methods 537, 537.1, 533, 1633](#). These methods demonstrate testing in the environmental space has been going on for nearly two decades, as EPA 537 V1.1 was finalized in 2009 with years of work prior to putting it together. [ASTM International](#), a group working towards the creation of voluntary consensus standards, also developed methods: D7979, D8421, D7968, and D8535.⁽¹⁻⁷⁾ These methods, similar to those developed by the EPA, are focused on soils, biosolids, but none were designed for testing food products. In the vacuum of food testing methods, [EPA 1633’s](#) inclusion of tissue samples led to the testing of farm animals using this method, but ultimately, it was not designed for testing human food products.

Methods developed by food protection agencies, notably [C-01.02](#) from the FDA and USDA’s [CLG - PFAS 2.03](#)⁽⁸⁻⁹⁾ have been single laboratory validated for PFAS in food. The FDA has recently published an update [C-010.03](#)⁽¹⁰⁾, a single laboratory validated method for 30 PFAS in 9 matrixes, expanding capabilities from the C-01.02.

AOAC International, an industry-led standardization body for food safety had called for methods to test 30 PFAS compounds from 11 food categories to address the lack of standardized food centric methods. The Standard Method Performance Requirements document [SMPR 2023.003](#) detailed the necessary specifications. We present Shimadzu’s response, [an accepted First Action method to AOAC’s OMA: a single LC-MS/MS method to address the needs of routine testing labs.](#)

Table 1: PFAS Target Analytes, Abbreviations, and CAS Numbers

Name	Abbreviation	CAS No.
Perfluorobutanoic acid	PFBA	375-22-4
Perfluoropentanoic acid	PFPeA	2706-90-3
Perfluorohexanoic acid	PFHxA	307-24-4
Perfluoroheptanoic acid	PFHpA	375-85-9
Perfluorooctanoic acid	PFOA	335-67-1
Perfluorononanoic acid	PFNA	375-95-1
Perfluorodecanoic acid	PFDA	335-76-2
Perfluoroundecanoic acid	PFUnA	2058-94-8
Perfluorododecanoic acid	PFDoA	307-55-1
Perfluorotridecanoic acid	PFTTrDA	72629-94-8
Perfluorotetradecanoic acid	PFTeDA	376-06-7
Perfluorobutanesulfonic acid	PFBS	375-73-5
Perfluoropentanesulfonic acid	PFPeS	2706-91-4
Perfluorohexanesulfonic acid	PFHxS	355-46-4
Perfluoroheptanesulfonic acid	PFHpS	375-92-8
Perfluorooctanesulfonic acid	PFOS	1763-23-1
Perfluorononanesulfonic acid	PFNS	68259-12-1
Perfluorodecanesulfonic acid	PFDS	335-77-3
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1
Perfluorododecanesulfonic acid	PFDoS	79780-39-5
Perfluorotridecanesulfonic acid	PFTTrDS	791563-89-8
Perfluorooctanesulfonamide	PFOSA	754-91-6
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9Cl-PF3ONS	756426-58-1
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11Cl-PF3OUdS	763051-92-9
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0

The list of PFAS in AOAC's SMPR is covered largely by the same standards used in EPA method 1633, with the addition of PFUnDS, PFTTrDS, and 10:2 FTS.

Instrument Setup

Two Shimadzu [LCMS-8060NX triple quadrupole LC/MS/MS](#) instruments with heated electrospray ionization (ESI) in negative mode, using [Lab Solutions Software](#), were optimized and coordinated in two different locations. These were used for the collection of all data to rapidly meet the requirements set in [SMPR 2023.003](#) in a collaborative effort. Methods submitted to AOAC undergo a review and approved methods are published in the [Official Methods of Analysis or OMA](#). First they are published as [First Action](#) where they have done a single lab validation. After a First Action acceptance, an inter-laboratory study is needed for submission for approval to a Final Action method. Setting up our initial instrumentation parameters in two locations and duplicating our settings sets a foundation for future validation work.

Mass spectrometry source and ionization conditions were first optimized by introducing the mixed PFAS standard into the [Shimadzu 8060NX MS](#) detector. In these experiments, the following were tested: collision induced dissociation (CID) gas pressure, desolvation line (DL) temperature, drying gas flow rate, heating block temperature, heating gas flow rate, interface temperature, interface voltage, and nebulizer gas flow rate. Small changes in these parameters can affect the sensitivity of the detector and the most optimal conditions can differ compound-to-compound. Determining the best setup for a given method provides the most sensitive results from the detector.

The results were assessed to maximize the response of the various PFAS, with particular focus given to the EU regulated PFAS: PFOA, PFHxS, PFNA, and PFOS.

These optimization experiments were conducted on both instruments at the two selected Shimadzu research centers in [Columbia, Maryland, USA](#) and in [Tokyo, Japan](#). Thirty-eight different variables were tested in each location on all 30 target PFAS analytes for a total of 1984 computations. At the Maryland location, the heating gas flow rate was also tested against the response of PFOA, PFHxS, PFNA, and PFOS as these targets required the most sensitive responses per the SMPR. Confirmation of these optimized conditions in two location allowed for duplicate systems to be run concurrently and share the research workload. [LabSolutions CS software](#) allows for real time sharing of data across the globe to facilitate rapid decision making and collaboration across teams.



LCMS-8060NX Triple Quadrupole Mass Spectrometer



Table 2: MRM Transitions, CE and Quantitation Isotopes

Analyte	Quantitation Ion	Qualifier Ion	Quantitation Ion CE	Qualifier Ion CE	Quantitation labeled Isotope
PFBA	212.9 > 168.8		10		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	9	22	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	9	16	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	11	17	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	10	16	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	10	18	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	13	17	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	19	26	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	13	25	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	12	28	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	33	26	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	39	30	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	44	33	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	49	39	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	54	39	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	54	43	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	45	45	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	52	52	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	54	48	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	54	53	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	31	25	¹³ C ₈ -PFOSA
9CI-PF3ONS	530.8 > 350.8	532.9 > 352.95	25	25	¹³ C ₈ -PFOS
11CI-PF3OUdS	631 > 451	632.9 > 452.95	27	28	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	7	27	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	12	29	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	31	31	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	36	36	¹³ C ₂ -6:2FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	26	46	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	27	45	¹³ C ₇ -PFUnA

Table 3: Isotopically Labeled Internal Standards

Analyte	Quantitation Ion	Qualifier Ion	Quantitation Ion CE	Qualifier Ion CE
$^{13}\text{C}_4$ -PFBA-ES	217 > 172		10	
$^{13}\text{C}_5$ -PFPeA-ES	268 > 222.9		8	
$^{13}\text{C}_5$ -PFHxA-ES	318 > 272.9		9	
$^{13}\text{C}_4$ -PFHpA-ES	367 > 321.9		9	
$^{13}\text{C}_8$ -PFOA-ES	421 > 375.9		9	
$^{13}\text{C}_9$ -PFNA-ES	472 > 426.9		10	
$^{13}\text{C}_6$ -PFDA-ES	519 > 473.8		11	
$^{13}\text{C}_7$ -PFUnA-ES	570 > 524.8		11	
$^{13}\text{C}_2$ -PFDoA-ES	615 > 569.8		11	
$^{13}\text{C}_2$ -PFTeDA-ES ^a	715 > 669.7		13	
$^{13}\text{C}_3$ -PFBS-ES	302 > 79.8		32	
$^{13}\text{C}_3$ -PFHXS-ES	402 > 79.8		42	
$^{13}\text{C}_8$ -PFOS-ES	507 > 79.8		55	
$^{13}\text{C}_8$ -PFOSA-ES	506 > 77.8		33	
$^{13}\text{C}_3$ -HFPO-DA-ES	287 > 168.8		7	
$^{13}\text{C}_2$ -4:2 FTS-ES	329 > 80.9	329 > 308.8	26	18
$^{13}\text{C}_2$ -6:2 FTS-ES ^a	429 > 80.9	429 > 408.8	35	23
$^{13}\text{C}_2$ -8:2 FTS-ES ^a	529 > 79.8	529 > 508.8	44	24

^a not used in SLV, but optional. $^{13}\text{C}_2$ -isotopes tend to suffer more interference particularly for higher mass compounds



PFAS Target Limits and Matrix Categories Set by SMPR 2023.003

Table 4: PFAS Target Limits Set by SMPR 2023.003

SMPR Category	Category 1 (ug/kg)	Category 2 (ug/kg)	Category 3 (ug/kg)
egg	≤0.3	≤3	≤3
seafood	≤0.3	≤3	≤3
meat & fish	≤0.1	≤1	≤1
offal	≤0.4	≤4	≤4
produce	≤0.01	≤1	≤0.1
coffee	≤0.3	≤3	≤3
milk	≤0.01	≤1	≤0.1
protein powder	≤0.08	≤1	≤0.8
fish oil	≤0.5	≤5	≤5
baby food	≤0.01	≤1	≤0.1
pet food	≤0.5	≤5	≤5

With AOAC being an international organization, existing EU regulations were taken into consideration. Four of the target PFAS are already regulated in certain foods in Europe: PFOS, PFOA, PFNA, and PFHxS. These are considered in category 1 of table 4. These are some of the most studied PFAS compounds with respect to exposure and testing, so they require lower limits of quantitation and stricter tolerances. They are regulated in the first 4 matrices of the table: eggs, seafood, meat, fish, and offal. Category 2 includes only PFBA and PFPeA, these are PFAS with only one specific MS/MS transition and have the highest limits of all the targeted compounds. Category 3 contains the remaining 24 PFAS targets which include perfluoroalkyl acids and their precursors. For the study, we used representative examples from each matrix category listed in the SMPR; the products validated for the study per category are shown in table 5.

Table 5: Matrix categories from the SMPR and corresponding products validated in this study

SMPR Matrix Category	Product validated in study
Produce	Carrots
Coffee	Organic brewed coffee
Milk (liquid)	Whole Milk (3.7% fat)
Dairy powders and plant-based protein powder	Organic plant-based protein powder
Eggs	Whole eggs
Seafood	Frozen shrimp
Fish meat and meat of terrestrial mammals	Raw tuna filet
Edible offal of terrestrial animals	Beef kidney
Fish oil	Fish oil capsules
Foods for infants and young children	Strained pumpkin and sweet potato
Pet food and animal feed	Bagged dry dog food

Liquid Chromatography Optimization

In order to ensure the best possible outcomes, the chromatographic conditions, including flow rate, gradient, injection volume were evaluated to find an optimal setup. The TQ 8060NX triple quadrupole LCMS was optimized across over 1900 parameters on two different units duplicated in two locations to ensure replication of the setup. The columns were optimized, as were the gradient conditions and flow rates. The final conditions chosen and used in the study are shown in Table 6. These conditions maximized separation of the PFAS analytes and their sensitivity while minimizing noise and systemic contamination.

Table 6: Chromatographic instrument settings

Instrument Parameter	Value
Mobile Phase A	2mM Ammonium Acetate in 5 / 95 = MeCN / Water
Mobile Phase B	MeCN
Analytical Column	Shimpack Scepter C18-120, 3µm, 2.1 x 100 mm
Delay Column	Shimpack Scepter C18-120, 3µm, 2.1 x 50 mm
Oven Temperature	40 °C

Time, min	Flow rate mL/min	Mobile Phase Composition	
		A %	B %
0	0.3	80	20
10	0.3	0	100
10.01 - 12	0.6	0	100
12.01 - 15	0.3	80	20



Chromatography Optimization and Interferences

An example chromatogram is shown from a sample of carrots which were used to demonstrate the extraction of PFAS from produce. All 30 target PFAS and the respective isotopically labelled internal standards were separated in a nine-minute chromatographic run. A Shimpack Scepter C18-120 delay column was used to isolate possible background contaminant PFAS compounds and avoid coelution with target analytes.




Figure 2: Chromatogram of 0.1 ng/g PFAS in a carrot, representing extraction from produce, with separation of all peaks in nine minutes

Part of optimization for chromatography involved the identification and separation of PFOA from potential cholic acid interferences, per the requirements set forth by AOAC INTERNATIONAL. Three were noted in SMPR 2023.003, Taurodeoxycholic acid (TCDA), Taurochenodeoxycholic acid (TCDCA) and Tauroursodeoxycholic acid (TUDCA).



Nexera HPLC, most used on LCMS-8060NX




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Sample Preparation and Analysis


Consistent sample preparation is key to obtaining the necessary repeatability to meet SMPR 2023.003. A general method was therefore developed for sample preparation involving cleanup with QuEChERS and SPE, then preconcentrating the samples before running the analysis on the Shimadzu LCMS-8060NX. Each matrix type was run using the same instrument method.

Knowing the challenges faced in routine labs, the goal in creating this method was to use the same LCMS-8060NX method, keep the workflow of sample preparation the same but simplify it wherever possible. We chose to omit steps on matrices where we found they were unnecessary, saving time, reagents, and consumable costs. For some matrices some components of the sample preparation did not improve results, so we investigated omitting these steps. We evaluated recovery and repeatability without the addition of the QuEChERS reagents and process, and without concentration after SPE extraction. In the case of milk, dispersive SPE proved more efficient than cartridge SPE.

We found that for offal and fish oil, a QuEChERS treatment was unnecessary, with better or equal recovery and repeatability when compared to the general method versus omitting this step. Similarly, concentration post-SPE cleanup was unnecessary for fish tissue, offal, fish oil, and animal feed or pet food.

The general workflow shown in chart 1 walks through the general procedure of: homogenization, QuEChERS where needed, followed by SPE for all samples, then concentration (again only where needed) before injection to an LC-MS/MS system. Where we could skip a step, we highlighted to do so, so we could streamline the procedure to save testing labs time and money.

Chart 1: General workflow of PFAS sampling



PFAS Confirmation Using QTOF LC/MS

The presence of PFAS compounds that do not have a confirmation ion, such as PFBA and PFPeA, must be confirmed if detected in a sample at or above the LOQ. Confirmation may be by using another column of a different stationary phase, or by using a High Resolution Mass Spectrometer (HRMS), such as the Shimadzu LCMS-9050 Quadrupole Time-of-Flight (QTOF) Liquid Chromatograph Mass Spectrometer. HRMS determines the atomic masses of organic molecules with a high degree of accuracy.

Second confirmation of identity for PFBA by Q-TOF




Figure 4: Confirmation by LC-QTOF for analytes without reference ion, example in egg for PFBA showing both the chromatogram as TIC and extracted ion chromatogram, followed by the mass spectrum compared to the theoretical reference.

Analysis using the LCMS-9050 begins by injecting the extracted sample into the LC operating with the same analytical conditions as used on the LC-MS/MS. Once the analytes elute from the column, they are ionized. The ions travel along the spectrometer's length, separating by their relative charge and mass and a ratio of charge-to-mass is recorded to five decimal places. This accurate mass from the unconfirmed target analyte in the sample can be compared to theoretical values for the mass of the target analyte or to the retention time and accurate mass of a calibration standard injected under the same conditions. Figure 4 shows an extracted ion chromatogram for the exact mass of PFBA (212.9792 m/z) in an egg matrix.



LCMS-9050 Quadrupole Time-of-Flight Liquid Chromatograph Mass Spectrometer

Analysis Results

We used the SMPR as guide for the required LOQ and tested in 3 concentration points: one below the estimated LOQ, and two between 2 – 100 times the estimated LOQ.

Figure 5 shows the recoveries and RSD results for the fish matrix samples, it is used as an example as fish is one of the matrices regulated in the EU. The RSD limit for fish for the 4 EU-regulated compounds was 20%, the scale on the right for RSD does not exceed 20% demonstrating results were well within the required limit. The recovery was required between 80-120%, and on the left scale for recovery rate, this was met for 5, 1 and 0.1 ug/kg level spiked samples.




Figure 5: Recoveries and RSD results for the fish matrix samples

The following page organizes the 30 PFAS included in the AOAC SMPR into groupings based on existing regulations in Europe for category 1, and differing requirements for PFBA and PFPeA than the remaining PFAS targets. EU-regulated matrices for PFAS include eggs, seafood, fish meat and meat of terrestrial animals, and offal. These food groups have stricter limits for PFOS, PFOA, PFNA and PFHxS, and make up the first quadrant of the table. Category 2 includes only PFBA and PFPeA, these are PFAS with only one specific MS/MS transition and have the highest limits of all the targeted compounds. Category 3 contains the remaining 24 PFAS targets which include perfluoroalkyl acids and their precursors.

Limits for each target in each category are followed by the limits of quantitation achieved. Each limit in category 1 was met, in some cases over an order of magnitude lower than the requested value, even as they were some of the lowest concentrations due to existing EU regulations. All limits in category 2 and 3 were met at the SMPR level or we could quantitate lower.

Table 7: Obtained LOQs and SMPR target limits for all matrices Values shown in ug/kg

PFAS	Egg	Seafood	Fish	Offal	Produce	Coffee	Milk	Protein powder	Fish oil	Baby food	Pet food
Limit (ppb) in SMPR per matrix for all Category A	≤0.3	≤0.3	≤0.1	≤0.4	≤0.01	≤0.3	≤0.01	≤0.08	≤0.5	≤0.01	≤0.5
PFOA	0.006	0.006	0.1	0.2	0.006	0.006	0.01	0.055	0.25	0.01	0.5
PFOS	0.055	0.055	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.01	0.5
PFNA	0.006	0.006	0.1	0.2	0.006	0.006	0.01	0.006	0.25	0.01	0.5
PFHxS	0.006	0.006	0.1	0.2	0.006	0.006	0.01	0.006	0.25	0.01	0.5
Limit (ppb) in SMPR per matrix for all Category B	≤3	≤3	≤1	≤4	≤1	≤3	≤1	≤1	≤5	≤1	≤5
PFBA	0.055	0.055	0.1	0.2	0.055	0.055	0.01	0.055	0.25	0.1	0.5
PFPeA	0.055	0.006	0.1	0.2	0.055	0.006	0.01	0.055	0.25	0.1	0.5
Limit (ppb) in SMPR per matrix for all Category C	≤3	≤3	≤1	≤4	≤0.1	≤3	≤0.1	≤0.8	≤5	≤0.1	≤5
PFHxA	0.055	0.006	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.1	0.5
PFHpA	0.006	0.055	0.1	0.2	0.006	0.006	0.01	0.055	0.25	0.01	0.5
PFDA	0.006	0.006	0.1	0.2	0.006	0.055	0.01	0.006	0.25	0.01	0.5
PFUnA	0.006	0.006	0.1	0.2	0.055	0.006	0.01	0.006	0.25	0.01	0.5
PFDoA	0.006	0.055	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.01	0.5
PFTTrDA	0.006	0.055	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.1	0.5
PFTeDA	0.055	0.055	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.1	0.5
PFBS	0.006	0.055	0.1	0.2	0.055	0.006	0.01	0.055	0.25	0.01	0.5
PFPeS	0.006	0.055	0.1	0.2	0.006	0.006	0.01	0.055	0.25	0.01	0.5
PFHpS	0.006	0.006	0.1	0.2	0.006	0.055	0.01	0.006	0.25	0.01	0.5
PFNS	0.006	0.006	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.01	0.5
PFDS	0.006	0.006	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.01	0.5
PFUnDS	0.055	0.006	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.1	0.5
PFDoS	0.055	0.006	0.1	0.2	0.006	0.055	0.01	0.055	0.25	0.01	0.5
PFTTrDS	0.055	0.055	0.1	0.2	0.055	0.55	0.01	0.55	0.25	0.1	0.5
PFOSA	0.55	0.055	0.1	0.2	0.006	0.55	0.01	0.55	0.25	0.1	0.5
9CI-PF3ONS	0.006	0.006	0.1	0.2	0.006	0.006	0.01	0.006	0.25	0.01	0.5
11CI-PF3OUdS	0.006	0.006	0.1	0.2	0.006	0.055	0.01	0.006	0.25	0.01	0.5
HFPO-DA	0.006	0.055	0.1	0.2	0.006	0.055	0.01	0.55	0.25	0.01	0.5
DONA	0.006	0.006	0.1	0.2	0.006	0.006	0.01	0.006	0.25	0.01	0.5
4:2 FTS	0.55	0.055	0.1	0.2	0.055	0.055	0.01	0.055	0.25	0.01	0.5
6:2 FTS	0.006	0.55	0.1	0.2	0.055	0.006	0.01	0.055	0.25	0.01	0.5
8:2 FTS	0.006	0.006	0.1	0.2	0.006	0.55	0.01	0.006	0.25	0.01	0.5
10:2 FTS	0.006	0.006	0.1	0.2	0.006	0.055	0.01	0.055	0.5	0.01	0.5

Results and Discussion

The same procedure from the test samples was used to process the calibration standards. A linear model provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Branched and linear isomers of PFHxS and PFOS were integrated together. Blank matrixes were analyzed in triplicate, depending on the matrix. Four different concentrations from below the LOQ required by the SMPR, and ranging up to 500 times the estimated LOQ, were also analyzed in triplicate to assess recovery and repeatability throughout the calibration and across matrices. Recovery and repeatability at the LOQ in egg matrix for each analyte in category A are given in table 8 to highlight existing EU regulations. Remaining data is appended.




Figure 6: LOQ peaks for PFHxS, and PFOA and their corresponding internal standards

Table 8: Recovery and repeatability for select PFAS at each concentration level in eggs.

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFOA	Blank	0			
	0.006	0.005	5.5	6.19	88.9
	0.055	0.054	2.41	2.46	97.9
	0.55	0.557	1.11	1.1	101.3
	5.5	6.338	2.33	2.02	115.2
PFNA	Blank	0			
	0.006	0.005	1.66	1.67	99.1
	0.055	0.056	2.05	2.01	102
	0.55	0.592	0.59	0.54	107.6
	5.5	6.122	8.8	7.91	111.3
PFOS	Blank	0.002			
	0.006	0.006	5.17	5.28	98.1
	0.055	0.051	2.3	2.49	92.4
	0.55	0.52	0.29	0.31	94.6
	5.5	5.563	1.78	1.76	101.2
PFHxS	Blank	-0.001			
	0.006	0.006	21.39	19.16	111.7
	0.055	0.051	1.16	1.24	93.7
	0.55	0.533	1.27	1.32	96.9
	5.5	5.779	1.15	1.09	105.1

The LOQs for all matrices and compounds were determined according to the SMPR, with retention time, recovery, repeatability, $S/N > 3$ for the qualifier ion and an ion ratio of $\pm 30\%$ all taken into consideration. Figure 6 includes examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards. The LOQ for each analyte was determined based on the specifications listed in the SMPR. The appended application notes detail the process and provide data for each matrix category.

Application News

LCMSTM-8060NX High Performance Liquid Chromatograph Mass Spectrometer
NexeraTM series High Performance Liquid Chromatograph

Determination of 30 PFAS in Foods for Infants and Young Children (Baby Food) by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Toshiya Matsubara¹, Dominika Gruszecka¹, Nozomi Maeshima², Manami Kobayashi²
¹ Shimadzu Scientific Instruments, Inc., ² Shimadzu Corporation

User Benefits

- ◆ Validated method for 30 PFAS in baby food meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in pumpkin and sweet potato puree baby food using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at four concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.1
Perfluoropentanoic acid	PFPeA	2706-90-3	0.1
Perfluorohexanoic acid	PFHxA	307-24-4	0.1
Perfluoroheptanoic acid	PFHpA	375-85-9	0.01
Perfluorooctanoic acid	PFOA	335-67-1	0.01
Perfluorononanoic acid	PFNA	375-95-1	0.01
Perfluorodecanoic acid	PFDA	335-76-2	0.01
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.01
Perfluorododecanoic acid	PFDoA	307-55-1	0.01

Perfluorotridecanoic acid	PFTeDA	72629-94-8	0.1
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.1
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.01
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.01
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.01
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.01
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.01
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.01
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.01
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.1
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.01
Perfluorotridecanesulfonic acid	PFTeDS	791563-89-8	0.1
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.1
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.01
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.01
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.01
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.01
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.01
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.01
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.01
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.01

■ Sample Preparation and Analysis Conditions

Samples were prepared by spiking blended baby food samples in triplicate at four different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using baby food test portions spiked with concentrations of 0.005, 0.01, 0.05, 0.10, 0.50 and 1.0 ng/g. Quantitation was carried out on additional baby food samples spiked in triplicate at 0.01, 0.10, 0.50 and 1.0 ng/g. Since standards were extracted in a baby food matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

Multiple jars of baby food were purchased, combined, and blended. Ten-gram portions were weighed, spiked with target analytes and internal standards, and 10 mL of acetonitrile was added. The samples were shaken for 1 minute and a QuEChERS packet was added. The sample was shaken for 1 minute, allowed to cool to room temperature and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. For greater sensitivity, the extract was concentrated to less than 1 mL, then adjusted to a final 1 mL in a methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.1 ng/g PFAS in a baby food matrix with separation of all peaks in nine minutes




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and at least three different concentrations ranging from the SMPR required LOQ to 100 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery (%)
PFBA	Blank	0.000			
	0.1	0.095	3.90	4.10	95.1
	0.5	0.506	0.38	0.37	101.2
	1	1.006	0.46	0.46	100.6
PFPeA	Blank	0.008			
	0.01	0.011	6.52	5.76	113.3
	0.1	0.101	0.91	0.90	100.9
	0.5	0.481	0.53	0.55	96.1
	1	0.976	0.40	0.41	97.6
PFHxA	Blank	ND			
	0.01	0.010	2.82	2.89	97.5
	0.1	0.099	1.85	1.86	99.4
	0.5	0.475	1.26	1.32	95.0
	1	0.977	1.37	1.41	97.7
PFHpA	Blank	ND			
	0.01	0.010	3.20	3.33	96.2
	0.1	0.104	2.60	2.49	104.1
	0.5	0.488	1.90	1.95	97.6
	1	0.990	0.35	0.35	99.0
PFOA	Blank	ND			
	0.01	0.010	3.04	3.07	99.3
	0.1	0.102	1.95	1.92	101.7
	0.5	0.495	0.95	0.96	99.0
	1	0.994	0.38	0.38	99.4
PFNA	Blank	0.000			
	0.01	0.011	5.44	5.13	106.2
	0.1	0.100	2.35	2.36	99.6
	0.5	0.480	2.24	2.33	96.1
	1	0.975	1.83	1.88	97.5
PFDA	Blank	ND			
	0.01	0.010	5.19	5.04	103.0
	0.1	0.102	1.17	1.15	102.0
	0.5	0.490	2.71	2.76	98.0
	1	1.000	0.52	0.52	100.0
PFUnA	Blank	ND			
	0.01	0.010	8.96	8.91	100.6
	0.1	0.104	0.97	0.94	103.6
	0.5	0.495	1.99	2.01	98.9
	1	0.995	3.20	3.22	99.5
PFDoA	Blank	ND			
	0.01	0.012	10.48	9.10	115.2
	0.1	0.103	6.81	6.63	102.8
	0.5	0.467	1.61	1.72	93.4
	1	0.957	3.90	4.08	95.7
PFTrDA	Blank	ND			
	0.1	0.103	1.04	1.01	103.0
	0.5	0.504	1.34	1.33	100.8
	1	1.024	3.89	3.80	102.4
PFTeDA	Blank	ND			
	0.1	0.117	5.89	5.03	117.3
	0.5	0.557	3.35	3.01	111.3
	1	1.139	3.42	3.00	113.9
PFBS	Blank	ND			
	0.01	0.010	15.82	15.43	102.5
	0.1	0.107	3.90	3.64	107.2
	0.5	0.503	2.70	2.68	100.7
	1	1.023	3.93	3.84	102.3
PFPeS	Blank	0.000			
	0.01	0.011	5.20	4.62	112.6
	0.1	0.104	2.75	2.63	104.4
	0.5	0.483	3.68	3.80	96.7
	1	0.969	1.64	1.70	96.9

PFHxS	Blank	0.001			
	0.01	0.011	7.99	7.45	107.4
	0.1	0.101	2.16	2.15	100.5
	0.5	0.457	2.40	2.62	91.5
	1	0.941	2.06	2.19	94.1
PFHpS	Blank	ND			
	0.01	0.009	2.86	3.10	92.2
	0.1	0.098	2.71	2.75	98.6
	0.5	0.461	4.58	4.97	92.1
	1	0.949	3.25	3.42	94.9
PFOS	Blank	ND			
	0.01	0.009	5.72	6.56	87.2
	0.1	0.099	3.74	3.78	99.1
	0.5	0.462	0.06	0.06	92.3
	1	0.943	1.66	1.76	94.3
PFNS	Blank	ND			
	0.01	0.010	4.62	4.59	100.8
	0.1	0.103	2.44	2.37	103.0
	0.5	0.495	1.33	1.34	99.1
	1	0.991	2.71	2.74	99.1
PFDS	Blank	ND			
	0.01	0.008	1.62	2.04	79.3
	0.1	0.106	6.18	5.83	105.9
	0.5	0.508	1.40	1.38	101.6
	1	1.020	3.31	3.24	102.0
PFUnDS	Blank	ND			
	0.1	0.098	3.71	3.76	98.5
	0.5	0.498	0.55	0.55	99.6
	1	1.005	1.80	1.79	100.5
	PFDoS	Blank	ND		
0.01		0.008	10.92	13.51	80.8
0.1		0.097	2.81	2.91	96.7
0.5		0.407	0.85	1.05	81.4
1		0.808	1.21	1.50	80.8
PFTrDS	Blank	ND			
	0.1	0.100	2.46	2.45	100.4
	0.5	0.504	5.22	5.17	100.8
	1	1.028	4.62	4.50	102.8
	PFOSA	Blank	ND		
0.01		0.010	10.68	10.71	99.7
0.1		0.101	4.27	4.24	100.9
0.5		0.478	1.86	1.94	95.6
1		0.961	0.25	0.26	96.1
9CI-PF3ONS	Blank	0.000			
	0.01	0.011	7.75	7.30	106.2
	0.1	0.104	3.58	3.44	104.0
	0.5	0.485	2.06	2.13	96.9
	1	0.987	2.03	2.06	98.7
11CI-PF3OUdS	Blank	ND			
	0.01	0.010	0.69	0.66	104.4
	0.1	0.099	3.76	3.79	99.2
	0.5	0.478	2.25	2.35	95.7
	1	0.941	4.44	4.72	94.1
HFPO-DA	Blank	ND			
	0.01	0.011	18.97	17.88	106.1
	0.1	0.108	1.46	1.35	107.7
	0.5	0.503	0.93	0.92	100.6
	1	1.009	1.21	1.19	100.9
DONA	Blank	ND			
	0.01	0.011	3.65	3.49	104.5
	0.1	0.103	2.15	2.09	103.0
	0.5	0.486	1.91	1.96	97.3
	1	0.995	0.40	0.41	99.5

4:2 FTS	Blank	0.000			
	0.01	0.012	11.06	9.41	117.5
	0.1	0.105	3.82	3.64	104.9
	0.5	0.492	2.05	2.09	98.2
	1	0.999	0.47	0.47	99.9
6:2 FTS	Blank	ND			
	0.01	0.011	7.86	7.17	109.6
	0.1	0.104	4.90	4.72	103.9
	0.5	0.492	1.23	1.25	98.5
	1	0.981	2.00	2.04	98.1
8:2 FTS	Blank	0.000			
	0.01	0.011	1.06	0.94	112.5
	0.1	0.101	4.33	4.29	100.9
	0.5	0.451	1.90	2.11	90.2
	1	0.895	3.04	3.40	89.5
10:2 FTS	Blank	ND			
	0.01	0.012	10.15	8.88	114.3
	0.1	0.104	3.84	3.69	104.0
	0.5	0.470	2.59	2.76	94.0
	1	0.911	2.92	3.20	91.1

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a baby food matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX

LCMS-8050RX

LCMS-8060RX

LCMS-2020

LCMS-2050

Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-163 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMS™-8060NX High Performance Liquid Chromatograph Mass Spectrometer
Nexera™ series High Performance Liquid Chromatograph

Determination of 30 PFAS in Egg by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Dominika Gruszecka¹, Toshiya Matsubara¹

¹ Shimadzu Scientific Instruments, Inc.

User Benefits

- ◆ Validated method for 30 PFAS in Egg meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in eggs using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at four concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.055
Perfluoropentanoic acid	PFPeA	2706-90-3	0.055
Perfluorohexanoic acid	PFHxA	307-24-4	0.055
Perfluoroheptanoic acid	PFHpA	375-85-9	0.0055
Perfluorooctanoic acid	PFOA	335-67-1	0.0055
Perfluorononanoic acid	PFNA	375-95-1	0.0055
Perfluorodecanoic acid	PFDA	335-76-2	0.0055
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.0055
Perfluorododecanoic acid	PFDoA	307-55-1	0.0055

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.0055
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.055
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.0055
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.0055
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.0055
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.0055
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.055
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.0055
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.0055
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.055
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.055
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.055
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.55
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.0055
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.0055
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.0055
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.0055
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.55
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.0055
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.0055
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.0055

■ Sample Preparation and Analysis Conditions

Large brown eggs were purchased locally. Samples were prepared cracking and removing the whole egg from the shell and grinding with dry ice for 30 seconds at 4000 rpm. The ground material was placed in a beaker. Test portions were spiked in triplicate at four different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram test portions spiked with concentrations of 0.001, 0.01, 0.1, 1.0 and 10.0 ng/g. Quantitation was carried out on additional whole egg samples spiked in triplicate at 0.0055, 0.055, 0.55 and 5.5 ng/g. Since standards were extracted in an egg matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram portions were weighed, spiked with target analytes and internal standards, and 10 mL of acetonitrile was added. The samples were vortexed for 1 minute and a QuEChERS packet was added. The sample was shaken for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. For greater sensitivity, the extract was concentrated to dryness and dissolved in 0.4 mL of a methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.55 ng/g PFAS in an egg matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and four different concentrations ranging from below the SMPR required LOQ to approximately 20 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that ¹³C₂ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3.

The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	0.002			
	0.055	0.051	1.38	1.48	93.3
	0.55	0.528	1.30	1.36	95.9
	5.5	5.664	0.26	0.26	103.0
PFPeA	Blank	0.004			
	0.0055	0.006	10.48	9.67	108.4
	0.055	0.050	4.96	5.42	91.5
	0.55	0.486	0.50	0.57	88.3
PFHxA	5.5	5.364	0.87	0.90	97.5
	Blank	0.004			
	0.0055	0.005	4.52	4.82	93.8
	0.055	0.055	0.62	0.62	100.3
PFHpA	0.55	0.550	1.37	1.36	100.0
	5.5	6.107	1.10	0.99	111.1
	Blank	0.000			
	0.0055	0.005	4.39	4.62	94.9
PFNA	0.055	0.053	1.04	1.09	95.6
	0.55	0.542	0.49	0.50	98.5
	5.5	6.271	0.92	0.80	114.0
	Blank	0.000			
PFOA	0.0055	0.005	5.50	6.19	88.9
	0.055	0.054	2.41	2.46	97.9
	0.55	0.557	1.11	1.10	101.3
	5.5	6.338	2.33	2.02	115.2
PFDA	Blank	0.000			
	0.0055	0.005	1.66	1.67	99.1
	0.055	0.056	2.05	2.01	102.0
	0.55	0.592	0.59	0.54	107.6
PFUnA	5.5	6.122	8.80	7.91	111.3
	Blank	ND			
	0.0055	0.005	3.66	3.91	93.7
	0.055	0.055	0.98	0.99	99.7
PFDoA	0.55	0.569	2.10	2.03	103.5
	5.5	6.085	7.05	6.37	110.7
	Blank	0.000			
	0.0055	0.005	3.20	3.35	95.5
PFTrDA	0.055	0.057	2.15	2.10	102.8
	0.55	0.587	0.85	0.80	106.8
	5.5	6.281	8.37	7.33	114.2
	Blank	ND			
PFTeDA	0.0055	0.004	12.06	16.81	71.7
	0.055	0.047	6.41	7.53	85.1
	0.55	0.504	1.88	2.05	91.6
	5.5	5.892	1.60	1.49	107.1
PFBS	Blank	0.000			
	0.0055	0.005	9.21	10.21	90.2
	0.055	0.047	3.12	3.63	85.8
	0.55	0.515	0.81	0.86	93.6
PFPeS	5.5	6.927	3.26	2.59	125.9
	Blank	0.002			
	0.0055	0.006	24.12	21.48	112.3
	0.055	0.040	0.61	0.83	73.5
PFBS	0.55	0.508	0.91	0.98	92.4
	5.5	7.220	3.10	2.36	131.3
	Blank	ND			
	0.0055	0.004	7.72	10.25	75.4
PFPeS	0.055	0.051	0.61	0.66	93.1
	0.55	0.534	2.17	2.23	97.1
	5.5	5.763	1.50	1.43	104.8
	Blank	0.001			
PFPeS	0.0055	0.005	17.67	18.11	97.6
	0.055	0.051	4.60	4.99	92.2
	0.55	0.516	0.67	0.71	93.9
	5.5	5.494	1.49	1.49	99.9

PFHxS	Blank	ND			
	0.0055	0.006	21.39	19.16	111.7
	0.055	0.051	1.16	1.24	93.7
	0.55	0.533	1.27	1.32	96.9
	5.5	5.779	1.15	1.09	105.1
PFHpS	Blank	0.001			
	0.0055	0.005	2.93	2.98	98.3
	0.055	0.057	3.84	3.73	103.0
	0.55	0.515	1.29	1.38	93.6
	5.5	5.058	3.91	4.25	92.0
PFOS	Blank	0.002			
	0.0055	0.006	5.17	5.28	98.1
	0.055	0.051	2.30	2.49	92.4
	0.55	0.520	0.29	0.31	94.6
	5.5	5.563	1.78	1.76	101.2
PFNS	Blank	0.000			
	0.0055	0.005	10.26	11.84	86.7
	0.055	0.048	1.71	1.94	88.1
	0.55	0.523	2.52	2.65	95.1
	5.5	6.736	3.30	2.69	122.5
PFDS	Blank	ND			
	0.0055	0.005	2.25	2.47	91.3
	0.055	0.050	2.42	2.67	90.5
	0.55	0.503	1.82	1.99	91.5
	5.5	5.541	3.26	3.23	100.8
PFUnDS	Blank	0.001			
	0.0055	0.005	4.47	5.03	88.9
	0.055	0.049	3.86	4.36	88.7
	0.55	0.536	2.55	2.62	97.4
	5.5	6.074	2.62	2.37	110.5
PFDoS	Blank	ND			
	0.0055	0.004	7.73	11.46	67.5
	0.055	0.048	1.70	1.94	87.7
	0.55	0.548	1.64	1.65	99.7
	5.5	5.444	3.69	3.72	99.0
PFTrDS	Blank	ND			
	0.055	0.045	1.86	2.25	82.5
	0.55	0.514	1.37	1.46	93.4
	5.5	6.457	3.92	3.34	117.4
	PFOSA	Blank	ND		
0.0055		0.005	10.12	12.42	81.5
0.055		0.054	2.96	3.00	98.7
0.55		0.539	1.73	1.77	98.1
5.5		5.676	2.21	2.14	103.2
9CI-PF3ONS	Blank	0.000			
	0.0055	0.006	6.82	6.82	100.0
	0.055	0.048	1.26	1.45	86.6
	0.55	0.507	1.59	1.73	92.2
	5.5	7.191	2.75	2.10	130.7
11CI-PF3OUdS	Blank	0.002			
	0.0055	0.006	4.92	4.23	116.2
	0.055	0.045	0.96	1.18	81.4
	0.55	0.494	2.01	2.23	89.9
	5.5	7.332	3.60	2.70	133.3
HFPO-DA	Blank	0.000			
	0.0055	0.005	10.01	10.88	92.0
	0.055	0.052	0.50	0.53	94.7
	0.55	0.533	1.25	1.29	96.9
	5.5	5.695	1.62	1.56	103.5
DONA	Blank	0.000			
	0.0055	0.005	3.26	3.24	100.6
	0.055	0.057	0.53	0.51	104.2
	0.55	0.583	0.31	0.29	106.1
	5.5	6.148	9.36	8.38	111.8
4:2 FTS	Blank	0.017			
	0.0055	0.005	3.67	4.43	82.7
	0.055	0.052	3.28	3.44	95.5
	0.55	0.548	2.03	2.04	99.6
	5.5	5.733	2.18	2.09	104.2

6:2 FTS	Blank	0.000			
	0.0055	0.005	14.57	15.09	96.6
	0.055	0.053	2.80	2.88	97.2
	0.55	0.538	1.99	2.03	97.9
	5.5	5.480	1.60	1.61	99.7
8:2 FTS	Blank	0.000			
	0.0055	0.004	2.76	3.47	79.5
	0.055	0.056	3.32	3.29	101.2
	0.55	0.613	2.81	2.52	111.6
	5.5	4.938	3.12	3.47	89.8
10:2 FTS	Blank	0.000			
	0.0055	0.005	10.90	13.79	79.1
	0.055	0.052	1.12	1.17	95.5
	0.55	0.596	2.71	2.50	108.3
	5.5	4.553	4.08	4.93	82.8

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in an egg matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX



LCMS-8050RX



LCMS-8060RX



LCMS-2020



LCMS-2050



Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-165 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMS™-8060NX High Performance Liquid Chromatograph Mass Spectrometer
Nexera™ series High Performance Liquid Chromatograph

Determination of 30 PFAS in Fish by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Toshiya Matsubara¹, Dominika Gruszecka¹, Nozomi Maeshima², Kota Ishioka², Manami Kobayashi²

¹ Shimadzu Scientific Instruments, Inc., ² Shimadzu Corporation

User Benefits

- ◆ Validated method for 30 PFAS in Fish meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in fish using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at three concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.1
Perfluoropentanoic acid	PFPeA	2706-90-3	0.1
Perfluorohexanoic acid	PFHxA	307-24-4	0.1
Perfluoroheptanoic acid	PFHpA	375-85-9	0.1
Perfluorooctanoic acid	PFOA	335-67-1	0.1
Perfluorononanoic acid	PFNA	375-95-1	0.1
Perfluorodecanoic acid	PFDA	335-76-2	0.1
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.1
Perfluorododecanoic acid	PFDoA	307-55-1	0.1

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.1
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.1
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.1
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.1
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.1
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.1
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.1
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.1
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.1
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.1
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.1
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.1
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.1
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.1
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.1
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.1
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.1
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.1
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.1
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.1
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.1

■ Sample Preparation and Analysis Conditions

Samples were prepared by dicing the edible portion of tuna filets and grinding with dry ice. The ground material was placed in a freezer overnight. Test portions were spiked in triplicate at three different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram test portions spiked with concentrations of 0.05, 0.1, 0.5, 1.0 and 5.0 ng/g. Quantitation was carried out on additional tuna samples spiked in triplicate at 0.1, 1.0, and 5.0 ng/g. Since standards were extracted in a tuna matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram portions were weighed, spiked with target analytes and internal standards, and 10 mL of acetonitrile was added. The samples were shaken by hand for 1 minute and a QuEChERS packet was added. The sample was shaken for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. A 0.5 mL aliquot was transferred to a LC vial and acidified with formic acid.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 1.0 ng/g PFAS in a fish matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and three different concentrations ranging from the SMPR required LOQ to 50 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that $^{13}\text{C}_2$ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3.

The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, $S/N > 3$ for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a $S/N > 10$. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	0.004			
	0.1	0.103	2.51	2.43	103.4
	1	0.970	0.21	0.21	97.0
	5	4.830	0.79	0.82	96.6
PFPeA	Blank	ND			
	0.1	0.096	2.11	2.20	96.2
	1	0.976	0.46	0.47	97.6
	5	4.920	0.61	0.62	98.4
PFHxA	Blank	0.002			
	0.1	0.099	1.64	1.66	98.8
	1	0.957	0.55	0.58	95.7
	5	4.883	1.72	1.77	97.6
PFHpA	Blank	ND			
	0.1	0.100	0.44	0.44	100.1
	1	0.980	1.72	1.76	98.0
	5	4.998	1.42	1.42	100.0
PFOA	Blank	ND			
	0.1	0.100	1.50	1.51	99.6
	1	0.960	1.43	1.49	96.0
	5	4.852	0.81	0.84	97.0
PFNA	Blank	ND			
	0.1	0.101	2.74	2.72	100.9
	1	0.984	1.11	1.12	98.4
	5	5.039	2.25	2.24	100.8
PFDA	Blank	0.000			
	0.1	0.098	3.36	3.42	98.2
	1	0.936	1.39	1.48	93.6
	5	4.800	0.57	0.59	96.0
PFUnA	Blank	0.002			
	0.1	0.100	1.65	1.65	99.8
	1	0.932	1.44	1.55	93.2
	5	4.896	1.63	1.66	97.9
PFDoA	Blank	ND			
	0.1	0.093	10.48	11.29	92.9
	1	0.939	2.17	2.31	93.9
	5	4.833	1.42	1.47	96.6
PFTrDA	Blank	0.004			
	0.1	0.107	8.56	8.05	106.4
	1	0.996	3.93	3.95	99.6
	5	4.889	4.60	4.71	97.8
PFTeDA	Blank	0.009			
	0.1	0.114	8.25	7.23	114.1
	1	1.029	12.47	12.12	102.9
	5	4.860	5.60	5.76	97.2
PFBS	Blank	ND			
	0.1	0.100	2.77	2.78	99.7
	1	0.945	1.76	1.86	94.5
	5	4.836	1.26	1.30	96.7
PFPeS	Blank	ND			
	0.1	0.092	9.10	9.85	92.4
	1	0.990	2.20	2.22	99.0
	5	4.857	1.96	2.01	97.2
PFHxS	Blank	0.002			
	0.1	0.105	6.51	6.21	104.9
	1	0.968	2.02	2.09	96.8
	5	4.838	3.11	3.21	96.7
PFHpS	Blank	0.001			
	0.1	0.099	3.85	3.88	99.1
	1	0.925	2.29	2.47	92.5
	5	4.613	0.35	0.38	92.2

PFOS	Blank	0.003			
	0.1	0.102	1.89	1.85	102.2
	1	0.968	0.60	0.62	96.8
	5	4.834	0.90	0.93	96.7
PFNS	Blank	ND			
	0.1	0.101	6.29	6.26	100.4
	1	0.929	2.48	2.68	92.9
	5	4.840	1.61	1.66	96.8
PFDS	Blank	ND			
	0.1	0.097	4.10	4.24	96.7
	1	0.982	2.63	2.68	98.2
	5	4.877	3.65	3.74	97.6
PFUnDS	Blank	0.000			
	0.1	0.106	13.25	12.58	105.3
	1	1.054	7.89	7.49	105.4
	5	4.759	6.56	6.90	95.2
PFDoS	Blank	0.005			
	0.1	0.114	16.66	14.59	114.1
	1	1.084	20.30	18.72	108.4
	5	4.694	2.70	2.88	93.9
PFTrDS	Blank	0.008			
	0.1	0.119	14.78	12.41	119.1
	1	1.113	24.60	22.11	111.3
	5	4.911	7.02	7.15	98.2
PFOSA	Blank	0.000			
	0.1	0.096	2.16	2.25	96.0
	1	0.945	0.97	1.03	94.5
	5	4.898	0.67	0.68	98.0
9CI-PF3ONS	Blank	ND			
	0.1	0.097	1.14	1.17	97.2
	1	0.965	0.44	0.45	96.5
	5	4.901	1.82	1.86	98.0
11CI-PF3OUdS	Blank	ND			
	0.1	0.097	6.63	6.82	97.3
	1	0.943	3.67	3.89	94.3
	5	4.796	3.06	3.20	95.9
HFPO-DA	Blank	0.004			
	0.1	0.102	6.03	5.91	101.9
	1	0.964	1.36	1.41	96.4
	5	4.778	0.21	0.22	95.5
DONA	Blank	0.001			
	0.1	0.100	1.42	1.43	99.6
	1	0.956	0.90	0.94	95.6
	5	4.974	1.10	1.11	99.5
4:2 FTS	Blank	0.000			
	0.1	0.099	3.73	3.77	98.9
	1	0.955	1.62	1.69	95.5
	5	4.718	1.86	1.97	94.4
6:2 FTS	Blank	0.001			
	0.1	0.098	2.91	2.98	97.7
	1	0.944	1.03	1.09	94.4
	5	4.664	1.00	1.07	93.3
8:2 FTS	Blank	0.010			
	0.1	0.107	4.31	4.02	107.2
	1	1.096	9.63	8.79	109.6
	5	5.041	19.12	18.97	100.8
10:2 FTS	Blank	0.008			
	0.1	0.110	5.70	5.21	109.4
	1	1.029	7.31	7.10	102.9
	5	4.028	8.96	11.12	80.6

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

■ Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a fish matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

■ Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX



LCMS-8050RX



LCMS-8060RX



LCMS-2020



LCMS-2050



Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-167 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMSTM-8060NX High Performance Liquid Chromatograph Mass Spectrometer
NexeraTM series High Performance Liquid Chromatograph

Determination of 30 PFAS in Protein Powder by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Dominika Gruszecka¹, Toshiya Matsubara¹

¹ Shimadzu Scientific Instruments, Inc.

User Benefits

- ◆ Validated method for 30 PFAS in Protein Powder meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in protein powder using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at four concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.055
Perfluoropentanoic acid	PFPeA	2706-90-3	0.055
Perfluorohexanoic acid	PFHxA	307-24-4	0.055
Perfluoroheptanoic acid	PFHpA	375-85-9	0.055
Perfluorooctanoic acid	PFOA	335-67-1	0.055
Perfluorononanoic acid	PFNA	375-95-1	0.0055
Perfluorodecanoic acid	PFDA	335-76-2	0.0055
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.0055
Perfluorododecanoic acid	PFDoA	307-55-1	0.055

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.055
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.055
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.055
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.055
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.0055
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.0055
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.055
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.055
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.055
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.055
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.055
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.55
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.55
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.0055
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.0055
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.55
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.0055
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.055
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.055
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.0055
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.055

■ Sample Preparation and Analysis Conditions

Organic, plant-based protein powder was purchased locally. Test portions were weighed directly from the store packaging and spiked in triplicate at four different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram protein powder test portions spiked with concentrations of 0.001, 0.01, 0.1, 1.0 and 10.0 ng/g. Quantitation was carried out on additional protein powder samples spiked in triplicate at 0.0055, 0.055, 0.55 and 5.5 ng/g. Since standards were extracted in a protein powder matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram portions were weighed, spiked with target analytes and internal standards, 10 mL of water and 10 mL of acetonitrile was added. The samples were shaken for 1 minute and a QuEChERS packet was added. The sample was shaken again for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 13 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. For greater sensitivity, the extract was concentrated to dryness and dissolved in 0.4 mL of a methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.055 ng/g PFAS in a protein powder matrix with separation of all peaks in nine minutes. The chromatogram shows other co-eluting peaks between 2.7 and 3.6 minutes, but they did not interfere.

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUDS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and at least three different concentrations ranging from below the SMPR required LOQ to approximately 70 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that ¹³C₂ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3.

The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	ND			
	0.055	0.071	6.83	5.31	128.7
	0.55	0.599	2.41	2.22	108.9
	5.5	6.077	0.70	0.64	110.5
PFPeA	Blank	ND			
	0.055	0.061	3.70	3.36	110.2
	0.55	0.598	0.29	0.27	108.6
	5.5	6.275	2.88	2.53	114.1
PFHxA	Blank	ND			
	0.055	0.049	2.50	2.81	89.1
	0.55	0.546	2.08	2.10	99.2
	5.5	5.758	1.32	1.26	104.7
PFHpA	Blank	0.005			
	0.055	0.057	1.46	1.41	103.3
	0.55	0.541	4.25	4.32	98.4
	5.5	5.339	2.40	2.47	97.1
PFOA	Blank	0.000			
	0.055	0.060	4.83	4.40	109.7
	0.55	0.614	2.34	2.10	111.6
	5.5	5.962	2.80	2.58	108.4
PFNA	Blank	ND			
	0.0055	0.006	10.70	10.53	101.7
	0.055	0.060	0.67	0.61	109.2
	0.55	0.604	2.45	2.23	109.8
	5.5	6.158	2.84	2.54	111.9
PFDA	Blank	ND			
	0.0055	0.005	4.86	5.09	95.6
	0.055	0.058	1.52	1.45	104.8
	0.55	0.589	3.72	3.48	107.1
	5.5	5.996	3.88	3.55	109.0
PFUnA	Blank	ND			
	0.0055	0.007	21.51	18.99	113.3
	0.055	0.059	6.09	5.73	106.3
	0.55	0.595	2.56	2.36	108.2
	5.5	5.655	4.91	4.78	102.8
PFDoA	Blank	ND			
	0.055	0.060	10.10	9.34	108.2
	0.55	0.597	6.61	6.09	108.6
	5.5	5.748	7.59	7.26	104.5
PFTrDA	Blank	ND			
	0.0055	0.004	12.72	17.77	71.6
	0.055	0.052	2.20	2.30	95.5
	0.55	0.561	2.15	2.10	102.0
	5.5	5.583	3.40	3.35	101.5
PFTeDA	Blank	0.000			
	0.0055	0.005	4.51	4.94	91.2
	0.055	0.049	5.97	6.69	89.3
	0.55	0.568	1.33	1.29	103.3
	5.5	5.432	4.41	4.46	98.8
PFBS	Blank	ND			
	0.055	0.062	11.70	10.33	113.2
	0.55	0.620	7.16	6.35	112.7
	5.5	5.952	6.67	6.16	108.2
PFPeS	Blank	0.001			
	0.055	0.063	3.61	3.14	114.8
	0.55	0.602	6.16	5.62	109.5
	5.5	6.168	0.86	0.77	112.2
PFHxS	Blank	ND			
	0.0055	0.006	6.85	6.02	113.8
	0.055	0.059	3.86	3.60	107.1
	0.55	0.587	1.17	1.10	106.7
	5.5	5.486	2.16	2.17	99.7
PFHpS	Blank	0.001			
	0.0055	0.006	15.90	14.87	106.9
	0.055	0.058	2.05	1.94	105.6
	0.55	0.563	3.63	3.55	102.3
	5.5	6.008	7.58	6.94	109.2

PFOS	Blank	ND			
	0.055	0.057	17.30	16.79	103.0
	0.55	0.603	4.60	4.20	109.6
	5.5	5.325	2.41	2.49	96.8
PFNS	Blank	ND			
	0.055	0.054	4.15	4.24	98.0
	0.55	0.610	8.84	7.97	110.9
	5.5	6.174	5.02	4.48	112.3
PFDS	Blank	0.003			
	0.055	0.056	18.37	18.08	101.6
	0.55	0.532	3.72	3.85	96.6
	5.5	5.808	6.93	6.56	105.6
PFUnDS	Blank	0.001			
	0.0055	0.005	10.13	12.93	78.4
	0.055	0.056	3.07	3.03	101.4
	0.55	0.561	3.32	3.25	101.9
	5.5	4.688	3.45	4.04	85.3
PFDoS	Blank	0.003			
	0.0055	0.006	7.39	6.50	113.7
	0.055	0.058	4.98	4.75	104.8
	0.55	0.523	3.60	3.79	95.0
	5.5	4.703	1.74	2.04	85.5
PFTrDS	Blank	ND			
	0.55	0.695	2.08	1.65	126.4
	5.5	6.531	6.31	5.32	118.8
PFOSA	Blank	ND			
	0.055	0.053	7.62	7.91	96.3
	0.55	0.614	3.36	3.01	111.7
	5.5	5.397	2.50	2.55	98.1
9CI-PF3ONS	Blank	0.000			
	0.0055	0.006	3.21	2.76	116.4
	0.055	0.063	2.40	2.10	114.3
	0.55	0.607	3.84	3.48	110.3
	5.5	5.687	1.52	1.47	103.4
11CI-PF3OUdS	Blank	ND			
	0.0055	0.004	9.37	12.62	74.2
	0.055	0.060	7.15	6.56	109.0
	0.55	0.557	0.90	0.89	101.3
	5.5	5.777	3.99	3.80	105.0
HFPO-DA	Blank	0.030			
	0.55	0.663	8.31	6.89	120.5
	5.5	6.199	8.09	7.18	112.7
DONA	Blank	ND			
	0.0055	0.005	5.20	5.91	88.0
	0.055	0.060	3.83	3.54	108.2
	0.55	0.546	4.27	4.30	99.3
	5.5	5.671	1.33	1.29	103.1
4:2 FTS	Blank	0.012			
	0.055	0.054	5.74	5.85	98.0
	0.55	0.589	2.02	1.89	107.1
	5.5	5.963	5.17	4.77	108.4
6:2 FTS	Blank	0.004			
	0.055	0.071	15.56	12.10	128.6
	0.55	0.582	4.10	3.87	105.9
	5.5	5.729	4.83	4.64	104.1
8:2 FTS	Blank	0.000			
	0.0055	0.005	7.43	7.54	98.6
	0.055	0.055	7.09	7.10	99.8
	0.55	0.553	5.34	5.32	100.5
	5.5	5.004	4.06	4.46	91.0
10:2 FTS	Blank	ND			
	0.055	0.055	10.56	10.61	99.5
	0.55	0.673	8.99	7.35	122.3
	5.5	5.601	11.11	10.91	101.8

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

■ Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a protein powder matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

■ Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX



LCMS-8050RX



LCMS-8060RX



LCMS-2020



LCMS-2050



Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-172 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMS™-8060NX High Performance Liquid Chromatograph Mass Spectrometer
Nexera™ series High Performance Liquid Chromatograph

Determination of 30 PFAS in Pet Food by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Toshiya Matsubara¹, Dominika Gruszecka¹, Nozomi Maeshima², Yui Higashi², Manami Kobayashi²

¹ Shimadzu Scientific Instruments, Inc., ² Shimadzu Corporation

User Benefits

- ◆ Validated method for 30 PFAS in Pet Food meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, Low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in pet food using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at three concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.5
Perfluoropentanoic acid	PFPeA	2706-90-3	0.5
Perfluorohexanoic acid	PFHxA	307-24-4	0.5
Perfluoroheptanoic acid	PFHpA	375-85-9	0.5
Perfluorooctanoic acid	PFOA	335-67-1	0.5
Perfluorononanoic acid	PFNA	375-95-1	0.5
Perfluorodecanoic acid	PFDA	335-76-2	0.5
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.5
Perfluorododecanoic acid	PFDoA	307-55-1	0.5

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.5
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.5
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.5
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.5
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.5
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.5
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.5
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.5
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.5
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.5
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.5
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.5
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.5
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.5
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.5
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.5
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.5
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.5
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.5
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.5
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.5

■ Sample Preparation and Analysis Conditions

Dog food was purchased locally. Sub samples were removed from the bag and crushed with dry ice. The ground material placed into containers and stored in a freezer overnight to allow the dry ice to sublimate away completely. Test portions were spiked in triplicate at three different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 2- gram test portions spiked with concentrations of 0.25, 0.5, 2.5, 5, 12.5, and 25 ng/g. Quantitation was carried out on additional dog food samples spiked in triplicate at 0.5, 5.0, and 25 ng/g. Since standards were extracted in a pet food matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

2-gram portions were weighed, spiked with target analytes and internal standards, 10 mL of water and 10 mL of acetonitrile was added. The samples were shaken by hand for 1 minute and a QuEChERS packet was added. The sample was shaken for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. A 0.5 mL aliquot of the eluted solution was transferred to a vial and acidified with formic acid.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 5 ng/g PFAS in a pet food matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve





Figure 5: PFNA Calibration Curve



Blank matrixes and three different concentrations ranging from the SMPR required LOQ to 50 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Figure 4: PFHxS Calibration Curve





Figure 6: PFOS Calibration Curve



The LOQs for all compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	ND			
	0.5	0.541	13.70	12.66	108.2
	5	4.956	1.69	1.71	99.1
	25	22.753	2.35	2.58	91.0
PFPeA	Blank	ND			
	0.5	0.514	0.65	0.63	102.8
	5	5.071	0.59	0.58	101.4
	25	23.462	0.29	0.31	93.9
PFHxA	Blank	ND			
	0.5	0.502	2.66	2.65	100.3
	5	5.088	0.25	0.25	101.8
	25	24.247	1.07	1.10	97.0
PFHpA	Blank	ND			
	0.5	0.559	3.74	3.35	111.8
	5	5.073	0.85	0.84	101.4
	25	23.973	0.36	0.38	95.9
PFOA	Blank	ND			
	0.5	0.525	2.82	2.68	105.0
	5	5.147	0.81	0.79	102.9
	25	24.564	0.42	0.42	98.3
PFNA	Blank	ND			
	0.5	0.513	0.76	0.74	102.8
	5	5.101	0.95	0.93	102.0
	25	24.739	0.64	0.65	99.0
PFDA	Blank	ND			
	0.5	0.525	1.19	1.14	105.1
	5	5.197	1.80	1.74	103.9
	25	24.522	0.76	0.78	98.1
PFUnA	Blank	ND			
	0.5	0.541	1.25	1.16	108.3
	5	5.141	2.02	1.97	102.8
	25	23.863	0.76	0.79	95.4
PFDoA	Blank	ND			
	0.5	0.482	4.82	5.00	96.4
	5	4.919	3.08	3.13	98.4
	25	23.291	1.75	1.88	93.2
PFTrDA	Blank	ND			
	0.5	0.490	6.38	6.51	97.9
	5	4.921	5.49	5.58	98.4
	25	24.337	0.87	0.90	97.4
PFTeDA	Blank	ND			
	0.5	0.508	1.01	0.99	101.7
	5	4.858	2.87	2.96	97.2
	25	24.158	1.93	2.00	96.6
PFBS	Blank	ND			
	0.5	0.533	11.52	10.80	106.7
	5	5.091	1.32	1.30	101.8
	25	22.671	0.89	0.98	90.7
PFPeS	Blank	ND			
	0.5	0.428	13.61	15.89	85.6
	5	4.429	0.35	0.39	88.6
	25	19.906	5.80	7.28	79.6
PFHxS	Blank	ND			
	0.5	0.539	2.23	2.07	107.7
	5	4.900	2.96	3.02	98.0
	25	22.645	2.36	2.60	90.6
PFHpS	Blank	ND			
	0.5	0.494	2.61	2.64	98.7
	5	5.026	0.82	0.81	100.5
	25	24.918	1.45	1.46	99.7

PFOS	Blank	ND			
	0.5	0.475	6.50	6.84	95.0
	5	4.930	3.29	3.34	98.6
	25	22.823	1.41	1.55	91.3
PFNS	Blank	ND			
	0.5	0.512	7.18	7.01	102.3
	5	4.926	2.71	2.75	98.5
	25	22.611	0.29	0.32	90.4
PFDS	Blank	ND			
	0.5	0.512	9.79	9.57	102.3
	5	4.921	1.68	1.71	98.4
	25	23.794	2.55	2.68	95.2
PFUnDS	Blank	ND			
	0.5	0.531	3.02	2.84	106.1
	5	5.145	1.75	1.70	102.9
	25	24.117	1.78	1.84	96.5
PFDoS	Blank	ND			
	0.5	0.502	5.59	5.57	100.5
	5	5.103	2.54	2.49	102.1
	25	23.501	3.05	3.25	94.0
PFTrDS	Blank	ND			
	0.5	0.511	5.99	5.87	102.1
	5	4.752	5.22	5.49	95.0
	25	23.836	1.86	1.95	95.3
PFOSA	Blank	ND			
	0.5	0.523	0.95	0.91	104.7
	5	5.139	0.62	0.61	102.8
	25	24.134	0.55	0.57	96.5
9CI-PF3ONS	Blank	ND			
	0.5	0.494	2.51	2.54	98.9
	5	5.095	1.51	1.49	101.9
	25	22.288	0.81	0.91	89.1
11CI-PF3OUdS	Blank	ND			
	0.5	0.501	1.86	1.86	100.1
	5	4.919	3.34	3.39	98.4
	25	23.572	1.67	1.78	94.3
HFPO-DA	Blank	ND			
	0.5	0.360	1.91	2.66	72.0
	5	5.044	1.04	1.03	100.9
	25	24.698	2.86	2.89	98.8
DONA	Blank	ND			
	0.5	0.493	6.51	6.59	98.7
	5	4.744	2.14	2.25	94.9
	25	21.780	7.33	8.41	87.1
4:2 FTS	Blank	ND			
	0.5	0.582	4.99	4.29	116.4
	5	5.084	1.85	1.82	101.6
	25	23.506	0.26	0.28	94.0
6:2 FTS	Blank	ND			
	0.5	0.545	3.78	3.47	109.0
	5	5.080	2.71	2.66	101.6
	25	24.619	6.15	6.24	98.5
8:2 FTS	Blank	ND			
	0.5	0.543	0.85	0.79	108.6
	5	5.626	3.15	2.80	112.5
	25	23.754	3.35	3.52	95.0
10:2 FTS	Blank	ND			
	0.5	0.502	10.12	10.08	100.4
	5	5.372	0.85	0.79	107.4
	25	22.176	3.46	3.90	88.7

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

■ Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a pet food matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

■ Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX

LCMS-8050RX

LCMS-8060RX

LCMS-2020

LCMS-2050

Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-170 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMSTM-8060NX High Performance Liquid Chromatograph Mass Spectrometer
NexeraTM series High Performance Liquid Chromatograph

Determination of 30 PFAS in Produce by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Toshiya Matsubara¹, Dominika Gruszecka¹, Nozomi Maeshima², Kota Ishioka², Manami Kobayashi²

¹ Shimadzu Scientific Instruments, Inc., ² Shimadzu Corporation

User Benefits

- ◆ Validated method for 30 PFAS in produce meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, Low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in carrots using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at four concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.055
Perfluoropentanoic acid	PFPeA	2706-90-3	0.055
Perfluorohexanoic acid	PFHxA	307-24-4	0.0055
Perfluoroheptanoic acid	PFHpA	375-85-9	0.0055
Perfluorooctanoic acid	PFOA	335-67-1	0.0055
Perfluorononanoic acid	PFNA	375-95-1	0.0055
Perfluorodecanoic acid	PFDA	335-76-2	0.0055
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.055
Perfluorododecanoic acid	PFDoA	307-55-1	0.0055

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.0055
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.0055
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.055
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.0055
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.0055
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.0055
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.0055
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.0055
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.0055
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.0055
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.0055
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.055
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.0055
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.0055
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.0055
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.0055
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.0055
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.055
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.055
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.0055
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.0055

■ Sample Preparation and Analysis Conditions

Samples were prepared by dicing the tuber portion of carrots and grinding with dry ice. The ground material was placed in a freezer overnight. Test portions were spiked in triplicate at four different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram carrot test portions spiked with concentrations of 0.001, 0.01, 0.1, 1.0 and 10.0 ng/g. Quantitation was carried out on additional carrot samples spiked in triplicate at 0.0055, 0.055, 0.55 and 5.5 ng/g. Since standards were extracted in a carrot matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram portions were weighed, spiked with target analytes and internal standards, and 10 mL of acetonitrile was added. The samples were vortexed for 1 minute and a QuEChERS packet was added. The sample was shaken for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. For greater sensitivity, the extract was concentrated to dryness and dissolved in 0.4 mL of a methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.055 ng/g PFAS in a carrot matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and at least three different concentrations ranging from below the SMPR required LOQ to approximately 500 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that $^{13}\text{C}_2$ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3.

The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, $S/N > 3$ for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a $S/N > 10$. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	0.000			
	0.055	0.047	1.74	2.03	85.5
	0.55	0.504	0.40	0.44	91.7
	5.5	5.485	0.32	0.32	99.7
PFPeA	Blank	ND			
	0.0055	0.005	9.98	11.62	85.9
	0.055	0.051	0.40	0.44	92.4
	0.55	0.519	1.91	2.02	94.3
PFHxA	Blank	0.001			
	0.0055	0.005	9.51	11.41	83.4
	0.055	0.057	3.03	2.94	103.0
	0.55	0.533	1.36	1.40	96.9
PFHpA	Blank	0.000			
	0.0055	0.005	2.80	3.20	87.5
	0.055	0.054	1.46	1.50	97.1
	0.55	0.558	0.87	0.86	101.4
PFOA	Blank	0.001			
	0.0055	0.006	3.55	3.30	107.3
	0.055	0.055	1.16	1.16	99.8
	0.55	0.560	0.65	0.64	101.7
PFNA	Blank	0.001			
	0.0055	0.006	1.45	1.38	104.8
	0.055	0.057	0.67	0.64	103.6
	0.55	0.592	2.73	2.54	107.7
PFDA	Blank	0.000			
	0.0055	0.005	14.07	15.03	93.6
	0.055	0.058	1.85	1.77	104.7
	0.55	0.606	0.85	0.78	110.1
PFUnA	Blank	0.000			
	0.0055	0.004	4.53	5.65	80.3
	0.055	0.052	1.79	1.89	94.8
	0.55	0.548	2.22	2.23	99.6
PFDoA	Blank	0.000			
	0.0055	0.005	6.79	7.44	91.2
	0.055	0.044	4.99	6.22	80.2
	0.55	0.483	16.54	18.85	87.8
PFTrDA	Blank	0.000			
	0.0055	0.005	7.05	8.34	84.6
	0.055	0.041	1.95	2.58	75.4
	0.55	0.477	11.69	13.46	86.9
PFTeDA	Blank	0.001			
	0.0055	0.005	5.88	6.00	97.9
	0.055	0.043	1.99	2.51	79.3
	0.55	0.484	10.52	11.96	88.0
PFBS	Blank	0.003			
	0.055	0.050	2.00	2.21	90.6
	0.55	0.520	0.56	0.59	94.5
	5.5	5.669	0.06	0.06	103.1
PFPeS	Blank	0.000			
	0.0055	0.005	1.13	1.16	97.5
	0.055	0.054	2.88	2.95	97.7
	0.55	0.504	1.66	1.82	91.6
PFHxS	Blank	0.000			
	0.0055	0.005	6.82	7.21	94.5
	0.055	0.050	1.82	1.98	91.8
	0.55	0.514	1.02	1.09	93.5
	5.5	5.182	10.95	11.62	94.2

PFHpS	Blank	0.002			
	0.0055	0.007	10.64	8.67	122.8
	0.055	0.060	3.50	3.22	108.8
	0.55	0.545	6.17	6.22	99.1
	5.5	6.013	0.38	0.35	109.3
PFOS	Blank	ND			
	0.0055	0.006	3.72	3.77	98.9
	0.055	0.050	2.49	2.75	90.4
	0.55	0.511	1.08	1.16	93.0
	5.5	5.627	0.15	0.15	102.3
PFNS	Blank	0.000			
	0.0055	0.005	18.78	20.65	90.9
	0.055	0.048	1.61	1.86	86.5
	0.55	0.488	12.40	13.98	88.7
	5.5	4.823	2.00	2.28	87.7
PFDS	Blank	0.000			
	0.0055	0.004	5.30	6.84	77.5
	0.055	0.046	3.56	4.24	84.1
	0.55	0.501	10.91	11.96	91.2
	5.5	5.495	3.06	3.06	99.9
PFUnDS	Blank	0.001			
	0.0055	0.005	0.86	0.91	95.0
	0.055	0.055	2.30	2.31	99.7
	0.55	0.594	1.35	1.25	107.9
	5.5	6.127	0.92	0.82	111.4
PFDoS	Blank	0.000			
	0.0055	0.006	18.94	16.15	117.3
	0.055	0.064	1.96	1.68	117.0
	0.55	0.638	13.06	11.26	116.0
	5.5	6.836	1.31	1.05	124.3
PFTrDS	Blank	0.000			
	0.0055	0.005	1.42	1.45	97.9
	0.055	0.051	0.89	0.95	93.2
	0.55	0.529	7.67	7.97	96.2
	5.5	5.722	4.76	4.58	104.0
PFOSA	Blank	ND			
	0.0055	0.005	4.88	5.54	88.1
	0.055	0.053	1.66	1.71	96.6
	0.55	0.528	1.40	1.46	96.0
	5.5	5.828	2.20	2.07	106.0
9CI-PF3ONS	Blank	0.000			
	0.0055	0.006	3.20	2.93	109.1
	0.055	0.059	1.42	1.33	106.2
	0.55	0.604	2.48	2.26	109.8
	5.5	6.603	0.15	0.13	120.0
11CI-PF3OUdS	Blank	0.000			
	0.0055	0.005	0.95	0.97	97.4
	0.055	0.050	1.57	1.73	91.0
	0.55	0.530	11.47	11.91	96.4
	5.5	6.030	2.40	2.19	109.6
HFPO-DA	Blank	0.000			
	0.0055	0.005	4.01	4.32	92.9
	0.055	0.052	1.37	1.43	95.4
	0.55	0.538	0.79	0.81	97.8
	5.5	5.755	0.75	0.72	104.6
DONA	Blank	0.000			
	0.0055	0.006	5.41	5.26	102.8
	0.055	0.056	0.30	0.30	101.5
	0.55	0.590	0.87	0.81	107.3
	5.5	5.373	0.51	0.53	97.7
4:2 FTS	Blank	ND			
	0.0055	0.006	11.31	11.28	100.2
	0.055	0.049	0.96	1.08	88.8
	0.55	0.544	0.21	0.21	98.9
	5.5	5.793	0.25	0.24	105.3

6:2 FTS	Blank	0.002			
	0.055	0.051	0.85	0.92	92.4
	0.55	0.493	0.67	0.74	89.7
	5.5	5.107	1.99	2.14	92.8
8:2 FTS	Blank	0.000			
	0.0055	0.006	3.26	2.78	117.2
	0.055	0.057	3.83	3.68	104.0
	0.55	0.649	7.42	6.29	118.1
	5.5	6.487	2.75	2.33	117.9
10:2 FTS	Blank	0.000			
	0.0055	0.005	3.15	3.39	93.0
	0.055	0.050	0.86	0.95	90.4
	0.55	0.512	11.03	11.85	93.1
	5.5	4.647	2.25	2.66	84.5

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a carrot matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX

LCMS-8050RX

LCMS-8060RX

LCMS-2020

LCMS-2050

Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-171 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMS™-8060NX High Performance Liquid Chromatograph Mass Spectrometer
Nexera™ series High Performance Liquid Chromatograph

Determination of 30 PFAS in Offal by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Dominika Gruszecka¹, Eishi Imoto¹, Toshiya Matsubara¹

¹ Shimadzu Scientific Instruments, Inc.

User Benefits

- ◆ Validated method for 30 PFAS in Offal meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in beef kidney using Acetonitrile and Solid Phase Extraction (SPE) followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at five concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.2
Perfluoropentanoic acid	PFPeA	2706-90-3	0.2
Perfluorohexanoic acid	PFHxA	307-24-4	0.2
Perfluoroheptanoic acid	PFHpA	375-85-9	0.2
Perfluorooctanoic acid	PFOA	335-67-1	0.2
Perfluorononanoic acid	PFNA	375-95-1	0.2
Perfluorodecanoic acid	PFDA	335-76-2	0.2
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.2
Perfluorododecanoic acid	PFDoA	307-55-1	0.2

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.2
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.2
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.2
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.2
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.2
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.2
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.2
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.2
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.2
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.2
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.2
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.2
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.2
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.2
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.2
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.2
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.2
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.2
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.2
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.2
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.2

■ Sample Preparation and Analysis Conditions

Samples were prepared by cutting locally purchased beef kidney into slices, freezing and grinding with dry ice. The ground material was placed in a freezer overnight. Test portions were spiked in triplicate at five different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram kidney test portions spiked with concentrations of 0.1, 0.5, 1.5, 5.0 and 15 ng/g. Quantitation was carried out on additional kidney samples spiked in triplicate at 0.2, 0.4, 1, 4, and 10 ng/g. Since standards were extracted in a kidney matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram test portions were weighed, spiked with target analytes and internal standards, and 10 mL of acetonitrile was added. The samples were vortexed for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with a basic methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.4 ng/g PFAS in a beef kidney matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHXS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and five different concentrations ranging from below the SMPR required LOQ to 25 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that ¹³C₂ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3.

The LOQs for all compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	0.002			
	0.2	0.212	3.23	3.04	106.3
	0.4	0.442	1.78	1.61	110.5
	1.0	1.131	5.12	4.52	113.1
	4.0	4.102	2.55	2.49	102.5
	10	10.121	3.84	3.80	101.2
PFPeA	Blank	0.010			
	0.2	0.213	3.58	3.35	106.7
	0.4	0.434	1.53	1.41	108.5
	1.0	1.068	4.36	4.08	106.8
	4.0	3.908	2.66	2.72	97.7
	10	9.583	2.29	2.38	95.9
PFHxA	Blank	0.003			
	0.2	0.216	1.08	1.00	107.9
	0.4	0.454	4.87	4.29	113.5
	1.0	1.195	6.03	5.04	119.5
	4.0	4.441	3.84	3.46	111.0
	10	10.797	4.18	3.87	108.0
PFHpA	Blank	0.007			
	0.2	0.208	4.47	4.30	103.8
	0.4	0.412	0.46	0.45	102.9
	1.0	1.014	5.75	5.67	101.4
	4.0	3.690	1.90	2.06	92.3
	10	9.120	3.12	3.42	91.2
PFOA	Blank	0.012			
	0.2	0.221	3.89	3.52	110.7
	0.4	0.453	2.65	2.34	113.3
	1.0	1.174	6.21	5.29	117.4
	4.0	4.331	2.69	2.48	108.3
	10	10.503	3.84	3.65	105.0
PFNA	Blank	0.007			
	0.2	0.205	3.75	3.67	102.1
	0.4	0.435	2.77	2.55	108.9
	1.0	1.139	7.18	6.30	113.9
	4.0	4.204	4.53	4.31	105.1
	10	10.421	2.91	2.79	104.2
PFDA	Blank	0.020			
	0.2	0.211	3.42	3.25	105.3
	0.4	0.433	0.55	0.51	108.1
	1.0	1.053	8.41	7.98	105.3
	4.0	3.936	2.52	2.56	98.4
	10	9.547	3.61	3.78	95.5
PFUnA	Blank	0.011			
	0.2	0.213	1.89	1.77	106.4
	0.4	0.433	4.47	4.12	108.3
	1.0	1.094	4.58	4.18	109.4
	4.0	4.070	2.89	2.84	101.7
	10	9.848	4.64	4.71	98.5
PFDoA	Blank	0.016			
	0.2	0.226	3.26	2.88	113.2
	0.4	0.462	3.03	2.62	115.5
	1.0	1.124	2.03	1.80	112.4
	4.0	4.365	2.26	2.07	109.1
	10	10.375	3.51	3.38	103.7
PFTrDA	Blank	0.001			
	0.2	0.214	2.86	2.66	107.3
	0.4	0.433	3.99	3.69	108.2
	1.0	0.973	7.74	7.95	97.3
	4.0	3.927	1.39	1.41	98.2
	10	9.297	2.55	2.74	92.9
PFTeDA	Blank	0.008			
	0.2	0.205	2.66	2.60	102.5
	0.4	0.427	0.99	0.92	106.8
	1.0	0.990	4.18	4.22	99.0
	4.0	3.880	0.98	1.02	97.0
	10	9.041	0.56	0.62	90.4
PFBS	Blank	0.011			
	0.2	0.194	11.88	12.27	96.8
	0.4	0.415	9.99	9.64	103.6
	1.0	1.035	4.94	4.78	103.5
	4.0	3.824	0.17	0.18	95.6
	10	9.594	7.39	7.71	95.9

PFPeS	Blank	0.025			
	0.2	0.225	6.70	5.95	112.8
	0.4	0.455	2.64	2.31	113.9
	1.0	1.042	3.84	3.69	104.2
	4.0	4.009	0.78	0.78	100.2
PFHxS	10	9.704	3.97	4.09	97.0
	Blank	0.026			
	0.2	0.203	6.37	6.29	101.3
	0.4	0.383	1.65	1.72	95.7
	1.0	0.962	5.37	5.58	96.2
PFHpS	4.0	3.953	4.07	4.12	98.8
	10	9.357	3.12	3.34	93.6
	Blank	0.025			
	0.2	0.216	6.23	5.77	107.9
	0.4	0.456	5.34	4.68	114.0
PFOS	1.0	1.110	7.64	6.88	111.0
	4.0	4.110	2.21	2.15	102.8
	10	9.886	4.35	4.40	98.9
	Blank	0.027			
	0.2	0.212	1.90	1.80	105.9
PFNS	0.4	0.429	1.97	1.84	107.1
	1.0	1.092	4.66	4.27	109.2
	4.0	4.071	5.24	5.15	101.8
	10	10.156	7.04	6.94	101.6
	Blank	0.015			
PFDS	0.2	0.223	2.50	2.24	111.4
	0.4	0.437	4.84	4.43	109.4
	1.0	1.142	7.23	6.33	114.2
	4.0	4.119	1.61	1.56	103.0
	10	9.991	5.09	5.09	99.9
PFUnDS	Blank	ND			
	0.2	0.214	6.66	6.23	106.8
	0.4	0.444	7.66	6.91	110.9
	1.0	1.122	11.22	10.00	112.2
	4.0	4.189	0.80	0.77	104.7
PFDoS	10	9.983	4.23	4.24	99.8
	Blank	0.030			
	0.2	0.239	9.94	8.31	119.6
	0.4	0.438	3.92	3.59	109.5
	1.0	1.005	5.98	5.95	100.5
PFTrDS	4.0	4.038	3.74	3.71	101.0
	10	9.503	5.35	5.63	95.1
	Blank	0.016			
	0.2	0.232	4.32	3.73	115.9
	0.4	0.459	2.71	2.36	114.7
PFOSA	1.0	1.081	0.83	0.77	108.1
	4.0	4.281	3.52	3.29	107.0
	10	9.839	0.96	0.98	98.4
	Blank	0.009			
	0.2	0.222	5.77	5.20	111.0
9CI-PF3ONS	0.4	0.475	12.42	10.46	118.7
	1.0	1.088	2.72	2.50	108.8
	4.0	4.264	4.76	4.47	106.6
	10	9.680	2.69	2.78	96.8
	Blank	ND			
11CI-PF3OUdS	0.2	0.232	24.78	21.34	116.1
	0.4	0.419	13.60	12.98	104.7
	1.0	1.110	6.80	6.12	111.0
	4.0	4.157	5.28	5.09	103.9
	10	10.422	6.27	6.02	104.2
9CI-PF3ONS	Blank	0.016			
	0.2	0.202	2.37	2.34	101.1
	0.4	0.428	5.16	4.83	106.9
	1.0	1.143	4.80	4.20	114.3
	4.0	4.194	2.87	2.74	104.9
11CI-PF3OUdS	10	10.461	6.84	6.54	104.6
	Blank	0.008			
	0.2	0.210	2.53	2.41	104.9
	0.4	0.458	8.34	7.29	114.5
	1.0	1.123	7.60	6.77	112.3
11CI-PF3OUdS	4.0	4.352	2.84	2.61	108.8
	10	10.038	3.95	3.94	100.4

HFPO-DA	Blank	0.016			
	0.2	0.207	3.86	3.72	103.6
	0.4	0.418	4.55	4.36	104.5
	1.0	1.084	6.85	6.32	108.4
	4.0	3.938	4.03	4.09	98.4
	10	9.601	3.20	3.34	96.0
DONA	Blank	0.011			
	0.2	0.211	1.14	1.08	105.2
	0.4	0.428	0.55	0.51	107.0
	1.0	1.042	6.69	6.43	104.2
	4.0	3.827	1.91	2.00	95.7
	10	9.590	2.75	2.86	95.9
4:2 FTS	Blank	0.016			
	0.2	0.213	5.91	5.54	106.7
	0.4	0.422	6.26	5.92	105.6
	1.0	1.056	6.91	6.54	105.6
	4.0	3.881	1.99	2.05	97.0
	10	9.579	3.63	3.79	95.8
6:2 FTS	Blank	0.028			
	0.2	0.218	15.62	14.33	109.1
	0.4	0.431	2.45	2.27	107.7
	1.0	1.149	5.08	4.42	114.9
	4.0	4.304	3.14	2.92	107.6
	10	10.209	6.68	6.55	102.1
8:2 FTS	Blank	0.022			
	0.2	0.164	2.83	3.46	81.8
	0.4	0.340	4.56	5.36	85.0
	1.0	0.952	1.18	1.24	95.2
	4.0	3.974	9.63	9.70	99.3
	10	9.327	3.07	3.29	93.3
10:2 FTS	Blank	0.016			
	0.2	0.217	9.19	8.47	108.5
	0.4	0.489	14.16	11.58	122.4
	1.0	1.251	8.26	6.61	125.1
	4.0	4.956	4.67	3.77	123.9
	10	11.777	6.18	5.25	117.8

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

■ Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in an offal matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

■ Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX



LCMS-8050RX



LCMS-8060RX



LCMS-2020



LCMS-2050



Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-169 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMSTM-8060NX High Performance Liquid Chromatograph Mass Spectrometer
NexeraTM series High Performance Liquid Chromatograph

Determination of 30 PFAS in Milk by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Toshiya Matsubara¹, Dominika Gruszecka¹, Nozomi Maeshima², Yui Higashi², Manami Kobayashi²

¹ Shimadzu Scientific Instruments, Inc., ² Shimadzu Corporation

User Benefits

- ◆ Validated method for 30 PFAS in Milk meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in milk using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at three concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.01
Perfluoropentanoic acid	PFPeA	2706-90-3	0.01
Perfluorohexanoic acid	PFHxA	307-24-4	0.01
Perfluoroheptanoic acid	PFHpA	375-85-9	0.01
Perfluorooctanoic acid	PFOA	335-67-1	0.01
Perfluorononanoic acid	PFNA	375-95-1	0.01
Perfluorodecanoic acid	PFDA	335-76-2	0.01
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.01
Perfluorododecanoic acid	PFDoA	307-55-1	0.01

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.01
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.01
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.01
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.01
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.01
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.01
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.01
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.01
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.01
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.01
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.01
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.01
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.01
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.01
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.01
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.01
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.01
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.01
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.01
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.01
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.01

■ Sample Preparation and Analysis Conditions

Whole milk (3.7 % fat) was purchased and sampled directly from the carton. Test portions were spiked in triplicate at three different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram test portions spiked with 0.005, 0.01, 0.05, 0.10, 0.50, and 1.0 ng/g of each target analyte. Quantitation was carried out on additional whole milk samples spiked in triplicate at 0.01, 0.10 and 1.0 ng/g. Since standards were extracted in milk matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram portions were weighed, spiked with target analytes and internal standards, and 150 µL of formic acid and 10 mL of acetonitrile was added. The samples were shaken by hand for 10 seconds and a QuEChERS packet was added. The sample was shaken again for 5 minutes and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube to which additional QuEChERS reagent was added. The sample was shaken for 10 seconds by hand, for 5 minutes on a shaker, and then centrifuged for 5 minutes at 4000 rpm. An aliquot was removed to another tube and concentrated to less than 1 ml under nitrogen, reconstituted to 1 ml in a methanol-water mixture, transferred to a 1.5 mL tube and centrifuged for 10 minutes at 15,000 rpm. The supernatant was transferred to an LC vial for analysis.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.1 ng/g PFAS in a milk matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and three different concentrations ranging from the SMPR required LOQ to 100 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that ¹³C₂ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3.

The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	0.000			
	0.01	0.012	8.20	7.02	116.8
	0.1	0.103	2.10	2.02	103.5
	1	0.973	0.57	0.58	97.3
PFPeA	Blank	0.000			
	0.01	0.011	2.63	2.45	107.1
	0.1	0.099	0.66	0.66	99.3
	1	0.947	0.38	0.40	94.7
PFHxA	Blank	ND			
	0.01	0.010	1.74	1.67	104.4
	0.1	0.096	0.93	0.97	96.1
	1	0.947	0.10	0.11	94.7
PFHpA	Blank	0.001			
	0.01	0.011	7.89	6.79	116.2
	0.1	0.099	0.45	0.45	99.4
	1	0.968	1.50	1.55	96.8
PFOA	Blank	ND			
	0.01	0.010	11.84	11.35	104.3
	0.1	0.098	2.89	2.93	98.5
	1	0.977	1.00	1.02	97.7
PFNA	Blank	0.000			
	0.01	0.011	0.96	0.86	112.2
	0.1	0.100	0.92	0.93	99.5
	1	0.976	0.61	0.62	97.6
PFDA	Blank	0.001			
	0.01	0.011	2.91	2.59	112.7
	0.1	0.104	0.42	0.40	104.2
	1	0.994	2.71	2.73	99.4
PFUnA	Blank	0.000			
	0.01	0.011	1.20	1.08	111.5
	0.1	0.101	0.84	0.83	100.9
	1	0.980	0.72	0.74	98.0
PFDoA	Blank	0.000			
	0.01	0.012	13.65	11.65	117.1
	0.1	0.101	1.54	1.52	101.3
	1	0.951	1.47	1.55	95.1
PFTrDA	Blank	0.000			
	0.01	0.011	2.40	2.24	106.9
	0.1	0.101	1.60	1.59	100.9
	1	0.994	0.68	0.68	99.4
PFTeDA	Blank	0.000			
	0.01	0.011	4.14	3.80	108.8
	0.1	0.098	0.67	0.68	98.2
	1	0.974	3.40	3.49	97.4
PFBS	Blank	ND			
	0.01	0.012	8.49	7.19	118.1
	0.1	0.102	1.65	1.62	101.6
	1	0.958	0.25	0.26	95.8
PFPeS	Blank	ND			
	0.01	0.010	11.26	11.41	98.7
	0.1	0.102	1.05	1.03	101.9
	1	0.953	1.62	1.70	95.3
PFHxS	Blank	0.000			
	0.01	0.010	10.50	10.76	97.6
	0.1	0.099	3.88	3.89	99.6
	1	0.956	1.36	1.42	95.6
PFHpS	Blank	ND			
	0.01	0.010	7.26	6.85	106.0
	0.1	0.105	2.27	2.17	104.4
	1	0.983	2.10	2.14	98.3

PFOS	Blank	ND			
	0.01	0.010	6.44	6.54	98.4
	0.1	0.099	2.57	2.59	99.2
	1	0.961	1.32	1.37	96.1
PFNS	Blank	ND			
	0.01	0.010	10.24	10.20	100.4
	0.1	0.102	2.86	2.81	101.8
	1	0.992	2.58	2.60	99.2
PFDS	Blank	ND			
	0.01	0.009	3.74	3.95	94.8
	0.1	0.103	6.13	5.95	103.0
	1	1.015	1.22	1.20	101.5
PFUnDS	Blank	ND			
	0.01	0.011	2.19	2.01	109.3
	0.1	0.105	4.78	4.56	105.0
	1	0.988	3.49	3.54	98.8
PFDoS	Blank	ND			
	0.01	0.010	8.88	9.15	97.0
	0.1	0.101	2.24	2.22	100.6
	1	0.963	2.27	2.36	96.3
PFTrDS	Blank	ND			
	0.01	0.010	5.12	4.83	106.0
	0.1	0.099	3.30	3.32	99.4
	1	0.990	2.08	2.10	99.0
PFOSA	Blank	ND			
	0.01	0.011	3.59	3.36	106.8
	0.1	0.099	1.83	1.85	99.0
	1	0.950	0.55	0.58	95.0
9CI-PF3ONS	Blank	0.000			
	0.01	0.010	4.42	4.31	102.5
	0.1	0.103	0.96	0.94	102.5
	1	0.971	3.15	3.25	97.1
11CI-PF3OUdS	Blank	ND			
	0.01	0.009	5.32	5.61	94.8
	0.1	0.101	1.10	1.09	100.8
	1	0.992	1.22	1.23	99.2
HFPO-DA	Blank	ND			
	0.01	0.010	7.16	7.19	99.6
	0.1	0.099	2.17	2.18	99.3
	1	0.932	0.80	0.86	93.2
DONA	Blank	0.000			
	0.01	0.010	0.80	0.77	103.9
	0.1	0.099	0.42	0.42	98.9
	1	0.961	2.12	2.21	96.1
4:2 FTS	Blank	ND			
	0.01	0.011	4.97	4.63	107.5
	0.1	0.102	2.75	2.69	102.1
	1	0.989	1.28	1.29	98.9
6:2 FTS	Blank	0.001			
	0.01	0.012	11.91	10.39	114.6
	0.1	0.106	5.06	4.80	105.5
	1	0.984	1.40	1.43	98.4
8:2 FTS	Blank	0.001			
	0.01	0.012	8.01	6.86	116.7
	0.1	0.097	1.53	1.57	97.1
	1	0.888	0.23	0.26	88.8
10:2 FTS	Blank	0.000			
	0.01	0.011	10.94	10.09	108.4
	0.1	0.095	3.27	3.43	95.4
	1	0.862	1.22	1.42	86.2

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a high fat milk matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX

LCMS-8050RX

LCMS-8060RX

LCMS-2020

LCMS-2050

Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-168 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMSTM-8060NX High Performance Liquid Chromatograph Mass Spectrometer
NexeraTM series High Performance Liquid Chromatograph

Determination of 30 PFAS in Fish Oil by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Dominika Gruszecka¹, Eishi Imoto¹, Toshiya Matsubara¹

¹ Shimadzu Scientific Instruments, Inc.

User Benefits

- ◆ Validated method for 30 PFAS in fish oil meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in fish oil followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxA, PFNA, and PFOS.

In this study, we spiked samples at five concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.25
Perfluoropentanoic acid	PFPeA	2706-90-3	0.25
Perfluorohexanoic acid	PFHxA	307-24-4	0.25
Perfluoroheptanoic acid	PFHpA	375-85-9	0.25
Perfluorooctanoic acid	PFOA	335-67-1	0.25
Perfluorononanoic acid	PFNA	375-95-1	0.25
Perfluorodecanoic acid	PFDA	335-76-2	0.25
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.25
Perfluorododecanoic acid	PFDoA	307-55-1	0.25

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.25
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.25
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.25
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.25
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.25
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.25
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.25
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.25
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.25
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.25
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.25
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.25
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.25
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.25
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.25
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.25
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.25
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.25
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.25
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.25
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.5

■ Sample Preparation and Analysis Conditions

Fish oil soft gels were purchased locally. Samples were prepared by cutting the tip of the capsule with scissors and squeezing the oil out into a 50 ml tube. Test portions were spiked in triplicate at five different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram fish oil test portions spiked with concentrations of 0.1, 0.5, 1.5, 5.0 and 15.0 ng/g. Quantitation was carried out on additional fish oil samples spiked in triplicate at 0.25, 0.5, 1.0, 5.0, and 10 ng/g. Since standards were extracted in a fish oil matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram portions were weighed, spiked with target analytes and internal standards, and 10 mL of acetonitrile was added. The samples were vortexed for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.5 ng/g PFAS in a fish oil matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and five different concentrations ranging from below the SMPR required LOQ to 20 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, $S/N > 3$ for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a $S/N > 10$. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	ND			
	0.25	0.254	6.17	6.07	101.7
	0.5	0.434	7.21	8.31	86.8
	1.0	1.017	3.90	3.83	101.7
	5.0	5.338	7.41	6.94	106.8
	10	10.414	3.17	3.04	104.1
PFPeA	Blank	ND			
	0.25	0.249	2.87	2.88	99.7
	0.5	0.447	5.37	6.00	89.5
	1.0	1.025	2.20	2.15	102.5
	5.0	5.247	8.12	7.74	104.9
	10	10.238	2.98	2.91	102.4
PFHxA	Blank	0.008			
	0.25	0.264	4.76	4.50	105.6
	0.5	0.483	4.25	4.40	96.6
	1.0	1.103	1.17	1.06	110.3
	5.0	5.537	6.70	6.05	110.7
	10	10.836	5.22	4.81	108.4
PFHpA	Blank	ND			
	0.25	0.266	3.07	2.87	106.7
	0.5	0.482	4.35	4.51	96.4
	1.0	1.089	2.63	2.42	108.9
	5.0	5.556	8.68	7.81	111.1
	10	10.898	4.01	3.68	109.0
PFOA	Blank	0.009			
	0.25	0.270	3.15	2.92	108.0
	0.5	0.470	2.78	2.96	94.0
	1.0	1.022	1.14	1.11	102.2
	5.0	5.292	5.17	4.88	105.8
	10	10.508	1.59	1.51	105.0
PFNA	Blank	ND			
	0.25	0.261	1.91	1.83	104.3
	0.5	0.473	2.11	2.23	94.7
	1.0	1.054	2.91	2.76	105.4
	5.0	5.453	3.67	3.36	109.0
	10	10.879	1.39	1.27	108.8
PFDA	Blank	0.006			
	0.25	0.265	6.96	6.57	105.9
	0.5	0.491	3.98	4.05	98.2
	1.0	1.054	0.81	0.77	105.4
	5.0	5.568	3.46	3.11	111.4
	10	10.804	1.37	1.26	108.1
PFUnA	Blank	0.010			
	0.25	0.259	6.56	6.34	103.4
	0.5	0.529	4.65	4.39	105.9
	1.0	1.040	3.08	2.96	104.0
	5.0	5.453	4.99	4.58	109.0
	10	11.149	1.21	1.08	111.5
PFDoA	Blank	0.011			
	0.25	0.252	6.11	6.06	100.8
	0.5	0.552	2.93	2.65	110.5
	1.0	1.076	7.91	7.35	107.6
	5.0	5.845	4.88	4.17	116.9
	10	11.709	1.20	1.02	117.1
PFTrDA	Blank	0.019			
	0.25	0.235	18.57	19.79	93.8
	0.5	0.544	9.18	8.43	108.9
	1.0	1.064	1.30	1.22	106.4
	5.0	5.459	4.67	4.27	109.2
	10	10.772	1.80	1.67	107.8
PFTeDA	Blank	0.028			
	0.25	0.202	19.65	24.26	81.0
	0.5	0.472	2.76	2.93	94.3
	1.0	0.994	7.50	7.54	99.4
	5.0	5.695	7.96	6.99	113.9
	10	11.144	2.97	2.67	111.4
PFBS	Blank	ND			
	0.25	0.240	3.96	4.13	96.0
	0.5	0.467	6.35	6.79	93.5
	1.0	0.920	0.46	0.50	92.0
	5.0	4.573	4.12	4.51	91.4
	10	9.796	2.10	2.15	98.0

PFPeS	Blank	ND			
	0.25	0.276	7.43	6.72	110.5
	0.5	0.471	3.21	3.41	94.1
	1.0	1.079	4.42	4.10	107.9
	5.0	5.282	2.68	2.53	105.6
	10	11.338	4.10	3.62	113.4
PFHxS	Blank	0.013			
	0.25	0.262	4.08	3.89	104.9
	0.5	0.482	1.97	2.04	96.4
	1.0	0.941	1.99	2.12	94.1
	5.0	4.740	4.58	4.83	94.8
	10	9.739	0.49	0.51	97.4
PFHpS	Blank	ND			
	0.25	0.266	3.21	3.02	106.3
	0.5	0.464	5.18	5.59	92.8
	1.0	1.066	2.45	2.30	106.6
	5.0	5.037	0.35	0.35	100.7
	10	11.747	3.06	2.61	117.5
PFOS	Blank	0.008			
	0.25	0.242	5.29	5.46	96.8
	0.5	0.464	3.63	3.90	92.9
	1.0	0.922	3.99	4.33	92.2
	5.0	4.614	4.50	4.88	92.3
	10	9.836	2.56	2.60	98.3
PFNS	Blank	ND			
	0.25	0.256	5.54	5.41	102.4
	0.5	0.488	4.03	4.14	97.5
	1.0	1.011	3.78	3.74	101.1
	5.0	4.918	0.85	0.86	98.4
	10	11.269	1.63	1.44	112.7
PFDS	Blank	0.004			
	0.25	0.261	2.69	2.58	104.2
	0.5	0.480	7.88	8.21	96.0
	1.0	0.986	3.78	3.84	98.6
	5.0	4.884	2.71	2.77	97.7
	10	11.074	2.04	1.84	110.8
PFUnDS	Blank	0.038			
	0.25	0.229	11.65	12.74	91.4
	0.5	0.471	3.69	3.92	94.2
	1.0	0.936	5.08	5.43	93.6
	5.0	4.868	5.07	5.21	97.3
	10	9.633	2.03	2.10	96.3
PFDoS	Blank	0.014			
	0.25	0.240	17.49	18.23	95.9
	0.5	0.463	3.94	4.26	92.5
	1.0	1.008	8.40	8.34	100.8
	5.0	5.083	3.24	3.19	101.7
	10	10.273	3.84	3.74	102.7
PFTrDS	Blank	0.013			
	0.25	0.217	10.40	11.99	86.8
	0.5	0.449	2.72	3.03	89.9
	1.0	1.017	8.78	8.64	101.7
	5.0	5.471	5.86	5.36	109.4
	10	11.963	3.60	3.01	119.6
PFOSA	Blank	ND			
	0.25	0.254	5.16	5.08	101.7
	0.5	0.497	7.54	7.59	99.4
	1.0	0.998	6.46	6.48	99.8
	5.0	5.054	3.82	3.78	101.1
	10	10.768	3.10	2.88	107.7
9CI-PF3ONS	Blank	0.007			
	0.25	0.261	4.09	3.91	104.5
	0.5	0.500	1.76	1.76	100.1
	1.0	0.976	3.38	3.46	97.6
	5.0	4.989	5.16	5.17	99.8
	10	10.804	1.33	1.23	108.0
11CI-PF3OUdS	Blank	0.004			
	0.25	0.276	0.75	0.68	110.5
	0.5	0.496	3.36	3.39	99.3
	1.0	1.036	9.19	8.87	103.6
	5.0	5.483	4.11	3.75	109.7
	10	12.079	4.51	3.73	120.8

HFPO-DA	Blank	0.010			
	0.25	0.251	4.55	4.53	100.4
	0.5	0.469	7.34	7.81	93.9
	1.0	0.876	3.07	3.50	87.6
	5.0	4.538	4.21	4.63	90.8
	10	9.461	3.98	4.21	94.6
DONA	Blank	ND			
	0.25	0.244	3.92	4.01	97.7
	0.5	0.466	4.58	4.92	93.1
	1.0	1.026	1.17	1.14	102.6
	5.0	5.390	9.33	8.65	107.8
	10	10.210	3.70	3.62	102.1
4:2 FTS	Blank	ND			
	0.25	0.250	10.10	10.10	99.9
	0.5	0.480	2.93	3.05	96.0
	1.0	1.107	7.81	7.05	110.7
	5.0	5.443	9.61	8.83	108.9
	10	10.788	7.37	6.83	107.9
6:2 FTS	Blank	ND			
	0.25	0.267	12.51	11.70	106.9
	0.5	0.511	13.05	12.78	102.1
	1.0	1.079	9.92	9.19	107.9
	5.0	5.342	8.95	8.38	106.8
	10	10.826	9.54	8.81	108.3
8:2 FTS	Blank	0.009			
	0.25	0.225	9.23	10.25	90.1
	0.5	0.483	0.38	0.39	96.6
	1.0	1.016	0.85	0.84	101.6
	5.0	5.945	10.53	8.86	118.9
	10	10.140	3.49	3.44	101.4
10:2 FTS	Blank	0.003			
	0.5	0.478	10.34	10.82	95.6
	1.0	1.060	8.61	8.13	106.0
	5.0	5.939	9.65	8.13	118.8
	10	11.502	2.70	2.35	115.0

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a fish oil matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX

LCMS-8050RX

LCMS-8060RX

LCMS-2020

LCMS-2050

Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-166 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMS™-8060NX High Performance Liquid Chromatograph Mass Spectrometer
Nexera™ series High Performance Liquid Chromatograph

Determination of 30 PFAS in Coffee by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Dominika Gruszecka¹, Toshiya Matsubara¹

¹ Shimadzu Scientific Instruments, Inc.

User Benefits

- ◆ Validated method for 30 PFAS in Coffee meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, with SPE cleanup
- ◆ Simple and rapid extraction

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in brewed coffee followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxA, PFNA, and PFOS.

In this study, we spiked samples at four concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.055
Perfluoropentanoic acid	PFPeA	2706-90-3	0.0055
Perfluorohexanoic acid	PFHxA	307-24-4	0.055
Perfluoroheptanoic acid	PFHpA	375-85-9	0.0055
Perfluorooctanoic acid	PFOA	335-67-1	0.0055
Perfluorononanoic acid	PFNA	375-95-1	0.0055
Perfluorodecanoic acid	PFDA	335-76-2	0.055
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.0055
Perfluorododecanoic acid	PFDoA	307-55-1	0.055

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.055
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.055
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.0055
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.0055
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.0055
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.055
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.055
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.055
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.055
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.055
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.055
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.55
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.55
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9CI-PF3ONS	756426-58-1	0.0055
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11CI-PF3OUdS	763051-92-9	0.055
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.055
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.0055
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.055
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.0055
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.55
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.055

■ Sample Preparation and Analysis Conditions

Organic brewed coffee was purchased locally. Samples were prepared by shaking the bottle and removing aliquots directly from the original container. Test portions were spiked in triplicate at four different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using 10-gram test portions spiked with concentrations of 0.001, 0.01, 0.1, 1.0 and 10.0 ng/g. Quantitation was carried out on additional samples spiked in triplicate at 0.0055, 0.055, 0.55 and 5.5 ng/g. Since standards were extracted in a brewed coffee matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

10-gram portions were weighed, spiked with target analytes and internal standards, and 10 mL of acetonitrile was added. The sample was shaken for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. For greater sensitivity, the extract was concentrated to dryness and dissolved in 0.4 mL of a methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.55 ng/g PFAS in a coffee matrix with separation of all peaks in nine minutes

MRM chromatogram




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₃ -PFBS
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₂ -PFDoA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve




Figure 4: PFHxS Calibration Curve




Figure 5: PFNA Calibration Curve





Figure 6: PFOS Calibration Curve



Blank matrixes and at four different concentrations ranging from below the SMPR required LOQ to approximately 20 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that ¹³C₂ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3.

The LOQs for all compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery%
PFBA	Blank	ND			
	0.055	0.049	4.19	4.75	88.1
	0.55	0.526	2.00	2.09	95.8
	5.5	5.553	0.31	0.30	101.0
PFPeA	Blank	0.001			
	0.0055	0.005	8.23	8.41	97.9
	0.055	0.051	3.42	3.67	93.1
	0.55	0.504	1.34	1.46	91.7
PFHxA	5.5	5.274	1.07	1.12	95.9
	Blank	0.001			
	0.0055	0.005	17.58	18.38	95.7
	0.055	0.054	1.08	1.11	97.2
PFHpA	0.55	0.534	3.22	3.32	97.2
	5.5	5.744	1.98	1.89	104.4
	Blank	0.001			
	0.0055	0.005	16.77	17.85	94.0
PFHpA	0.055	0.053	1.46	1.53	96.0
	0.55	0.530	2.95	3.07	96.3
	5.5	5.549	2.10	2.08	100.9
	Blank	ND			
PFOA	0.0055	0.005	2.36	2.62	89.9
	0.055	0.053	0.55	0.58	95.7
	0.55	0.555	0.67	0.66	100.9
	5.5	5.798	2.46	2.33	105.4
PFNA	Blank	0.000			
	0.0055	0.005	2.91	3.48	83.4
	0.055	0.052	4.07	4.26	95.4
	0.55	0.523	4.73	4.97	95.2
PFDA	5.5	5.581	3.29	3.24	101.5
	Blank	0.000			
	0.0055	0.004	14.51	17.41	83.4
	0.055	0.051	3.33	3.56	93.6
PFUnA	0.55	0.521	2.63	2.78	94.7
	5.5	5.445	2.10	2.12	99.0
	Blank	0.001			
	0.0055	0.005	10.28	11.12	92.4
PFDoA	0.055	0.052	3.02	3.22	93.8
	0.55	0.520	0.78	0.82	94.6
	5.5	5.420	4.17	4.23	98.5
	Blank	0.001			
PFTrDA	0.055	0.056	11.37	11.21	101.4
	0.55	0.504	1.40	1.53	91.7
	5.5	5.587	3.54	3.48	101.6
	Blank	0.001			
PFTeDA	0.0055	0.005	13.20	14.63	90.2
	0.055	0.043	3.55	4.53	78.2
	0.55	0.450	1.87	2.29	81.9
	5.5	5.272	3.37	3.51	95.9
PFBS	Blank	ND			
	0.055	0.037	2.40	3.59	66.9
	0.55	0.440	2.75	3.43	80.1
	5.5	5.363	4.01	4.12	97.5
PFPeS	Blank	0.000			
	0.0055	0.005	10.02	10.89	92.0
	0.055	0.051	1.92	2.07	92.7
	0.55	0.515	1.65	1.76	93.5
PFHxS	5.5	5.359	0.78	0.80	97.4
	Blank	0.000			
	0.0055	0.006	9.29	9.22	100.7
	0.055	0.060	2.99	2.75	108.9
PFHpS	0.55	0.547	1.80	1.81	99.5
	5.5	5.546	2.99	2.96	100.9
	Blank	ND			
	0.0055	0.006	12.64	12.35	102.3
PFHpS	0.055	0.054	1.14	1.16	97.8
	0.55	0.538	1.75	1.79	97.8
	5.5	5.752	1.80	1.72	104.6
	Blank	0.001			
PFHpS	0.055	0.052	6.88	7.25	94.9
	0.55	0.484	1.62	1.84	88.1
	5.5	4.930	1.72	1.91	89.6
	Blank	0.001			

PFOS	Blank	0.000			
	0.0055	0.005	14.22	15.79	90.1
	0.055	0.048	1.11	1.26	88.1
	0.55	0.493	2.30	2.57	89.7
PFNS	5.5	5.395	7.45	7.59	98.1
	Blank	0.001			
	0.0055	0.005	18.55	20.76	89.3
	0.055	0.056	1.30	1.28	101.6
PFDS	0.55	0.573	4.06	3.89	104.2
	5.5	5.419	4.92	4.99	98.5
	Blank	0.007			
	0.055	0.051	8.80	9.47	93.0
PFUnDS	0.55	0.497	2.55	2.83	90.3
	5.5	5.096	4.15	4.48	92.7
	Blank	0.005			
	0.055	0.047	10.89	12.89	84.5
PFDoS	0.55	0.401	4.99	6.85	72.9
	5.5	5.537	9.92	9.85	100.7
	Blank	0.003			
	0.055	0.047	13.08	15.26	85.7
PFTrDS	0.55	0.478	5.59	6.44	86.8
	5.5	4.933	1.55	1.73	89.7
	Blank	0.001			
	0.055	0.038	6.92	10.10	68.5
PFOSA	0.55	0.437	3.70	4.65	79.5
	5.5	5.021	3.90	4.27	91.3
	Blank	ND			
	0.055	0.046	12.00	14.31	83.9
9CI-PF3ONS	0.55	0.605	7.45	6.77	110.0
	5.5	4.924	3.98	4.45	89.5
	Blank	0.000			
	0.0055	0.006	13.50	13.00	103.8
11CI-PF3OUdS	0.055	0.054	0.32	0.33	98.0
	0.55	0.538	1.85	1.89	97.8
	5.5	5.888	7.98	7.45	107.0
	Blank	0.003			
HFPO-DA	0.0055	0.007	13.02	11.27	115.5
	0.055	0.047	9.48	11.20	84.6
	0.55	0.451	1.20	1.46	82.0
	5.5	4.916	3.70	4.14	89.4
DONA	Blank	ND			
	0.055	0.046	2.80	3.33	84.3
	0.55	0.550	3.03	3.03	99.9
	5.5	5.770	2.77	2.64	104.9
4:2 FTS	Blank	0.000			
	0.0055	0.006	3.27	3.18	102.8
	0.055	0.053	3.57	3.69	96.7
	0.55	0.547	2.87	2.89	99.4
6:2 FTS	5.5	5.721	1.55	1.49	104.0
	Blank	0.005			
	0.055	0.051	6.44	6.97	92.5
	0.55	0.557	1.34	1.33	101.2
8:2 FTS	5.5	5.534	1.07	1.06	100.6
	Blank	0.001			
	0.0055	0.005	6.71	7.33	91.5
	0.055	0.056	16.30	16.11	101.2
10:2 FTS	0.55	0.559	10.64	10.46	101.7
	5.5	5.406	4.45	4.53	98.3
	Blank	0.000			
	0.0055	0.006	19.62	20.35	96.4
4:2 FTS	0.055	0.054	8.94	9.00	99.3
	0.55	0.609	2.99	2.69	110.9
	5.5	7.143	6.46	4.97	129.9
	Blank	0.001			
6:2 FTS	0.055	0.043	10.17	12.82	79.4
	0.55	0.443	3.81	4.73	80.4
	5.5	5.042	6.27	6.84	91.7

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

■ Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a coffee matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

■ Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX



LCMS-8050RX



LCMS-8060RX



LCMS-2020



LCMS-2050



Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-164 First Edition: September 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

LCMS™-8060NX High Performance Liquid Chromatograph Mass Spectrometer
Nexera™ series High Performance Liquid Chromatograph

Determination of 30 PFAS in Seafood by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS/MS)

William Lipps¹, Toshiya Matsubara¹, Dominika Gruszecka¹

¹ Shimadzu Scientific Instruments, Inc.

User Benefits

- ◆ Validated method for 30 PFAS in seafood (shrimp) meeting all criteria of AOAC SMPR 2023.003
- ◆ High precision, excellent recovery, low Limit of Quantification (LOQ)
- ◆ Simple and rapid extraction using QuEChERS

Introduction

Per- and polyfluorinated substances (PFAS) are a diverse group of man-made chemicals used in numerous products since the 1950s. PFAS can enter the food supply by contact in environmentally contaminated areas, during food processing, or exposure to packaging. Because PFAS have been linked to serious health effects, accurate methodology is needed. In this application news, we describe a single laboratory validation study with a rapid extraction of low concentrations of 30 PFAS in shrimp using the QuEChERS technique followed by analysis using the Shimadzu Nexera Liquid Chromatograph coupled to a Shimadzu LCMS-8060NX triple quadrupole mass spectrometer (Figure 1).

We optimized the chromatography and instrument operating parameters to achieve excellent peak shape, separation, and sensitivity. Sensitivity was improved for 4 PFAS, PFOA, PFHxS, PFNA, and PFOS.

In this study, we spiked samples at three concentrations in triplicate. For greater accuracy, standards were matrix-matched and extracted and spikes were quantified using the isotope dilution technique. Recovery and precision were compared to the requirements of AOAC SMPR 2023.003. In addition, we determined the Limit of Quantitation (LOQ) as the lowest concentration meeting accuracy and precision, ion ratio, retention time, and signal-to-noise ratio criteria of the qualifier ion. All recovery, precision, and LOQ's met the acceptance criteria of the SMPR. The target analytes, their acronym, chemical abstract number, and experimentally determined LOQ are shown in Table 1.



Figure 1: Nexera™ and LCMS™-8060NX. The ion focus design improves signal intensity with higher gas flows and higher effective temperatures.

Table 1: PFAS Analytes, Acronyms, CAS No. and Method LOQ

Analyte name	Acronym	CAS No.	LOQ (ppb)
Perfluorobutanoic acid	PFBA	375-22-4	0.055
Perfluoropentanoic acid	PFPeA	2706-90-3	0.0055
Perfluorohexanoic acid	PFHxA	307-24-4	0.0055
Perfluoroheptanoic acid	PFHpA	375-85-9	0.055
Perfluorooctanoic acid	PFOA	335-67-1	0.0055
Perfluorononanoic acid	PFNA	375-95-1	0.0055
Perfluorodecanoic acid	PFDA	335-76-2	0.0055
Perfluoroundecanoic acid	PFUnA	2058-94-8	0.0055
Perfluorododecanoic acid	PFDoA	307-55-1	0.055

Perfluorotridecanoic acid	PFTrDA	72629-94-8	0.055
Perfluorotetradecanoic acid	PFTeDA	376-06-7	0.055
Perfluorobutanesulfonic acid	PFBS	375-73-5	0.055
Perfluoropentanesulfonic acid	PFPeS	2706-91-4	0.055
Perfluorohexanesulfonic acid	PFHxS	355-46-4	0.0055
Perfluoroheptanesulfonic acid	PFHpS	375-92-8	0.0055
Perfluorooctanesulfonic acid	PFOS	1763-23-1	0.055
Perfluorononanesulfonic acid	PFNS	68259-12-1	0.0055
Perfluorodecanesulfonic acid	PFDS	335-77-3	0.0055
Perfluoroundecanesulfonic acid	PFUnDS	749786-16-1	0.0055
Perfluorododecanesulfonic acid	PFDoS	79780-39-5	0.0055
Perfluorotridecanesulfonic acid	PFTrDS	791563-89-8	0.055
Perfluorooctanesulfonamide	PFOSA	754-91-6	0.055
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9Cl-PF3ONS	756426-58-1	0.0055
11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11Cl-PF3OUdS	763051-92-9	0.0055
Hexafluoropropylene oxide dimer acid	HFPO-DA	13252-13-6	0.055
4,8-Dioxa-3H-perfluorononanoic acid	DONA	919005-14-4	0.0055
1H,1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2 FTS	757124-72-4	0.055
1H,1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2 FTS	27619-97-2	0.55
1H,1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2 FTS	39108-34-4	0.0055
1H,1H, 2H, 2H-Perfluorododecane sulfonic acid	10:2 FTS	120226-60-0	0.0055

■ Sample Preparation and Analysis Conditions

Samples were prepared by spiking prepared shrimp samples in triplicate at three different concentrations with 30 native PFAS (Table 1) and 16 isotopically labeled internal standards. Calibration curves for use in the quantitative analysis were prepared using shrimp test portions spiked with concentrations of 0.001, 0.01, 0.10, 1.0, and 10.0 ng/g. Quantitation was carried out on additional shrimp samples spiked in triplicate at 0.0055, 0.055, 0.55, and 5.5 ng/g. Since standards were extracted in a shrimp matrix and carried through the same procedure, the final concentration of each PFAS in the sample can be calculated directly from the curve.

Frozen shrimp were thawed at room temperature and then crushed in a grinder with dry ice for 30 seconds at 4000 rpm. The ground sample was stored in a freezer overnight to remove all the dry ice. Ten-gram portions were weighed, and 10 mL of acetonitrile was added. The samples were vortexed for 1 minute and a QuEChERS packet was added. The sample was shaken for 1 minute and then centrifuged for 5 minutes at 4000 rpm. An aliquot of the acetonitrile layer was transferred to a tube and diluted 5 times with PFAS-free reagent water. The sample was then passed through a weak anion exchange (WAX) Solid Phase Extraction (SPE) cartridge and the PFAS were eluted with basic methanol. For greater sensitivity, 2 mL of the extract was concentrated in 400 µL of a methanol-water mixture.

A volume suitable to obtain the required sensitivity of the extract was injected onto a UHPLC system (Shimadzu Nexera). Adequate separation of all compounds was achieved in nine minutes. (Figure 2 chromatogram shows the separation of all peaks).

For this study, Shimadzu evaluated 1984 different instrument settings, and 6 different column and gradient combinations, to achieve excellent peak shape and resolution between peaks, as well as to maximize the signal-to-noise ratio of PFOA, PFHxS, PFNA, and PFOS. Mass spectrometry was performed on a Shimadzu LCMS-8060NX with heated electrospray ionization operated in negative mode. Specific compound MRM transitions and associated internal standards are listed in Table 2. Chromatography was adjusted to provide sufficient separation of PFOA from potential cholic acid interferences, and to provide baseline resolution of branched and linear isomers (Figure 3).




Figure 2: Chromatogram of 0.55 ng/g PFAS in shrimp matrix with separation of all peaks in nine minutes




Figure 3: Separation of PFOS from Cholic Acids and baseline resolution between PFOS branched and linear isomers

Table 2: MRM Transitions and Internal Standard Associations

Analyte	Quantitation Ion	Qualifier Ion	Internal Standard
PFBA	212.9 > 168.8		¹³ C ₄ -PFBA
PFPeA	262.9 > 218.8		¹³ C ₅ -PFPeA
PFHxA	313 > 268.8	313 > 118.9	¹³ C ₅ -PFHxA
PFHpA	362.9 > 318.9	362.9 > 168.8	¹³ C ₄ -PFHpA
PFOA	412.9 > 368.9	412.9 > 168.8	¹³ C ₈ -PFOA
PFNA	462.9 > 418.5	462.9 > 218.6	¹³ C ₉ -PFNA
PFDA	512.9 > 468.95	512.9 > 268.55	¹³ C ₆ -PFDA
PFUnA	562.9 > 518.95	562.9 > 269	¹³ C ₇ -PFUnA
PFDoA	612.9 > 268.6	612.9 > 168.6	¹³ C ₇ -PFUnA
PFTTrDA	663 > 618.6	663 > 168.5	¹³ C ₇ -PFUnA
PFTeDA	713 > 669	713 > 168.5	¹³ C ₇ -PFUnA
PFBS	298.8 > 79.8	298.8 > 98.9	¹³ C ₅ -PFHxA
PFPeS	348.8 > 79.8	348.8 > 98.9	¹³ C ₄ -PFHpA
PFHxS	398.8 > 79.8	398.8 > 98.9	¹³ C ₃ -PFHxS
PFHpS	448.8 > 79.8	448.8 > 99	¹³ C ₆ -PFDA
PFOS	498.8 > 79.8	498.8 > 98.9	¹³ C ₈ -PFOS
PFNS	548.9 > 79.8	548.8 > 99	¹³ C ₇ -PFUnA
PFDS	598.8 > 99	598.8 > 79.8	¹³ C ₇ -PFUnA
PFUnDS	649 > 99	649 > 80	¹³ C ₈ -PFOS
PFDoS	699 > 80	699 > 99	¹³ C ₂ -PFDoA
PFTTrDS	749 > 80	749 > 279.6	¹³ C ₇ -PFUnA
PFOSA	498 > 78	498 > 477.95	¹³ C ₈ -PFOSA
9Cl-PF3ONS	530.8 > 350.8	532.9 > 352.95	¹³ C ₈ -PFOS
11Cl-PF3OUdS	631 > 451	632.9 > 452.95	¹³ C ₇ -PFUnA
HFPO-DA	284.8 > 168.8	284.8 > 118.9	¹³ C ₃ -HFPO-DA
DONA	376.9 > 250.8	376.9 > 84.9	¹³ C ₄ -PFHpA
4:2 FTS	326.9 > 80.9	326.9 > 306.8	¹³ C ₂ -4:2 FTS
6:2 FTS	426.9 > 80.8	426.9 > 406.9	¹³ C ₂ -6:2 FTS
8:2 FTS	526.9 > 506.45	526.9 > 80.55	¹³ C ₃ -PFHxS
10:2 FTS	626.9 > 606.7	626.9 > 80.9	¹³ C ₇ -PFUnA

Quantitative Analysis

Calibration standards were processed the same as samples. A linear model not forced through zero isotopic dilution calibration in matrix-matched standards provided the best fit and best recoveries of analytes. Residuals of each point in the curve were $\pm 25\%$ of the expected value. Calibration curves for PFOA, PFHxS, PFNA, and PFOS are shown in Figures 3 – 6 respectively. Branched and linear isomers of PFHxS and PFOS were integrated together.

Figure 3: PFOA Calibration Curve





Figure 5: PFNA Calibration Curve



Blank matrixes and at least three different concentrations ranging from below the SMPR required LOQ to approximately 20 times the estimated LOQ were analyzed in triplicate. Recovery and repeatability for each analyte at each concentration are given in Table 3. The LOQ for each analyte was estimated by spiking at concentrations at, or below, the required LOQs listed in SMPR-2023_003. The spiked samples were analyzed in triplicate and the mean and repeatability standard deviation were calculated. Then, the standard deviation was divided by the mean and multiplied by 100% to calculate the repeatability percent relative standard deviation (RSD).

Exact labeled analogs were used as isotope dilution standards, except in a few cases where we found that ¹³C₂ suffered matrix interferences. In these cases, another non-interfering isotope was chosen. The analytes tested and the associated calibration isotopes are shown in Table 3. During the single laboratory validation, the shrimp matrix contained an unknown interference coeluting with the ¹³C₃-PFBS isotope so ¹³C₅-PFHxA was used instead.

Figure 4: PFHxS Calibration Curve





Figure 6: PFOS Calibration Curve



The LOQs for all matrixes and compounds were determined using an Excel worksheet that compared each of the requirements of the SMPR including retention time, recovery, repeatability, S/N > 3 for the qualifier ion and an ion ratio of $\pm 30\%$. PFBA, PFPeA, and PFOSA LOQ were set at the minimum concentration, meeting recovery and repeatability requirements and a S/N > 10. The lowest concentration to meet all the requirements of the SMPR was set as the LOQ. Figure 7 shows examples of the LOQ spike for PFHxS, PFNA, PFOA, and PFOS and their corresponding internal standards.

Table 3: Recovery and repeatability for each analyte at each spike concentration

Analyte	Spike conc. (ppb)	Average conc. (ppb)	Standard Deviation	%RSD	Average Recovery (%)
PFBA	Blank	0.007			
	0.055	0.056	0.98	0.98	100.3
	0.55	0.518	0.32	0.34	94.1
	5.5	5.478	0.35	0.35	99.6
PFPeA	Blank	ND			
	0.0055	0.005	3.52	3.81	92.6
	0.055	0.054	2.17	2.23	97.6
	0.55	0.517	0.67	0.71	94.1
PFHxA	Blank	0.001			
	0.0055	0.005	9.08	9.50	95.6
	0.055	0.057	1.48	1.43	103.7
	0.55	0.538	0.66	0.67	97.8
PFHpA	Blank	0.001			
	0.0055	0.005	3.90	4.06	96.1
	0.055	0.051	0.75	0.82	92.0
	0.55	0.518	0.26	0.28	94.2
PFOA	Blank	0.001			
	0.0055	0.005	3.90	4.06	96.1
	0.055	0.051	0.75	0.82	92.0
	0.55	0.518	0.26	0.28	94.2
PFNA	Blank	0.000			
	0.0055	0.006	7.32	6.77	108.1
	0.055	0.056	0.68	0.67	102.1
	0.55	0.583	0.55	0.52	106.0
PFDA	Blank	0.000			
	0.0055	0.006	4.45	4.38	101.6
	0.055	0.056	0.72	0.71	101.9
	0.55	0.572	0.86	0.83	104.0
PFUnA	Blank	ND			
	0.0055	0.007	6.07	4.92	123.3
	0.055	0.056	2.29	2.27	100.9
	0.55	0.579	1.03	0.98	105.2
PFDoA	Blank	0.001			
	0.055	0.054	6.52	6.67	97.7
	0.55	0.528	3.51	3.66	96.0
	5.5	5.863	3.10	2.91	106.6
PFTTrDA	Blank	ND			
	0.055	0.063	2.31	2.01	115.0
	0.55	0.586	3.89	3.65	106.6
	5.5	6.215	2.56	2.26	113.0
PFTeDA	Blank	0.001			
	0.055	0.055	5.22	5.23	99.7
	0.55	0.555	4.27	4.23	100.8
	5.5	6.313	3.03	2.64	114.8
PFBS	Blank	0.001			
	0.055	0.051	5.55	5.93	93.6
	0.55	0.527	2.83	2.95	95.8
	5.5	5.940	2.31	2.14	108.0
PFPeS	Blank	0.002			
	0.0055	0.004	15.88	20.67	76.8
	0.055	0.055	3.92	3.89	100.7
	0.55	0.598	3.29	3.02	108.8
	5.5	5.713	0.71	0.68	103.9

PFHxS	Blank	ND			
	0.0055	0.004	14.59	17.62	82.8
	0.055	0.052	2.32	2.45	94.9
	0.55	0.535	1.06	1.09	97.2
PFHpS	5.5	6.070	0.25	0.23	110.4
	Blank	0.001			
	0.0055	0.006	20.96	18.27	114.8
	0.055	0.057	9.51	9.14	104.0
PFOS	0.55	0.590	8.60	8.02	107.2
	5.5	6.223	4.30	3.80	113.1
	Blank	0.001			
	0.055	0.053	5.12	5.28	96.9
PFNS	0.55	0.517	1.77	1.88	94.0
	5.5	5.700	1.03	0.99	103.6
	Blank	0.001			
	0.0055	0.006	5.17	4.42	117.0
PFDS	0.055	0.052	3.68	3.86	95.1
	0.55	0.526	3.20	3.35	95.6
	5.5	6.117	4.19	3.76	111.2
	Blank	0.000			
PFUnDS	0.0055	0.006	9.46	8.11	116.6
	0.055	0.056	0.87	0.85	102.6
	0.55	0.569	4.82	4.67	103.4
	5.5	6.148	4.56	4.08	111.8
PFDoS	Blank	ND			
	0.0055	0.004	15.95	19.84	80.4
	0.055	0.056	5.75	5.66	101.5
	0.55	0.594	0.61	0.56	107.9
PFTrDS	5.5	5.614	0.85	0.84	102.1
	Blank	0.000			
	0.0055	0.007	14.15	11.90	118.9
	0.055	0.056	5.90	5.80	101.7
PFOSA	0.55	0.607	4.68	4.24	110.3
	5.5	5.949	2.90	2.68	108.2
	Blank	0.000			
	0.0055	0.005	8.10	7.94	102.0
9CI-PF3ONS	0.055	0.057	4.16	4.03	103.2
	0.55	0.581	5.25	4.97	105.7
	5.5	6.228	3.96	3.50	113.2
	Blank	0.001			
11CI-PF3OUdS	0.0055	0.006	9.45	8.91	106.1
	0.055	0.055	3.80	3.81	99.8
	0.55	0.556	0.78	0.77	101.1
	5.5	5.973	1.75	1.61	108.6
HFPO-DA	Blank	0.000			
	0.0055	0.004	7.22	8.57	84.3
	0.055	0.054	1.65	1.67	98.6
	0.55	0.570	0.78	0.75	103.6
PF3ONS	5.5	5.896	2.62	2.44	107.2
	Blank	0.001			
	0.0055	0.006	5.66	5.06	112.0
	0.055	0.054	1.81	1.85	97.8
PF3OUdS	0.55	0.546	3.03	3.06	99.2
	5.5	6.022	3.61	3.30	109.5
	Blank	0.002			
	0.0055	0.005	6.55	6.59	99.4
PF3OUs	0.055	0.053	2.17	2.24	96.9
	0.55	0.517	3.93	4.18	94.0
	5.5	5.653	0.95	0.93	102.8
	Blank	0.002			

DONA	Blank	0.000			
	0.0055	0.005	2.04	2.07	98.7
	0.055	0.060	3.96	3.59	110.1
	0.55	0.619	2.37	2.10	112.6
	5.5	5.714	1.95	1.87	103.9
4:2 FTS	Blank	ND			
	0.055	0.052	1.92	2.05	93.6
	0.55	0.535	0.35	0.36	97.2
	5.5	5.663	0.80	0.78	103.0
6:2 FTS	Blank	0.018			
	0.055	0.048	4.12	4.66	88.4
	0.55	0.495	1.60	1.78	90.0
	5.5	5.290	1.26	1.31	96.2
8:2 FTS	Blank	0.000			
	0.0055	0.005	5.61	6.54	85.8
	0.055	0.060	1.83	1.67	110.1
	0.55	0.609	3.59	3.25	110.6
	5.5	5.790	0.83	0.79	105.3
10:2 FTS	Blank	0.000			
	0.0055	0.006	13.17	12.02	109.5
	0.055	0.065	2.75	2.34	117.8
	0.55	0.641	2.69	2.31	116.6
	5.5	6.102	4.75	4.28	110.9

ND = average results less than zero




Figure 7: LOQ Peaks with Internal Standards

Conclusion

The Shimadzu LCMS-8060NX Triple Quadrupole Mass Spectrometer coupled with a Shimadzu Nexera Liquid Chromatograph was used in a single laboratory study to measure 30 PFAS compounds in a shrimp matrix and compared to criteria set by AOAC SMPR 2023.003. Chromatography conditions and the mass spectrometer were optimized to achieve excellent separation of all analytes, baseline resolution between linear and branched isomers, and a two-minute separation between PFOS and potentially interfering cholic acids.

Precision and recovery and the experimentally determined LOQ are well within the requirements of the SMPR.

Reference

- 1) AOAC SMPR 2023.003

UAFMS

ULTRA FAST MASS SPECTROMETRY



LCMS-8045RX



LCMS-8050RX



LCMS-8060RX



LCMS-2020



LCMS-2050



Q-TOF LCMS-9030/9050

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. Established in 1975, Shimadzu Scientific Instruments (SSI), the American subsidiary of Shimadzu Corporation, provides a comprehensive range of analytical solutions to laboratories throughout North, Central, and parts of South America. SSI maintains a network of ten regional offices strategically located across the United States, with experienced technical specialists, service and sales engineers situated throughout the country, as well as applications laboratories on both coasts.

For information about Shimadzu Scientific Instruments and to contact your local office, please visit our Web site at www.ssi.shimadzu.com

02-SSI-LCMS-157 First Edition: May 2024, Revised July 2024



Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

References

1. Shoemaker, J., Grimmet, P.E, and Boutin, P.K. Determination of Selected Perfluorinated Alkyl Acids in Drinking Water by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS) Method 537 Ver. 1.1. National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency; Cincinnati OH, 2009.
2. Shoemaker, J. and Dan Tettenhorst. Method 537.1: Determination of Selected Per- and Polyfluorinated Alkyl Substances in Drinking Water by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS). U.S. Environmental Protection Agency, Office of Research and Development, National Center for Environmental Assessment, Washington, DC, 2018.
3. U.S. Environmental Protection Agency, Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous, Solid, Biosolids, and Tissue Samples by LC-MS/MS, Method 1633, Office of Water, Office of Science and Technology Engineering and Analysis Division 1200 Pennsylvania Avenue, NW Washington, DC 20460, January 2024
4. ASTM International, D7979, Standard Test Method for Determination of Per- and Polyfluoroalkyl Substances in Water, Sludge, Influent, Effluent and Wastewater by Liquid Chromatography Tandem Mass Spectrometry (LC/MS/MS), 2017, DOI: 10.1520/D7979-17
5. ASTM International, D8421, Standard Test Method for Determination of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous Matrices by Co-solvation followed by Liquid Chromatography Tandem Mass Spectrometry (LC/MS/MS), 2021, 10.1520/D8421-21
6. ASTM International, D7968a, Standard Test Method for Determination of Polyfluorinated Compounds in Soil by Liquid Chromatography Tandem Mass Spectrometry (LC/MS/MS), 2017, DOI: 10.1520/D7968-17A
7. ASTM International, D8535, Standard Test Method for Determination of Per- and Polyfluoroalkyl Substances (PFAS) in Soil/Biosolid Matrices by Solvent Extraction, Filtering, and Followed by Liquid Chromatography Tandem Mass Spectrometry (LC/MS/MS), 2023, DOI: 10.1520/D8535-23
8. Genualdi, S. and deJager, L., Method C-010.02, Determination of 16 Per and Polyfluoroalkyl Substances (PFAS) in Processed Food using Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS), U.S Food and Drug Administration, 2021
9. United States Department of Agriculture, CLG - PFAS 2.03, Screening, Determination and Confirmation of PFAS by UPLC-MS-MS, Food Safety and Inspection Service, Office of Public Health Science, Washington DC, 2021
10. Genualdi, S and deJager L., Method C-010.03, Determination of 30 Per and Polyfluoroalkyl Substances (PFAS) in Food and Feed using Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS), U.S. Food and Drug Administration, 2024
11. <https://www.aoac.org/resources/aoac-smpr-2023-003> [accessed 11 Mar 2024]



Shimadzu Corporation

www.shimadzu.com/an/

Shimadzu Scientific Instruments

7102 Riverwood Drive, Columbia, Maryland 21046, U.S.A.
Phone: 800-477-1227/410-381-1227, Fax: 410-381-1222

www.ssi.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.