

Introduction

Labs that test cannabis and hemp typically work in a high throughput environment and receive samples that include plant material, extracts, tinctures, inhaled products, and a wide variety of edibles. Typically, these samples are batched by sample type and then run as batches in a microwave digestion system. Previously, the only way to batch all of the sample types listed above required a high pressure device which are very expensive.

CEM has developed a method to prepare all of these samples in a single batch in a MARS 6 microwave digestion system equipped with a simple, three-part MARSXpress vessel. Up to forty mixed samples can be prepared in a single batch. It provides the highest throughput and most cost effective solution available to laboratories. For this study, eleven different sample types, of all varieties, were digested in a single batch, in duplicate. All together 24 samples and blanks were run simultaneously in a single batch for maximum throughput.

Experimental Conditions

Samples: Eleven samples were chosen to simulate the wide variety of samples that are tested in the cannabis lab. These were sourced from local hemp dispensaries. The samples are shown in this photograph and include: hemp flower, hemp oil, MCT oil, topical cream, crude CBD oil, beef jerky, peanut butter, ghee (butter), granola bar, gummy snacks, hard candy.





Figure 1. Cannabis and hemp containing samples

Table 1. MARS 6 Cannabis One Touch Method Conditions

Parameter	Setting
Power	1800 W
Ramp Time	20 min
Hold Time	10 min
Temperature	210 °C



Sample Preparation: 0.5 grams of each sample were digested in a MARS 6 with MARSXpress Plus vessels in order to meet the sampling requirements for California heavy metals testing. An acid mixture of 9 mL nitric acid and 1 mL HCl was then added and the vessels were sealed. A pre-programmed Cannabis method was selected on the MARS 6 touchscreen interface. This method was developed to heat the samples to 210 °C and hold them for a specified time in order to achieve a complete digestion of all sample types. The samples were then diluted to 50 mL volume with DI water and analyzed via ICP-MS.



MARS 6 with MARSXpress Plus vessels

Experimental Conditions

Table 2. ICP-MS 2030 Operating Conditions

Parameter	Setting
Radio Frequency Power	1.20 kW
Sampling Depth	5.0 mm
Plasma Gas	8.0 L/min
Auxiliary Gas	1.10 L/min
Carrier Gas	0.70 L/min
Mix Gas	0.00 L/min
Cell Gas	6.0 mL/min
Cell Voltage	-21 V
Energy Filter	7.0 V
Chamber Temperature	5 °C



Shimadzu ICP-S-2030M

Results and Discussion

Figure 3 illustrates the ramping of the varied sample types. Note that during the ramp certain samples are higher in temperatures than others. These samples are likely higher in sugars or carbohydrates as these bonds break first. The MARS 6 with iWave technology and advanced software makes continuous adjustments to modulate power in order to maintain the proper heating conditions. The robustness of the high throughput MARSXpress Plus vessel also contributes to the ability to run mixed samples in one batch. In Figure 4 the samples are in the hold phase and the temperatures of the varied samples have equilibrated. Figure 5 shows a subset of samples that have been completely digested and diluted to volume.




Figure 3. Temperature of all vessels during ramp




Figure 4. Temperature of all vessels during hold



Figure 5. Diluted samples

ICP-MS Analysis

Table 3. Instrument Measurement Parameters

Element	Mass	Internal Standard	Cell Gas	Cal. Range (ppb)	Cal. R ²	Detection Limit (ppb)	Integ. Time (sec)	No. of Scan	Repeat No.	Spiked Concn. (ppb) ^a
As	75	Sc (45)	On	0.6 - 10	0.99946	0.00818	2.0	10	3	5
Cd	111	In (115)	Off	0.6 - 10	0.99992	0.00131	2.0	10	3	5
Hg	200	Bi (209)	Off	0.6 - 10	0.99978	0.00226	2.0	10	3	5
Pb	208	Bi (209)	Off	0.6 - 10	0.99986	3.14E-4	2.0	10	3	5

a. Spiked concentration is the spiked concentrations of different elements in the final measurement solutions after dilution.

Table 4. ICP-MS Results

		⁷⁵ As	¹¹¹ Cd	²⁰⁰ Hg	²⁰⁸ Pb		⁷⁵ As	¹¹¹ Cd	²⁰⁰ Hg	²⁰⁸ Pb	
Blank	Mean value	n.d.	n.d.	n.d.	n.d.	Lotion	Mean value	n.d.	n.d.	n.d.	
	RSD (n = 3)	---	---	---	---		RSD (n = 3)	---	---	---	---
Fortified blank	Mean value	5.23	5.23	5.1	5.08	Fortified lotion	Mean value	4.96	5.09	5.16	5.03
	RSD (n = 3)	2.03	0.76	3.57	1.59		RSD (n = 3)	11.08	1.21	3.12	2.62
Recovery (%)		105	105	102	102	Recovery (%)		99	102	103	101
Hard candy	Mean value	n.d.	n.d.	n.d.	0.0116	Froggy	Mean value	n.d.	n.d.	n.d.	0.0407
	RSD (n = 3)	---	---	---	1.79		RSD (n = 3)	---	---	---	0.66
Fortified hard candy	Mean value	5.12	4.92	5.05	5.13	Fortified froggy	Mean value	4.76	4.9	4.83	5.07
	RSD (n = 3)	3.51	1.68	0.19	1.74		RSD (n = 3)	0.95	0.74	1.46	1.92
Recovery (%)		102	98	101	102	Recovery (%)		95	98	97	101
Granola bar	Mean value	n.d.	0.125	n.d.	0.0186	Hemp flower	Mean value	0.0231	n.d.	n.d.	0.163
	RSD (n = 3)	---	1.63	---	3.02		RSD (n = 3)	6.07	---	---	1.09
Fortified granola bar	Mean value	4.55	5.07	5.09	5.06	Fortified hemp flower	Mean value	4.81	5.05	5.11	5.15
	RSD (n = 3)	13.53	1.61	1.58	1.15		RSD (n = 3)	9.68	1.22	0.6	1.48
Recovery (%)		91	99	102	101	Recovery (%)		96	101	102	100
MCT oil	Mean value	n.d.	n.d.	n.d.	n.d.	Beef jerky	Mean value	n.d.	0.0482	n.d.	0.057
	RSD (n = 3)	---	---	---	---		RSD (n = 3)	---	0.77	---	1.65
Fortified MCT oil	Mean value	5.39	5.04	5.13	4.94	Fortified beef jerky	Mean value	4.81	5.07	5.06	5.22
	RSD (n = 3)	1.36	1.54	1.24	1.89		RSD (n = 3)	7.54	1.17	1.25	1.86
Recovery (%)		108	101	103	99	Recovery (%)		96	100	101	103
Ghee	Mean value	n.d.	n.d.	n.d.	n.d.	Peanut butter	Mean value	n.d.	0.111	n.d.	0.0511
	RSD (n = 3)	---	---	---	---		RSD (n = 3)	---	3.92	---	1.32
Fortified ghee	Mean value	5.46	5.05	5.15	5	Fortified peanut butter	Mean value	5.22	5.28	5.13	5.11
	RSD (n = 3)	3.8	2.01	1.19	1.73		RSD (n = 3)	5.32	1.08	1.45	2.93
Recovery (%)		109	101	103	100	Recovery (%)		104	103	103	101
Hemp oil	Mean value	n.d.	n.d.	n.d.	0.0723	Conc. CBD oil	Mean value	n.d.	n.d.	n.d.	0.0662
	RSD (n = 3)	---	---	---	1.6		RSD (n = 3)	---	---	---	0.78
Fortified hemp oil	Mean value	5.27	5.18	5.2	5.09	Fortified Conc. CBD oil	Mean value	5.41	5.17	4.93	5.05
	RSD (n = 3)	11.91	1.53	2.1	1.5		RSD (n = 3)	10.43	0.7	2.03	2.74
Recovery (%)		105	104	104	100	Recovery (%)		108	103	99	100

Conclusion

- The MARS 6 was able to completely digest a wide variety of cannabis samples in a single batch
- Up to 24 samples can be run at a time in a single batch
- The MARSXpress Plus vessel provides the combination of high throughput and robustness necessary for this work
- The Glass Expansion Niagara Rapid Rinse accessory significantly reduced our analysis time
- The Shimadzu ICP-MS 2030 provided superior analysis and showed excellent recoveries and RSD's