

# Determination of $\Delta^9$ -tetrahydrocannabinol and two of its metabolites in whole blood, plasma and urine by UHPLC-MS/MS using QuEChERS sample preparation

**ASMS 2014** ThP600

Sylvain DULAURENT<sup>1</sup>, Mikaël LEVI<sup>2</sup>, Jean-michel GAULIER<sup>1</sup>,  
Pierre MARQUET<sup>1,3</sup> and Stéphane MOREAU<sup>2</sup>

<sup>1</sup> CHU Limoges, Department of Pharmacology and Toxicology,  
Unit of clinical and forensic toxicology, Limoges, France ;

<sup>2</sup> Shimadzu France SAS, Le Lizard 2, Boulevard Salvador  
Allende, 77448 Marne la Vallée Cedex 2

<sup>3</sup> Univ Limoges, Limoges, France

# Determination of $\Delta^9$ -tetrahydrocannabinol and two of its metabolites in whole blood, plasma and urine by UHPLC-MS/MS using QuEChERS sample preparation

## Introduction

In France, as in other countries, cannabis is the most widely used illicit drug. In forensic as well as in clinical contexts,  $\Delta^9$ -tetrahydrocannabinol (THC), the main active compound of cannabis, and two of its metabolites [11-hydroxy- $\Delta^9$ -tetrahydrocannabinol (11-OH-THC) and 11-nor- $\Delta^9$ -tetrahydrocannabinol-9-carboxylic acid (THC-COOH)] are regularly investigated in biological fluids for example in Driving Under the Influence of Drug context (DUID) (figure 1).

Historically, the concentrations of these compounds were determined using a time-consuming extraction procedure

and GC-MS. The use of LC-MS/MS for this application is relatively recent, due to the low response of these compounds in LC-MS/MS while low limits of quantification need to be reached. Recently, on-line Solid-Phase-Extraction coupled with UHPLC-MS/MS was described, but in our hands it gave rise to significant carry-over after highly concentrated samples. We propose here a highly sensitive UHPLC-MS/MS method with straightforward QuEChERS sample preparation (acronym for Quick, Easy, Cheap, Effective, Rugged and Safe).




Figure 1: Structures of THC and two of its metabolites

## Methods and Materials

Isotopically labeled internal standards (one for each target compound in order to improve method precision and accuracy) at 10 ng/mL in acetonitrile, were added to 100  $\mu$ L of sample (urine, whole blood or plasma) together with 50 mg of QuEChERS salts ( $\text{MgSO}_4$ /NaCl/Sodium

citrate dehydrate/Sodium citrate sesquihydrate) and 200  $\mu$ L of acetonitrile. Then the mixture was shaken and centrifuged for 10 min at 12,300 g. Finally, 15  $\mu$ L of the upper layer were injected in the UHPLC-MS-MS system. The whole acquisition method lasted 3.4 min.

# Determination of $\Delta^9$ -tetrahydrocannabinol and two of its metabolites in whole blood, plasma and urine by UHPLC-MS/MS using QuEChERS sample preparation

## UHPLC conditions (Nexera MP system)

Column	: Kinetex C18 50x2.1 mm 2.6 $\mu$ m (Phenomenex)
Mobile phase A	: 5mM ammonium acetate in water
B	: CH <sub>3</sub> CN
Flow rate	: 0.6 mL/min
Time program	: B conc. 20% (0-0.25 min) - 90% (1.75-2.40 min) - 20% (2.40-3.40 min)
Column temperature	: 50 °C

## MS conditions (LCMS-8040)

Ionization	: ESI, negative MRM mode		
Ion source temperatures	: Desolvation line: 300°C Heater Block: 500°C		
Gases	: Nebulization: 2.5 L/min Drying: 10 L/min		
MRM Transitions:			
	Compound	MRM	Dwell time (msec)
	THC	313.10>245.25 (Quan)	60
		313.10>191.20 (Qual)	60
		313.10>203.20 (Qual)	60
	THC-D <sub>3</sub>	316.10>248.30 (Quan)	5
		316.10>194.20 (Qual)	5
	11-OH-THC	329.20>311.30 (Quan)	45
		329.20>268.25 (Qual)	45
		329.20>173.20 (Qual)	45
	11-OH-THC-D <sub>3</sub>	332.30>314.40 (Quan)	5
		332.30>271.25 (Qual)	5
	THC-COOH	343.20>245.30 (Quan)	50
		343.20>325.15 (Qual)	50
		343.20>191.15 (Qual)	50
		343.20>299.20 (Qual)	50
	THC-COOH-D <sub>3</sub>	346.20>302.25 (Quan)	5
		346.20>248.30 (Qual)	5
Pause time	: 3 msec		
Loop time	: 0.4 sec (minimum 20 points per peak for each MRM transition)		

# Determination of $\Delta^9$ -tetrahydrocannabinol and two of its metabolites in whole blood, plasma and urine by UHPLC-MS/MS using QuEChERS sample preparation

## Results

### Chromatographic conditions

A typical chromatogram of the 6 compounds is presented in figure 1.




Figure 1: Chromatogram obtained after an injection of a 15  $\mu$ L whole blood extract spiked at 50  $\mu$ g/L

### Extraction conditions

As described by Anastassiades et al. J. AOAC Int 86 (2003) 412-31, the combination of acetonitrile and QuEChERS salts allowed the extraction/partitioning of compounds of interest from matrix. This extraction/partitioning process is not only

obtained with whole blood and plasma-serum where deproteinization occurred and allowed phase separation, but also with urine as presented in figure 2.




Figure 2: influence of QuEChERS salts on urine extraction/partitioning: A: acetonitrile with urine sample lead to one phase / B: acetonitrile, QuEChERS salts and urine lead to 2 phases.

# Determination of $\Delta^9$ -tetrahydrocannabinol and two of its metabolites in whole blood, plasma and urine by UHPLC-MS/MS using QuEChERS sample preparation

## Validation data

One challenge for the determination of cannabinoids in blood using LC-MS/MS is the low quantification limits that need to be reached. The French Society of Analytical Toxicology proposed 0.5  $\mu\text{g/L}$  for THC et 11-OH-THC and 2.0  $\mu\text{g/L}$  for THC-COOH. With the current application, the

lower limit of quantification was fixed at 0.5  $\mu\text{g/L}$  for the three compounds (3.75 pg on column). The corresponding extract ion chromatograms at this concentration are presented in figure 3.




Figure 3: Chromatogram obtained after an injection of a 15  $\mu\text{L}$  whole blood extract spiked at 0.5  $\mu\text{g/L}$  (lower limit of quantification).

The upper limit of quantification was set at 100  $\mu\text{g/L}$ . Calibration graphs of the cannabinoids-to-internal standard peak-area ratios of the quantification transition versus

expected cannabinoids concentration were constructed using a quadratic with 1/x weighting regression analysis (figure 4).




Figure 4: Calibration curves of the three cannabinoids

Contrary to what was already observed with on-line Solid-Phase-Extraction no carry-over effect was noted using the present method, even when blank samples were

injected after patient urine samples with concentrations exceeding 2000  $\mu\text{g/L}$  for THC-COOH.

## Determination of $\Delta^9$ -tetrahydrocannabinol and two of its metabolites in whole blood, plasma and urine by UHPLC-MS/MS using QuEChERS sample preparation

### Conclusions

- Quick sample preparation based on QuEChERS salts extraction/partitioning, almost as short as on-line Solid Phase Extraction.
- Low limit of quantification compatible with determination of DUID.
- No carry over effect noticed.

First Edition: June, 2014