SHIMADZU

Climbing the oligonucleotide ladder toward rapid and wide-ranging oligonucleotide analysis using MALDI-MS.

Dominique B. Figueroa¹, Vikki Johnson¹, M. Nazim Boutaghou², Jordan M. Frost¹ 1 Shimadzu Scientific Instruments, Carlsbad, CA, USA. 2 Shimadzu Scientific Instruments, Columbia, MD, USA.

Results obtained from gel electrophoresis workflows can be acquired in seconds harnessing MALDI-MS.

1. Introduction

Biopharmaceutical precision and continue to technologies medicine life science design, transform drug clinical diagnostics, and research. demanding mass spectrometry (MS) fast, high-throughput techniques for oligonucleotide analysis requiring mass sequence confirmation. Matrixdesorption/ionization laser assisted (MALDI-MS), provides a departure from toxic, time-consuming verification using and ethidium electrophoresis MALDI-MS techniques are bromide. amenable to rapid sample preparation mass verification in and high-throughput studies, able to provide results in seconds.



2-1. Sample Preparation

Oligonucleotide length standards were obtained as the 10/60 ladder from Integrated DNA Technologies (San Diego, CA). One µL volume was spotted onto Fleximass Stainless Steel MALDI Target Plates (Shimadzu Scientific Instruments, Columbia, MD). The matrix 0.5M 2',4',6'-trihydroxyacetophenone (THAP) was dissolved in 50% water, 50% acetonitrile. and Diammonium citrate was added to chelate sodium ions.

2-2. Instrumentation

Benchtop MALDI-8020 technology was used to detect the prepared oligonucleotide mixture (Figure 2).

3. Results

3-1. Oligonucleotide molecular ladder weights

Nucleic acid sequences for each oligonucleotide length standards were obtained. Average [M+H]+ values were



Figure An 1. oligonucleotide ladder used for oligonucleotide gel electrophoresis.




calculated using the Mongo Oligo Mass Calculator v2.08 from The RNA Institute, University at Albany, State University of New York. All oligonucleotides consist of single-stranded DNA and were unmodified with a hydroxyl group at both the 5' and 3' ends (Table 1)


Length	Sequence	Average [M+H]+
10mer	ATC GCG GAT T	3044.053
15mer	GCT GCG ACG AGG CTG	4635.066
20mer	ATC GCG GAT TAG CAC TAC GT	6118.052
25mer	ATC TCG GAT TAG CAC TAC GCA TCG G	7643.039
30mer	ATC GCG GAT TAG CAC TAC GCA TCG GTT ACA	9192.051
40mer	ATC GCG GAT TAG CAC TAC GCA TCG GTT ACA AAC GAG TAC C	12275.063
50mer	ATC GCG GAT TAG CAC TAC GCA TCG GTT ACA AAC GAG GAC CTG ATG CAC TT	15380.074
60mer	ATC GCG GAT TAG CAC TAC GCA TCG GTT ACA AAC GAG GAC CTG ATG CAC TTT GAC AGC ATG	18494.098

Table 1. Oligonucleotide sequences and expected m/z values to be detected.

MALDI-8020 Specifications			
Laser	355 nm, solid state, 200 Hz		
Mass Range	<i>m/z</i> = 1 – 500,000		
Resolution	>5000		
Sensitivity	250 fmol protein 250 amol peptide		
Ionization	Positive (+)		

Figure 2. Benchtop MALDI-8020 instrument and specifications.





MALDI-8020 Acqu	MA	
Mass Range	2000-20000	Baseline
Laser Power	150	Smoothi
Accumulation Rate	30 shots @ 200 Hz	Smoothi
Pulsed Extract (<i>m/z</i>)	12360	Peak Wi
Blanking Mass (<i>m/z</i>)	1500	Peak De

Figure 4. Data acquisition and processing parameters used to generate mass spectra above.

Research Use Only (RUO). Not for use in diagnostic procedures.