

### ASMS 2016 MP-466

Shuichi Kawana<sup>1</sup>, Yumi Unno<sup>1</sup>, Yukihiko Kudo<sup>1</sup>, Takero Sakai<sup>1</sup>, Takashi Kobayashi<sup>2</sup>, Shin Nishiumi<sup>2</sup>, Masaru Yoshida<sup>2</sup>, Noriyuki Ojima<sup>1</sup> <sup>1</sup> Shimadzu Corporation. 1, Nishinokyo-Kuwabaracho Nakagyo-ku, Kyoto 604–8511, Japan, <sup>2</sup> Kobe University Graduate School of Medicine. Research Build. A6-010, 7-5-1, Kusunoki-cho, Chu-o-ku, Kobe, Hyogo 650-0017, Japan

PO-CON1626E



### Introduction

Gas-Chromatography Tandem Mass Spectrometry (GC-MS/MS) is one of highly suitable techniques for metabolome analysis because of the high separation, reproducible retention times and sensitive selective mass detection. However, metabolites are commonly derivatized before GC-MS/MS analysis. In case of TMS derivatization, the samples should be analyzed within 24 hours after derivatization, because the TMS derivatives will deteriorate after 24 hours. Moreover, metabolome analysis requires the measurement of more than 100 samples. Therefore, an effective analysis procedure and a reduction of exposure of operators to toxic reagents are required. In order to improve analysis accuracy, efficiency and safety, an Automated TMS Derivatization GC-MS/MS System was developed. In this study, the system is evaluated for analysis of metabolites in human plasma.

### Methods and Materials

#### Sample Preparation

Reference; Nishiumi S et. al., Metabolomics, 2010 Nov;6(4):518-528



#### 

#### Analysis of Plasma Metabolites Using Gas-Chromatography Tandem Mass Spectrometry System with Automated TMS Derivatization

Analytical Conditions				
GC				
Inj. Temp. Column oven Temp. Linear velocity Split Ratio Injection volume	: 250°C : 60°C (2.00 min) → (15°C/min) → 330°C (3.00 min) : Total 23 min : 39.0 cm/sec : 30 : 1 µL			
Column				
BPX-5 (30m x 0.25mm l	.D. df=0.25µm, SGE)			
MS				
Interface Temp. Ion source Temp. Data acquisition	: 280°C : 200°C : MRM (475 compounds 950 transitions Avg. dwell time : 3.5ms Min. dwell time : 1.0ms			
Database				
Smart Metaboites Datab	base (Shimadzu)			

### Automated TMS Derivatization GC-MS/MS system



Figure 1 Automated TMS Derivatization GC-MS/MS system

#### Analysis Workflow

🕀 SHIMADZU

1) Conventional Manual Derivatization Method



2) Automated TMS Derivatization GC-MS/MS System



The developed system can derivatize plural samples simultaneously and analyze the samples immediately after TMS derivatization.

## Result

# Analysis of Plasma metabolites using Automated TMS derivatization GC-MS/MS System

The 179 metabolites were detected in pooled human plasma by the Automated TMS Derivatization GC-MS/MS System. To validate the system, 7 replicate samples were analyzed and the relative standard deviation (RSD) value in the 179 metabolites was calculated. For 2-Isopropylmalic acid as an internal standard, RSD% of peak area was 8.47%.



Figure 2 MRM Chromatograms of metabolites in pooled human plasma

For the detected metabolites, the peak area of each ion was calculated and normalized to the peak area of 2-isopropylmalic acid as an internal standard. The RSD% of 134 metabolites were less than 20%. The result was almost the same as the conventional manual derivatization method.

%RSD (n=7)	Number of Compounds		
- 5.00%	43		
5.01 - 10.00%	42		
10.01 - 15.00%	27		
15.01 - 20.00%	22		
≤ 20.00%	134		
> 20.00%	45		
Total	179		

	Table 1	Reproducibility	of relative	peak area	in pooled	human plasma
--	---------	-----------------	-------------	-----------	-----------	--------------

### Signal intensities with lapse of time

We have measured the samples which passed from 0 to 24 hours after TMS derivatization. The signal intensities of some metabolites like Lysine, Tyrosine, Kynurenine and Tryptophan were decreased gradually during 24 hours. As time passed after TMS derivatization, the TMS derivatives ratio in some compounds were changed by the reaction

progress. On the other hand, Automated TMS Derivatization system can analyze samples immediately after TMS derivatization. The samples which passed from 0 to 24 hours after the methoximation were analyzed. As the result, the signal intensities in these compounds were not changed even with the lapse of time.



Figure 3 Mass Chromatograms of Lysine-4TMS, Tyrosine-3TMS, Kyunurenine-3TMS and Tryptophan-3TMS. Blue lines are 0 hour and red lines are 24 hours after setting the samples on the auto-sampler.



### Conclusions

The Automated TMS Derivatization GC-MS/MS System was developed and evaluated for analysis of metabolites in human plasma. The novel system improves the accuracy, efficiency and safety for metabolite analysis in human plasma compared to the conventional manual derivatization method. This system contributes to the dissemination of metabolomics studies.

A part of this study was financially supported by Medical Research and Development Programs Focused on Technology Transfer: Development of Advanced Measurement and Analysis Systems of Japan Agency for Medial Research and Development.

Disclaimer: GCMS-TQ8040 is intended for Research Use Only (RUO). Not for use in diagnostic procedures.



Shimadzu Corporation

www.shimadzu.com/an/

#### For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "@". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names of ther the networks and trade names of the na

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

First Edition: May, 2016