## 

# **Chromatographic and Mass Spectrometry Approaches to Supporting Biologics Development**

<u>Mridul Mandal</u>; Brian Domanski<sup>1</sup>; Stephen Kurzyniec<sup>1</sup>; Stephen Swartz<sup>1</sup> <sup>1</sup>Shimadzu Scientific Instruments, Columbia, MD

### **Overview**


Accurate characterization of biologics is vital to the development of biotherapeutics.


- It is essential to know the properties of biotherapeutic products to optimize bioprocess production, product formulation, and dose.
- The latest innovations from Shimadzu include LC, LC/MS, mini and conventional MALDI, and in-house as well as collaborator informatics to support the complex analysis of biologics.

### **1. Introduction**

Biopharmaceuticals are being developed using biotechnology techniques such as genome editing and cell fusion. These drugs are getting attention worldwide as they are expected to be effective in treating various chronic diseases including cancers. Antibody drugs, such as monoclonal antibodies (mAbs) and antibody-drug conjugates (ADCs), are known to offer high therapeutic efficacy and reduced side effects due to their specificity and affinity for target molecules. However, because these biopharmaceuticals are manufactured using animal cells, the challenge is to ensure their structural homogeneity is not encountered in small molecule pharmaceuticals manufactured by chemical synthesis. Accordingly, biopharmaceuticals require appropriate quality controls at every production step. For example, ICH-Q6B (1), (2), which proposes specifications and test procedures for biological products, stipulates that productrelated impurities should be separated and/or their percentage levels determined in the manufacturing process. In most cases, these analyses are performed using liquid chromatography

### 2. Results









Figure 9: N-linked glycans are attached to proteins at Asn residues on the Fc region.






### 2. Methods

A variety of instruments have been used for analysis, including UHPLC, LCMS triple quad, LCMS QTOF, and MALDI. A brief summary of the methods is given below:


#### LC Method


Adsorption of an analyte by a wetted surface of a UHPLC instrument poses some critical challenges when analyzing biomolecules. In part of the study, the Nexera XS inert system was used to separate biomolecules by combining the elevated pressure tolerance of a UHPLC system with a flow path that was completely inert, as there were no wetted metal surfaces in the sample flow path, assuring ultra-high corrosion resistance.



**Figure 3:** Mixed standard solutions of AMP, ADP, and ATP (50  $\mu$ g/mL) measured in the "metal-based".

#### **Drug-to-Antibody Ratio by Nexera XS inert**





**Figure 4:** Mixed standard solutions of AMP, ADP, and ATP (50  $\mu$ g/mL) measured in the "metal-free".




Figure 7: Deconvoluted mass spectrum of NIST mAb.




Figure 10: Fluorescence chromatogram (upper) and extracted ion chromatogram (lower) for procainamide labeled RNaseB Glycans.




Figure 11: Detection of trypsin fragments containing CDR1-3 using standard monoclonal antibodies subjected to trypsin digestion (Data Analysis by Skyline).

#### LCMS-9030 QTOF Method

Innovative TOF architecture allows ultra-fast, precise pulsing of ions into the flight tube and back to the detector. Thus, highspeed data acquisition is compatible with high-throughput laboratories. Utilizing the Nexera Mikros, the Micro-ESI 9030 interface allows microflow analysis to be performed on the Q-TOF LCMS-9030. Using the system, recombinant human IgG1\* NIST mAbs, Trastuzumab, and Bevacizumab were evaluated as model biotherapeutics.

#### LCMS-8060 (Triple Quad) Method


LCMS-8060NX redefines LC-MS/MS detection limits and quantitative capabilities. A new ion focusing technique in the 8060NX transforms trace quantitative detection by redefining "sensitivity". In addition, the fastest polarity switching available in the industry enables switching from positive to negative modes in 5 msec. With data acquisition occurring so quickly, it is possible to perform triggered MS/MS scans and MRM quantification using the same acquisition method.

#### MALDI Method

Shimadzu's MALDI-8030 is part of a long line of MALDI-TOF instrumentation. This dual-polarity, benchtop MALDI-TOF mass spectrometer provides outstanding performance in a compact footprint, making it the ideal choice for today's high-quality laboratories. It can deliver the performance required for analytical method development and quality control applications.

Figure 6: NIST mAb Intact TIC chromatogram (top, left), deconvoluted spectrum (top right), and

Figure 13: An unparalleled convenient and rapid workflow provided by the nSMOL Antibody BA Kit dramatically improves the productivity and robustness of LCMS mAb bioanalysis.






Figure 12: Calibration Data for the Transition terfere with analysis 531.27 > 603.32 on a Linear Scaling .

### 4. Conclusions

The comprise innovator topics biopharmaceuticals, biosimilars, antibodydrug conjugates, and keeping up with changes in the regulatory environment for development of biologics.

In addition to full protein, subunit, peptide fragment, and N-Glycan analysis, these workflows can analyze mAbs, proven (ADC), antibody drug conjugates oligonucleotides, proteins, and peptides.

A newly developed bioinert UHPLC system

#### • Data processing

mass spectrum (bottom).

Chromatographic data were collected and analyzed by Accurate characterization of monoclonal antibodies is essential to development of LabSolutions software. LC-MS data were taken by LabSolutions biotherapeutics. Thorough understanding of biotherapeutic properties aids in the optimization and analyzed by Protein Metrics and Skyline. The MALDI-MS data of bioprocess production, product formulation, and product dosage. LCMS-9030 QTOF and were collected by MALDI Solutions and analyzed with Protein Nexera XS inert have been used to characterize the recombinant human NIST mAb and Bevacizumab Metrics software

Fc S Subject to the server of	effectively eliminates metal ion interaction with targets such as proteins and nucleic
<u>MRM peptide</u> Identified in both digestion Identified only in-solution digestion	acids, preventing peak tailing and low
Figure 14: nSMOL coupled with LC-MS bioanalysis	recovery of target components. As a result,
for monitoring the Fc-fusion biopharmaceuticals	data can be obtained with high sensitivity
Etanercept and Abatacept in human serum.	and sharp peaks even with low
	concentrations of sample.

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.