


Trace level determination of AsH₃, PH₃, COS and H₂S in Propylene by using GC-MS

Prashant Dattatray Hase¹; Sanket Anand Chiplunkar¹; Dheeraj Handique¹; Nitish Suryawanshi¹; Aseem Wagle¹; Rahul Dwivedi¹; Jitendra Kelkar¹; Pratap Rasam¹;Satyendra Singh²; ¹Shimadzu Analytical (India) Pvt.Ltd., 1 A/B Rushabh Chambers, Makwana Road, Marol, Andheri (E), Mumbai-400059, Maharashtra, India.²Shimadzu Analytical (India) Pvt Ltd., Plot No-2, Sethi Chamber, M O R Land, New Rajender Nagar, New Delhi, Delhi 110060.

. Introduction

Arsine (AsH₃), Phosphine (PH₃), Carbonyl Sulfide (COS), and Hydrogen Sulfide (H₂S) are compounds present as impurities (Figure 1) in petroleum products and their derivatives. They contribute to the corrosion of facilities and the inhibition of catalytic reactions, leading to the deterioration of product quality. High purity propylene is used as feedstock to produce polypropylene and the quality of the monomer is critical for successful polymerization. The presence of trace amounts of AsH₃, PH₃, COS, and H₂S can have detrimental effects on the catalyst and product yield. To maximize the lifetime of the catalyst, propylene producers need to monitor AsH₃ PH₃ COS and H₂S at the lowest possible concentration levels.

Shimadzu GCMS-QP2020 NX (Figure 2) single quadrupole system with vaporizer and Gas sampling valve (GSV) was used to simultaneously determine AsH₃, PH₃, COS and H₂S impurities in propylene at trace levels.



Arsine (AsH₃)

Phosphine (PH₃) Carbonyl Sulphide (COS)

Hydrogen Sulphide

Figure 1. Structure of AsH₃, PH₃, COS and H₂S

2. Materials and methods

2-1. Procedure

Individual standard of about 500 ppb for AsH₃, PH₃, COS and H₂S in balance Helium was procured for the study. The standard was further diluted to prepare the linearity concentration levels as mentioned in Table 1.

Table 1. Linearity standards preparation

Sr. no	Linearity Conc. (Mol-ppb)	AsH ₃	PH ₃	cos	H ₂ S
1	5	√	✓	✓	NA
2	10	√	√	✓	NA
3	25	√	√	✓	NA
4	50	√	√	√	\checkmark
5	100	✓	✓	✓	√
6	250	NA	NA	NA	√
7	500	NA	NA	NA	√

✓ Applicable. NA – Not Applicable

The liquified petroleum sample was passed through the vaporizer to eliminate the sample condensation effects. Then the sample was introduced into Shimadzu GCMS-QP2020 NX single quadrupole system automatically through GSV equipped with 1mL loop. Sulfinert® treated parts were used to ensure inert sample path for trace level analysis.

Plot column was used for simultaneous determination of these impurities in propylene matrix. The concentrations of targeted contaminants are determined against multi-point external standard calibration curve. The standard and samples were analyzed under identical conditions and injection volumes.



Figure 2. Shimadzu GCMS-QP2020 NX

2-2. Analytical Conditions

The analysis was carried out on Shimadzu GCMS-QP2020 NX with GSV and liquid sample vaporizer as per the conditions given in Table 2.

Table 2. Analytical conditions

Chromatographic parame	<u>eters</u>				
• Column	: Plot (50 m x 0.32 mm l.D. x 10.0 μ m)				
 Injection Temp 	: 80 °C				
 Injection Mode 	: Split				
 Split Ratio 	: 1:5				
 Carrier Gas 	: Helium				
 Flow Control Mode 	: Column flow				
 Column Flow 	: 1.5mL/min				
 Total Program Time 	: 35.00 min				
 Column Temp. Program 	: Rate (°C /min)	Temperature (°C)	Hold time (min)		
		40.0	10.00		
	1.50	50.0	3.00		
	10.00	150.0	5.33		
 GSV temperature 	: 80 °C				
 Vaporizer temperature 	: 100 °C				
Mass Spectrometry parar					
 Ion Source Temp. 	: 150.0 °C				
 Interface Temp. 	: 150.0 °C				

: El (Electron Ionization)

: 34

3. Results

Ionization Mode

Acquisition Mode

AsH₃ m/z

COS m/z

• PH₃ m/z

H₂S m/z

The analysis of of AsH₃, PH₃, COS, and H₂S was carried out as per the method shown in Table 2 and Selective Ion Monitoring (SIM) chromatogram of 50 ppb standard mixture acquired using this method is shown in Figure 3




Figure 3. SIM chromatogram for 50 mol-ppb standard

The Relative Standard deviation (RSD) for three replicate injections of AsH₃, PH₃, COS at 5 mol-ppb and H₂S at 50 mol-ppb was found to be less than 5%. Linearity was plotted at concentration level ranging from 5 to 100 mol-ppb for AsH₃, PH₃, COS and from 50 to 500 mol-ppb for H₂S. For all target components, the regression coefficient $(r^2) \ge 0.995$ was obtained (refer Figure 3 to 6). The limit of detection (LOD) and limit of quantitation (LOQ) were determined based on S/N ratio. Statistical evaluation of the data showed that the method was selective, sensitive and reproducible (refer Table 3).




Figure 3. Calibration curve, linearity overlay graph and precision overlay at 50mol-ppb for H₂S




Figure 4. Calibration curve, linearity overlay graph and precision overlay at 5mol-ppb for PH3




Figure 5. Calibration curve, linearity overlay graph and precision overlay at 5mol-ppb for AsH₃




Figure 6. Calibration curve, linearity overlay graph and precision overlay at 5mol-ppb for COS

Table 3. Analytical Results

		Linearity		Precision		*LOQ	*LOD
Name	m/z	range (mol-ppb)	r ²	Conc. (mol-ppb)	%RSD (n=3)	(**S/N=10) (mol-ppb)	(**S/N=3) (mol-ppb)
PH ₃	34	5 – 100	0.9994	5	4.6	3.6ppb	1.2ppb
H ₂ S	34	50 – 500	0.9967	50	4.9	18ppb	6ppb
ASH ₃	76	5 – 100	0.9985	5	3.1	1.2ppb	4.1 ppb
COS	60	5 – 100	0.9983	5	4.0	0.1ppb	0.4 ppb

^{*}LOD/LOQ are determined based on the values obtained from GCMS Solutions software.

4. Conclusion

The GCMS-QP2020 NX is an advanced gas chromatograph mass spectrometer which provides highest sensitivity in its class through its unique design and superior noise reduction technology, thereby enabling ultra trace level determination of low molecular weight compounds like AsH₃, PH₃, COS and H₂S in propylene sample.

^{**} Noise calculation method used is ASTM with interval of 0.5mins.

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic