

Characterization of Cyclopropenoic Fatty Acid Containing Triacylglycerols of Monguba Oil by MALDI-TOF and High-Energy CID

MALDI-TOF/TOF Mass Spectrometry

Simona Salivo¹; Alexandra Pereira Rodrigues²; Glaucia Maria Pastore²; Matthew E. Openshaw¹; Gerald Stübiger³

¹Shimadzu, Manchester, UK; ²University of Campinas, Brazil; ³Medical University of Vienna, Austria

Overview

We present the analysis of lipids from Monguba (*Pachira aquatica Aubl.*) seeds by MALDI-(TOF)/TOF MS. To the best of our knowledge, this is the first time MALDI-TOF/TOF has been used for the structural analysis and confirmation of triacylglycerols containing cyclopropenoic fatty acid e.g., steric acid and derivatives.

1. Introduction

Monguba (*Pachira aquatica Aubl.*) is a tree native to tropical regions, which is traditionally used in folk medicine. Its seeds are rich in oil containing saturated, monounsaturated and polyunsaturated fatty acids ranging from 46.67–84.87%, 6.62–39.30% and 5.24–11.81%, respectively. Beside these common fatty acids (FAs), the oil is characterized by varying amounts of cyclopropenoic fatty acids (CPFAs), i.e., steric acid, dihydrosteric acid, 2-hydroxy-steric acid (Figure 1), which are of great interest because of their anti-cancer potential. On the other hand, due to the toxicity for humans, their contents are of importance for potential food applications of the oil. Our work is intended to provide a method for the rapid detection and structural analysis of triacylglycerols (TAGs) of Monguba oil using a MALDI-TOF-(MS)/MS approach, based on a previous gas chromatography study.




Figure 2. Analysis workflow of native Monguba (*Pachira aquatica Aubl.*) oil.

3. Results

3-1. MALDI-MS analyses

The TAG profiles of Monguba oil were obtained on the MALDI-8030 in positive ion mode. Figure 3 shows the TAG profiles of the un-heated (red trace) and heated (blue trace) oils. TAG identification was carried out by HE-CID MALDI-MS/MS analyses (section 3-2).




Figure 3. MALDI-MS spectrum of Monguba oil: unheated (MO1; red trace) and heated (MO2; blue trace). TAG identification was carried out by high-energy CID MALDI-MS/MS. 16:0 = palmitic; 18:0 = stearic; 18:1 = oleic; 18:2 = linoleic; 19:1cy = steric; 19:0cy = dihydrosteric; 19:1cy-2OH = 2-hydroxy-steric.

3-2. High-energy CID MALDI-MS/MS analyses

High-energy-CID MALDI-MS/MS analyses were conducted on all visible, high- and low-abundant TAG peaks, with particular focus on those species containing CPFAs. MS/MS spectra were carefully interpreted for the presence of fragments indicative of: 1) the *sn*-positions of the fatty acids in the TAGs; 2) monoacylglycerols (MGs) and diacylglycerols (DAGs); 3) charge-remote fragments (CRFs) for the presence and position of the cyclopropene ring and double bond(s). Figure 4 shows the HE-CID MS/MS spectrum of *m/z* 883, which resulted in the combination of TG(16:0/19:1cy-2OH/16:0) and TG(16:0/18:1/18:0). Figure 4A shows the DAG fragment region, where all FAs composing the two TAGs were confirmed. Figure 4B shows the region of the MS/MS spectrum of the free FAs ($[M + Na]^+$), MGs ($[M] + H^+$), and diagnostic ions for the *sn*-1/3 and *sn*-2 positions. The fragments at *m/z* 359 and 331 indicate that the 19:1cy-2OH and 18:1 FAs are in the *sn*-2 position, respectively. Figure 4C shows the CRF region for a full, detailed structural characterisation of the FA chains. For the 2-hydroxy-steric acid (19:1cy-2OH), the diagnostic CRF at *m/z* 717 was successfully detected (blue series) along with the 66 Da unit, indicating the cyclopropene ring is present. For the oleic acid (18:1), the diagnostic CRF at *m/z* 729 was successfully detected (pink series), and the 40 Da unit, characteristic of a double bond, is present. CRFs of palmitic (16:0) and stearic (18:0) acids can also be observed (green and cyan series, respectively). Figure 5 shows the CRFs region of the HE-CID MS/MS spectrum of the TAG at *m/z* 867, which was identified as TG(16:0/19:1cy/16:0). Once again, the CRFs diagnostic of the steric acid were found (*m/z* 659, 687, 701; purple series).




Figure 4. High-energy CID MALDI-MS/MS spectrum of *m/z* 883 [TG(16:0/19:1cy-2OH/16:0) + TG(16:0/18:1/18:0)]. A) Diacylglycerol region. B) Free fatty acids, monoacylglycerols [], *sn*-position fragments region. C) Charge remote fragment region, confirming the presence and position of the cyclopropene ring and double bond (19:1cy-2OH and 18:1, respectively).

Conclusions

This work demonstrates the powerful combination of a benchtop linear MALDI-TOF instrument and a high-performance MALDI-TOF/TOF mass spectrometer for the rapid profiling and detailed characterisation of the TAGs in native Monguba oil. The high-sensitivity, high-resolution MS/MS capability of the MALDI-7090 allowed: 1) the FA composition of the TAGs to be obtained; 2) the presence of CPFAs, such as steric, dihydrosteric and 2-hydroxy-steric acid to be confirmed; 3) positional information of the FAs in the TAG to be obtained, as well as a full characterisation of the FA chains including cyclopropene ring and double bond positions based on the detection of CRFs. The differences in the TAG profiles of the un-heated and heated oils will be the subject of a future investigation.