# SHIMADZU

# Peptide Mapping of a Monoclonal Antibody Using an Integrated Protein Digestion LCMS Platform (Perfinity-QTOF LCMS-9030)

Xiaomeng (Kate) Xia<sup>1</sup>, Evelyn Wang<sup>1</sup>, Stephen Kurzyniec<sup>1</sup>, Tairo Ogura<sup>1</sup>, Yoshiyuki Okamura<sup>1</sup>, M. Nazim Boutaghou<sup>1</sup> 1 Shimadzu Scientific Instruments, Columbia, MD 21046 U.S.A.

### . Overview

The combination of the Perfinity, an integrated online protein digestion platform, and a Q-TOF provides an efficient and reproducible protein digestion protocol for peptide mapping analysis.

### 2. Introduction

Typical peptide-mapping procedures, especially for protein digestion, consist of many laborious steps, which can vary due to the differences among protocols or operators. As such they are susceptible to changes that affect reproducibility, reduce assay sensitivity, and significantly increase analysis time. An automated, new method is developed to overcome these drawbacks and provide a simple, efficient, and reliable alternative. We utilized Perfinity workstation – an automated on-column protein digestion platform – directly coupled with a LCMS-QTOF mass spectrometer for high throughput and high sequence coverage peptide mapping.



### 3. Methods

Peptide mapping was performed using the NISTmAb standard (reference material 8671). Three types of enzymatic columns – trypsin column, nonreduction/alkylation (NoRA) trypsin column, and Lys-C column – were used for protein digestion on Perfinity. A C18 desalt column was used to clean up the samples after digestion and a Restek Ultra C18 (2.7µm, 150mm x 2.1mm) column was used for reserved phase separation.

All data was obtained on a Q-TOF mass spectrometer, LCMS-9030. Data dependent acquisition (DDA) with three dependent events was used for data acquisition. For one set of experiments, NISTmAb was diluted to 1 mg/mL with water then directly injected into Perfinity for a 6-minute on-column protein digestion. For another set of experiments, NISTmAb was first reduced by adding 50 mM dithiothreitol (DTT), then alkylated with 500 mM iodoacetamide prior to the injection into Perfinity.

Protein Metrics software was used for comprehensive peptide identification of the digested protein.



#### NISTabb Light Chain

4. Results



NISTmAb peptide coverage. Black: trypsin column without Figure 1: reduction/alkylation (RA). Blue: NoRA trypsin column without RA. Red: trypsin column with RA. Green: NoRA trypsin column with RA. Magenta: Lys-C column with RA.

	Sequence coverage (%)						
NISTmAb	Trypsin column w/o RA, 50 °C	NoRA trypsin column w/o RA, 70 °C	Trypsin column with RA, 50 °C	NoRA trypsin column with RA, 70 °C	Lys-C column with RA, 50 °C		
Heavy chain	30.21	37.33	98.67	100	94.42		
Light chain	44.13	50.23	100	100	95.31		

attributes post translational modifications (PTMs), such as oxidation, glycosylation, and deamidation, was identified and quantified by relative abundance. Figure 2 illustrates the MS/MS spectra comparison of the native (precursor at m/z = 418.2207, +2) and Met-oxidized peptides (precursor at m/z =426.2182, +2) in the heavy chain 252-258. The differences in m/z (+15.99 Da in the oxidized form) for the y4–y5 fragment ions (green boxes) clearly distinguish the modified forms from the native species, and the m/z of y1-y3 ions remain constant, indicating that the Met-4 in the peptide is the location of oxidation. Table 2 summarizes the quantification results of the critical PTMs identified with different digestion conditions.

- samples.



Figure 2: PTM (methionine oxidation) analysis. MS/MS spectra of Met-oxidized and native peptides (heavy chain peptide 252-258). Top: oxidation at Met 255 in heavy chain. Bottom: native peptide. Green boxes show confirmed fragment ions. Disclaimer: All content contained herein resulted solely from Shimadzu, and no conflict of interest exists. The products and applications are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

## **MP 608**

**Critical PTMs Identification and Quantification.** A set of quality

 Complete conversion of N-terminal glutamine (Gln) to pyroglutamate (Pyro-Glu) was observed under all conditions

• The loss of lysine from the C-terminus of the heavy chain, in the range of 87-90%, was comparable at different digestion conditions.

• Methionine oxidation at heavy chain M255, and asparagine deamidation at heavy chain N387 and N392 are usually monitored for NISTmAb. Similar levels for the two modifications were detected for all digested

• High abundance of glycosylation at the heavy chain N300 was observed for the reduced and alkylated samples where the main glycans were complex biantennary oligosaccharides containing from 0 or 1 nonreducing galactose with fucose attached to the reducing end of Nacetylglucosamine (A2G0F, A2G1F).

- The sequence was not covered for the non-reduced/alkylated samples, so the glycosylation modification was not detected.
- Overall, the reported PTM values are very close among experiments and show low amounts of sample preparation-induced modifications.

### Modification

HC Q1+GIn -> PyroGlu HC M255 + Oxid HC ~N392/N387 + Deam HC K450 Lys loss HC N300+A2G0F

HC N300+A2G1F

HC N300 unglycosylated

### **5.** Conclusions

- using Shimadzu's unique solution: Perfinity-LCMS-9030 system. separation, and accurate detection upgraded the process of peptide mapping from time-consuming and tedious to an easy, reproducible
- We demonstrate here a fast, reproducible peptide mapping workflow • The combination of automated sample preparation, reversed-phase and automated.
- The analysis of NISTmAb RM 8671 protein produced outstanding sequence coverage with high confidence due to the QTOF mass accuracy. Reduction and alkylation steps were necessary in order to improve sequence coverage. Even with the added steps, the total analysis time was still less than two hours.
- The most targeted modifications were successfully identified and quantified relatively. We confirmed that PTM induced by sample preparation were low and insignificant.



	Relative abundance (%)							
	Trypsin column w/o RA, 50 °C	NoRA trypsin column w/o RA, 70 °C	Trypsin column with RA, 50 °C	NoRA trypsin column with RA, 70 °C	Lys-C column with RA, 50 °C			
	100	100	100	100	100			
	2.39	3.01	3.63	4.97	2.53			
	1.24	2.89	1.21	1.95	2.54			
S	90.55	89.75	87.73	88.79	90.23			
F	NA	NA	40.14	38.36	37.52			
F	NA	NA	36.65	34.11	33.56			
	NA	NA	1.53	1.28	1.28			

**Table 2:** Summary of critical PTMs identified and quantified for
 NISTmAb with different digestion conditions