SHIMADZU

Characterization of Unsaturated Fatty Acids in Negative OAD-MS/MS using LCMS-9050

Hidenori Takahashi, Mami Okamoto, Yuta Miyazaki, Yohei Arao, Natsuyo Asano Shimadzu Corporation, Kyoto, Japan

1. Introduction to OAD-MS/MS

- ◆ While low-energy CID-MS/MS is one of the most effective fragmentation techniques for structural analysis, it may not be ideally suited for the analysis of certain isomers.
- Several novel fragmentation techniques have been proposed to complement low-energy CID-MS/MS.

Table 1. Example of proposed novel fragmentation techniques

Electron-based fragmentation


EIEIO, ECD(Electron Capture Dissociation) by Zubarev et al. (1996)

Anion-based fragmentation

ETD(Electron Transfer Dissociation) by Syka et al. (2004)

Photon-based fragmentation

IRMPD (Infrared), UVPD (Ultraviolet), BRID (Blank body infrared)



2. Lipid MS/MS: CID vs. OAD vs. HAD Comparison

◆ CID-MS/MS does not provide detailed structural information within carbon chains. Instead, CID selectively cleaves labile polar head groups (*m*/*z* 184).




Fig. 2-1. Typical CID spectrum of a model lipid of PC (18:1).

◆ OAD-MS/MS clearly provides C=C positional information. Atomic oxygen selectively oxidizes and cleaves at C=C.




Fig. 2-2. Typical OAD spectrum of a model lipid of PC (18:1).

• HAD-MS/MS provides sequential structural information within carbon chains. The sensitivity of HAD is lower than OAD for C=C position assignment, since fragment ions become too complicated in a complex sample mixture.





Fig. 2-3. Typical HAD spectrum of a model lipid of PC (18:1).

3. Negative OAD-MS/MS in lipid analysis

Table 2. Analytical Condition


HPLC condition System: Shimadzu Nexera X3 Flow rate: 0.3 mL/min Mobile phase A: 1:1:3 (v/v/v) ACN:MeOH: water with ammonium acetate (5mM) and 10nM EDTA Mobile phase B: 100% IPA with ammonium acetate (5mM) and 10nM EDTA. **Column:** Acquity UPLC Peptide BEH C18 (50 × 2.1mm; 1.7µm; Waters, USA) Gradient: Omin 0% (B); 1min 0% (B); 5min 40% (B); 7.5min 64% (B); 12min 64% (B); 12.5min 82.5% (B); 19 min 85% (B); 20min 95% (B); 20.1min 0% (B); and 25min 0% (B).


Sample: Lipid extract from a mouse's brain with the addition of several commercially available fatty acid standards. The method for preparing the lipid extract adheres to the procedure described in Uchino.H et al. Commun Chem. 5, 162 (2022).



- ♦ OAD was successfully integrated into the
- ◆ OAD has demonstrated its effectiveness as a
- Without any need for derivatization, OAD

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.





MS condition System: Shimadzu LCMS-9050 Mode: OAD Polarity: Negative ESI CE: +10 V MS/MS: Auto-MS/MS Top 10