BHIMADZU

Analyses of biopharmaceuticals using high-resolution multi-turn TOF-MS system

Yusuke Tateishi¹, Hiroyuki Miura¹, Hiroko Morinaga¹, Koichi Kimura¹, Tetsuo Iida¹, Junna Nakazono¹, Masaru Nishiguchi¹, Osamu Furuhashi¹, Daisuke Okumura¹, Yuki Yamaguchi², Susumu Uchiyama² 1 Shimadzu Corporation, Kyoto, Japan, 2 Osaka University, Osaka, Japan

1. Introduction

In recent years, the development of new biopharmaceuticals such as next-generation antibodies and viral vectors has been promoted [1]. In the field of quality control of these biopharmaceuticals, analyses of structural changes such as chemical modification have important roles. For this purpose, we have been developing a highresolution multi-turn (MT) TOF-MS. In this presentation, we demonstrate the analyses of oxidized monoclonal antibodies (mAb) and the adeno-associated virus (AAV) capsid proteins.

2. MT TOF-MS system

Our MT TOF-MS consists of rotationally symmetric sector electrodes. Ions fly along a unique 3D open-loop orbit (Figure 1). It achieves nearly 50-m flight path within the compact size of Φ500 x 240 mm [2]. High resolution of over 150k has been demonstrated by simulation (Figure 2).




Figure 2 Simulated (a) ion trajectories and (b) TOF spectrum.

-40 -32 -24 -16 -8 0 8 16 24 32 40

Figure 3 shows the overview of the MT TOF-MS system. The system consists of high-performance liquid chromatography (HPLC) system, UV detector, electrospray ionization (ESI) source, quadrupole (QP) mass filter, collision cell (CC), linear ion trap (LIT), and MT TOF-MS. lons passed through QP and CC are temporary accumulated and compressed in LIT. Then these ions are extracted and injected toward the MT TOF-MS.





Figure 3 Overview of MT TOF-MS system.

A mass spectrum of myoglobin for a performance test is shown in Figure 4. We have achieved mass resolution of over 200k.



3. Analyses of biopharmaceuticals

Analyses of oxidized mAb subunits and the AAV viral proteins (VPs) were conducted.

3-1. Sample Preparations

Oxidized mAbs

For oxidation, 70% tBHP was added to NISTmAb (NIST8671, Sigma-Aldrich) and incubated at room temperature. Then, the whole amount of samples was applied to MicroSpin G-25 columns and the buffer was exchanged to PBS (pH7.4). FabRICATOR Enzyme and IgGZERO Enzyme (Genovis AB) were added and digested at 37°C for 60 min. PBS and 1M TCEP hydrochloride solution was added and incubated at 37°C for 60 min to reduce them to subunits. 6% TFA was added to stop the reaction.

AAV-VPs

10% acetic acid was added to AAV6-empty and incubated at room temperature for 15 min to fragment into VP1, VP2 and VP3 [3].

3-2. Analytical Conditions

LC-UV-MT TOF-MS system has been constructed, with which quantitative analysis by LC-UV and high-resolution mass spectrometry by MT TOF-MS were achieved.

HPLC conditions (Shimadzu HPLC)

	Oxidized mAbs
Column	Restek C4 5um, 150 x 2.1
Mobile Phase A	0.1%FA in H ₂ O
Mobile Phase B	0.1%FA in ACN
Flow Rate	200 uL/min
Injection volume	20 uL
UV	_

MS conditions (MT TOF-MS)

	Oxidized mAbs
Ionization	ESI
Neb gas	3 L/min
Heat gas	10 L/min
Dry gas	10 L/min
I.F. voltage	+4 kV
DL temp.	250 °C
I.F. temp.	300 °C
H.B. temp.	400 °C

AAV-VPs

AAV-VPs

1.5 L/min

5 L/min

10 L/min

+3.5 kV

200 °C

200 °C

300 °C

ESI

Waters ACQUITY BEH C4 1.7um, 150 x 2.1 0.1%DFA in H2O 0.1%DFA in ACN 200 uL/min 50 uL 280 nm










Changes of oxidation products as a function of oxidation time were confirmed with over 100k resolution (Figure 5).



LC chromatogram of VPs and mass spectrum of VP3 were obtained. Phosphorylated VP3 was observed in mass spectrum (Figure 6). 50000

Figure 6 (a) TIC of AAV6 and (b) Mass spectrum of VP3 and its phosphorylation product.




3-3. Results Analysis of oxidized mAb subunits

Fc/2 ~24100 Da(+17, m/z 1421)

Figure 5 Changes in oxidant formation of mAb subunits.

Analysis of AAV-VPs



4. Conclusions

 We demonstrated the high-resolution mass analyses of oxidized mAbs and the AAV capsid proteins with MT TOF-MS.

References

[1] M. Amano et al., Anal. Chem., 86, 7536-7543 (2014)

[2] Y. Tateishi et al., The 70th ASMS Conference, MP-019 (2022)

[3] H. Oyama et al., HUMAN GENE THERAPY, vol 32, number 21-22, 1403–1416 (2021)

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

The authors declare no competing financial interest.