

Rapid Characterization of Alkaloids using Probe ESI Q-TOF LCMS-9050 in OAD-MS/MS

Hidenori Takahashi, Kaoru Nakagawa, Mami Okamoto, Yuta Miyazaki, Yohei Arao, Tetsuo Iida Shimadzu Corporation, Kyoto, Japan

1. Introduction to OAD-MS/MS

- ♦ While low-energy CID-MS/MS is one of the most effective fragmentation techniques for structural analysis, it may not be ideally suited for the analysis of certain isomers.
- Several novel fragmentation techniques have been proposed to complement low-energy CID-MS/MS.

Table 1. Example of proposed novel fragmentation techniques

Electron-based fragmentation

EIEIO, ECD(Electron Capture Dissociation) by Zubarev et al. (1996)

Anion-based fragmentation

ETD(Electron Transfer Dissociation) by Syka et al. (2004)

Photon-based fragmentation

IRMPD (Infrared), UVPD (Ultraviolet), BRID (Blank body infrared)



2. Introduction to a direct ionization of PESI

◆ Probe Electro Spray Ionization (PESI) is one of the direct ionization techniques. Fig. 2 shows the scheme of the PESI system.



Fig. 2. Schematic diagram of the DPiMS[™] QT system.

3. PESI-OAD Synergistic Lipid Analysis

Check WP019

◆ OAD-MS/MS reveals double-bond (C=C) positions not accessible with CID-MS/MS. Atomic oxygen selectively oxidizes and cleaves at C=C.



Fig. 3. Typical OAD spectrum of a model lipid of PC (18:1).

WP 227

4. PESI-OAD Synergistic Alkaloids Analysis

+ Atomic oxygen selectively oxidizes carbon atoms adjacent to nitrogen in nitrogen-containing heterocycles, leading to OADspecific fragmentation distinct from CID, as shown in Fig. 4 (*m*/*z* 450.2127).



- ◆ OAD provides unique structural information on nitrogen-containing heterocycles, which is distinct from CID, particularly observed in alkaloids.

• OAD specific fragment ions observed in alkaloids offer the potential for rapid characterization of alkaloids. Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

