## Spatial Mapping of Lipids and Elements by Mass Microscopy and Integration with LA-ICP-MS in the **Diabetic Mouse Pancreata**

Koji Okuda<sup>1</sup>; Ei-ichi Matsuo<sup>1</sup>, Jeffrey Dahl<sup>2</sup>; Toshiya Matsubara<sup>2</sup>; Shinichi Yamaguchi<sup>1</sup>; Jannine I. Gamayot<sup>3</sup>; Shiuhwei Chen<sup>3</sup>; Christine M. Kusminski<sup>3</sup>; Philipp. E. Scherer<sup>3</sup>; <u>Ruth Gordillo<sup>3</sup></u> 1. Shimadzu Corporation, Kyoto, Japan; 2. Shimadzu Scientific Instruments, Columbia, MD; 3. University of Texas Southwestern Medical Center, Dallas, TX,

## Introduction.

Pancreas endocrine function critically relies on the islets of Langerhans that secrete insulin and glucagon as well as other hormones. Previous pancreas LC-MS/MS lipidomics studies required islets isolation by a lengthy process involving tissue digestion and manual selection, potentially compromising metabolite basal levels [1]. MALDI-TOF imaging allows for in situ global assessment of lipids in pancreata cryosections. Our methodology was applied to a mouse model of diabetes, the adipocytespecific-doxycycline-inducible mitochondrial ferritin (FtMT) overexpression, which displays massive beta-cell hyperplasia as its most striking phenotype [2].



Figure 1. A) Total pancreatic islet area (mm<sup>2</sup>) [2]. B) Insulin IF staining of pancreata and (left) and H&E staining of whole pancreas from WT and FtMT-Adip mice [2]. C) Mouse pancreas false color mass spectrometry imaging (25 μm spatial resolution) of SM 34:1;O<sub>2</sub> m/z 741.527 [M+K]<sup>+</sup>. Transgenic mouse pancreas displays considerably larger size islets. D) Unsupervised spatial segmentation calculation of MSI data (15 clusters) allows for differentiations of 🚺 µm spatial resolution. Bottom panel (cyan) negative ion mode analysis of WT pancreas cryosection islets and blood vessel regions. E) DAPI fluorescence staining of red marked area in D.



Imaging Mass Microscope iMScope<sup>™</sup> QT

## Methods and Instrumentation.

Frozen mouse pancreata was cryosectioned into 10 µm-thick sections. Tissue micro-sections were mounted on In<sub>2</sub>O<sub>3</sub>-SnO<sub>2</sub> (Millipore-Sigma, Burlington, MA). Automated matrix deposition system was achieved with an iMLayer<sup>™</sup> (Shimadzu Corporation, Japan). 2-mercaptobenzothiazole (MBT, positive ion mode) and 1,5diaminonapthalene (DAN, negative ion mode) were used as MALDI matrix reagents. Mass spectrometry imaging (MSI) analysis was performed on an iMScope<sup>™</sup> QT mass microscope (Shimadzu Corporation). Spatial elemental analysis was performed on an imageBIO266 (Elemental Scientific Lasers, USA) coupled to an ICPMS-2030 (Shimadzu Corporation). Imaging data processing was performed using IMAGEREVEAL<sup>TM</sup> MS software package (Shimadzu Corporation).

References. 1. Ye et al. J clin Invest 2018, 128(3):1178-1189 2. Kusminski et al. Diabetes 2020 69(3):313-330 3. MetaboAnalyst 5.0 www.metaboanalyst.ca

Acknowledgements. Shimadzu Corporation (JP) and Shimadzu Scientific Instruments (USA) for the collaborative efforts in mass spectrometry technology resources and useful discussion. Dr. José A. Fernández at the University of the Basque Country for useful discussion.





Figure 3. A) Mouse WT pancreas false color imaging mass spectrometry (7 µm spatial resolution) of SM 34:1;O2 m/z 741.527 [M+K]<sup>+</sup> **B)** Post mass spectrometry analysis glucagon IF.



iMScope, iMLayer, and IMAGEREVEAL are trademarks of Shimadzu Corporation in Japan and/or other countries. For Research Use Only. Not for use in diagnostic procedures. This publication may contain references to products that are not available in your country. Please contact us to check th availability of these products in your country.

Figure 5. A) LA-ICP-MS imaging analysis of mouse pancreas at 50 µm spatial resolution. Zn accumulates in the islets of Langerhans, Fe in blood vessels, while P is ubiquitous. B) Overlay of elemental images of Zn (green) and Fe (red) in control and FtMT overexpressing pancreas (analysis at 5μm spatial resolution). TG mice cryosection displays islets hyperplasia. High resolution optical images are also displayed. C) Elements relative abundance in islets and exocrine regions.



**UTSouthwestern** 

 $\bigoplus$ 

Medical Center<sub>®</sub>

SHIMADZU

Excellence in Science