## SHIMADZU

## Identifying Double-Bond-Positions of Phospholipids in Mouse Liver by Using Simultaneous Positive/Negative Ion Switching Analysis of LCMS-9050 and OAD-MS/MS

Yohei ARAO<sup>1</sup>, <u>Yuta MIYAZAKI<sup>1</sup></u>, Hidenori TAKAHASHI<sup>1</sup> (1) Shimadzu Corporation, Kyoto, Japan


### **1. Introduction**

- Lipids, a class of biomolecules, play a significant role in the physiological system. Non-targeted lipidomics by High Resolution Accurate Mass Spectrometer (HRAM) is currently one of the major analytical workflows for the clinical research such as biomarker discovery as well as the basic study. However, it has been difficult to identify double-bond positions in lipids using conventional methods.
- A brand-new oxygen attachment dissociation (OAD) provides diagnostic fragment ions enabling the assignment of double-bond positions. Identifying the double-bond positions by OAD-MS/MS is expected to contribute for understanding of the lipid metabolism. In this study, it was focused on the fragmentation derived from phosphatidylcholine (PC) by using OAD-MS/MS and identified the chemical structures.

#### 2. Methods

All experiments were performed using an ESI Q-TOF LCMS-9050 (Shimadzu, Kyoto, Japan). OH and O radicals generated by a compact microwave-driven radical source was introduced into the collision cell (Q2) through the quartz tube to obtain OAD-MS/MS spectrum <sup>[1][2]</sup>. Shimadzu Nexera<sup>™</sup> LC system was used for the separation of the lipids in the extracts of mouse liver.

Extracts were subjected to analysis with LCMS-9050, which allows simultaneous positive/negative ion switching analysis. The positive ion mode was used to generate OAD-MS/MS based spectra to determine the positions of the double-bonds, and spectra based on collision induced dissociation (CID) was generated to confirm the lipid classes and fatty acids.



The 160 mg of mouse-liver was homogenized using a multi-beads shocker with a metal cone at 1500  $\times$  g for 15 s, and 800  $\mu$ L of MeOH was added to the homogenate. After the solvent was homogenized again under the same conditions, an appropriate amount of MeOH (30 mg tissue weight per 200 µL) was transferred to a 2 mL glass tube. After the solvent was removed up to 95% of the total volume of the sample, 200  $\mu$ L of CHCl<sub>3</sub> was added to the solvent and vortexed for 10 s. After the solvent was incubated for 1 h on ice, 20 µL of ultra pure water was added, and the mixture was vortexed for 10 s. After 10 min incubation on ice, the solvent was centrifuged at 2000  $\times$  g for 10 min at 20 °C, and the supernatant was transferred to an LC-MS vial<sup>[2]</sup>.

Η
S
FI
G
M
M
C

6.0e5

5.0e5 -

### 3. Results

Two peaks at  $m/z = 782.5695 [M+H]^+$  from positive ion analysis mode and at 826.5604 [M+HCOO]<sup>-</sup> from negative ion analysis mode were detected from a mouse liver sample(Fig. 2) when Data Dependent Acquisition(DDA) with positive/negative ion switching analysis was performed. It was estimated as PC(36:4).

#### 2.2. Sample Preparation & Analytical Condition


**Table 1**Analytical Condition

**IPLC** condition ystem: Shimadzu Nexera X3 low rate: 0.3 mL/min Gradient Elution lobile phase A: 20mM Ammonium formate Nobile phase B ACN/IPA = 1/1 (v/v) Column: Shim-pack Scepter Claris C18-120 (1.9 μm, 2.1. mml.D. x 100 mmL.)

MS condition System: Shimadzu LCMS-9050 Mode: CID/OAD Polarity: Pos/Neg switching CE: 25V/ - 35V (pos/neg)


### 3.1. Detailed Processing of MS/MS spectra Based on CID

MS/MS spectra of PC(36:4) acquired with CID-based negative ion analysis were shown in Fig. 3.



**Fig. 3** MS/MS spectra of PC(36:4) acquired with CID-based negative ion analysis

**Table 2** Combination list of fatty acids in PC(36:4)







Carbon chain lengths and the number of double-bonds were estimated from CID-based negative ion analysis, and double-bond positions were identified with OAD-based positive ion analysis.

### **3.2. Detailed Processing of MS/MS spectra Based on OAD**


MS/MS spectra of PC(36:4) acquired with OAD-based negative ion analysis were shown in Fig. 4 and 5.







# **WP 520**




MS/MS spectra of PC(36:4) acquired with OAD-based positive ion analysis

**Fig. 5** Enlarged view of MS/MS spectra based on OAD. Left) peak A, Right) peak B

### 3.3. Identification of double bond positions

A combination of CID/OAD and positive/negative switching analysis by Q-TOF was used to identify lipids in mouse liver. The result was shown in Table 3. Three compounds were found from two peaks(A and B) detected at the same m/z.




#### 4. Conclusion

- identified.

### 5. Reference

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.



**Table 3** Structure identification result of PC(36:4)

✓ By using the combination of OAD and pos/neg switching analysis, the chemical structure of lipids in a mouse liver were

 $\checkmark$  It was suggested that alpha-linolenic acid and gammalinolenic acid are able to be distinguished by OAD.

[1] Takahashi.H et al. *Mass Spectrometry*. 2019, S0080. [2] Uchino.H et al. *Commun Chem.* 5, 162 (2022)