## SHIMADZU

## Synthetic modified oligonucleotides analysis using a matrix-assisted laser desorption/ionization digital-ion-trap mass spectrometer (MALDI-DIT-MS)

OYuko Fukuyama<sup>1</sup>, Hideharu Shichi, Masaki Murase, Yoshihiro Yamada, Sadanori Sekiya, Shinichi Iwamoto, Koichi Tanaka

### **1. Overview**

- In the analysis of oligonucleotide therapeutics, the development of analytical techniques using LC-MS is accelerating, while the need for simpler and faster analytical techniques is also increasing.
- We report here a simple and rapid analytical method of synthetic modified oligonucleotides using MALDI-digital ion trap-MS (MALDI-DIT-MS).

## 2. Methods

#### 2-1. Analytes (model oligonucleotide samples)

• **Single-stranded DNA**: Mipomersen (GeneDesign Inc.), MW 7177

5' MG-MC-MC-MU-MC-dA-dG-dT-dC-dT-dG-dC-dT-dC-MG-MC-MA-MC-MC 3

A, C, G, and U represent adenosine, cytidine, guanosine, and uridine deoxyribonucleotides.; M: 2'-O-(2-methoxyethyl) nucleotide d: 2'-deoxynucleotide; Substitution at 5-position of cytosine and uracil base with a methyl group.

• Single-stranded RNA: Patisiran, sense strand (GeneDesign Inc.), MW 6764 5' G-Um-A-A-Cm-Cm-A-A-G-A-G-Um-A-Um-Um-Cm-Cm-A-Um-dT-dT 3'

• **Double-stranded RNA**: Vutrisiran (MedChemExpress LLC),

MW 16345 (sense strand: 8789; antisense strand: 7558) 5' Ums-Gms-Gm-Gm-Am-Um-Uf -Um-Cf - Af - Uf - Gm-Um-Am-Am-Cm-Cm-Am-Am-Gm- Am-R1 3' 3' Cms-Ums-Am - Cm - Cm-Cm-Um -Af- Am-Af -Gm-Um-Am-Cm- Af - Um-Um-Gf -Gm-Um-Um-sCf-sUm 5'



A and G represent adenosine and guanosine ribonucleotides. Af: 2'-fluoroadenosine; Am: 2'-O-methyladenosine; Cf: 2'-fluorocytidine; Cm: 2'-O-methylcytidine; Gf: 2'-fluoroguanosine; Gm: 2'-O-methylguanosine; Uf: 2'-fluorouridine; Um: uracil 2'-O-methyluridine; dT: thymidine deoxyribonucleotide; R1: triantennary GalNAc, s: Phosphodiester bonds between nucleotides are substituted with phosphorothioate bonds.

### **2-2.** Preparation

- 1. An aqueous solution of the analyte and a matrix solution: 40 mg/mL 50% ACN/water (v/v) with diammonium hydrogen citrate were prepared.
- 2. After mixing the analyte solution and the matrix solution at 1:1 (v/v), 1  $\mu$ L was dropped on the sample plate.

### **2-3. MALDI-MS measurement**

- Instrument: MALDImini<sup>TM</sup>-1 (Shimadzu Corp.) (Fig. 1)
- Measurement: Raster scan, Positive ion mode



Koichi Tanaka Mass Spectrometry Research Laboratory, Shimadzu Corporation, Kyoto, Japan.



and/or other countries. The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in

[1] McLuckey SA, J. Am. Soc. Mass. Spectrom., 1992, 3, 60-70.

diagnostic procedures.

 $\mathbf{w}_2 \mathbf{x}_2 \mathbf{y}_2 \mathbf{z}_2 \mathbf{w}_1 \mathbf{x}_1 \mathbf{y}_1 \mathbf{z}_1$ 

Nucleic acid fragment ion name [1]

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

# **ThP-571**

Single-stranded RNA was detected as [M+H]<sup>+</sup> with high sensitivity under the analytical conditions in Table 3.

Table 3. MW analysis conditions		
for single-stranded	I RNA	_
Instrumental:		
Dynode voltage (V)	7000	
Detector voltage (V)	1300	
RF delay ( <i>µ</i> s)	25	-
Others:		
Matrix (ex.)	3-HPA/2,4,6-THAP	for MV
ACD (concentration)	70 mM	-

Base sequence information including the modification was obtained under the analytical conditions in Table 4. In this case, d/wions were mainly detected.

 
 Table 4.
 Sequence analysis
conditions

#### Instrumental:

		_
Dynode voltage (V)	8000	High
Detector voltage (V)	1800	High
RF delay (μs)	17	Low
Others:		
Matrix (ex.)	2,4-DHAP	for seq
ACD (concentration)	70 mM	

**Fig. 6** Mass spectra of patisiran, sense strand (10 pmol/well)

We confirmed that molecular weight analysis, base sequence analysis including modification sites, and terminal GalNAc modification analysis for single-stranded DNA/RNA and doublestranded RNA are possible with a single MALDI-DIT-MS instrument.



