

MS Imaging Spectrum Annotation Tool using Accurate Mass and Isotope Pattern Analysis

Jeff Dahl and Toshiya Matsubara
Shimadzu Scientific Instruments, Columbia, Maryland

Introduction

MS Imaging reveals the spatial distribution of important biomolecules to discover the details of disease processes, find drug distribution patterns, or characterize new materials. Real-world samples can result in complex mass spectra which are time-consuming to process manually. Searching for targets may turn up false matches, even having less than one part-per-million mass error. Using isotope pattern can exclude false matches and speed up processing. We developed a new MS spectrum annotation tool using accurate mass and isotope pattern analysis to quickly match targets to data.

Methods

Tissue samples were prepared by vapor-depositing 1,5-Diamino-naphthalene onto mouse pancreas affixed to conductive glass slides using a Shimadzu iMLayer™ sublimation unit. Then tissues were analyzed by MALDI MS imaging using a Shimadzu IMScope QT™, composed of an ionization unit connected to a QToF mass spectrometer. Spectra were exported from the acquisition software and analyzed by a newly developed data processing tool.

Setting	Value
Pitch X	30.0 µm
Pitch Y	30.0 µm
Polarity	Negative
Scan Range (m/z)	500 to 1000
Number of Pixels	8752
Laser Irradiation Number	80 shots
Laser Repetition Frequency	2000 Hz
Laser Diameter	Approximately 35 µm
Laser Intensity	65%
DL Temperature	250 deg C

Table 1. Imaging MS settings optimized for mouse pancreas tissue, DAN coated, in negative MS mode.




Figure 1. MS images of mouse pancreas: TIC, representative lipids with acceptable matching criteria, and optical image. Alpha cells (round clusters, one example marked α) are visible in the optical and selected MS images.

Results

Mass spectra were imported into the newly developed tool along with a list of targets of interest which might be expected to be detected in the sample. The settings in the tool allow to select ionization mode, adducts, mass error and intensity tolerances, and other match settings. For each target with a matching signal in the MS spectrum, the isotope pattern match, monoisotopic ion (A-number) and other characteristics were calculated and shown in a summary table.




Figure 2. (Left) Concept of A-number assignment. (Right) Workflow of data processing using the tool.

1. Import target list (name + formula)
• Same list format as ImageReveal
2. Import spectrum data
• Same as spectrum export format of ImageReveal
3. Choose settings like adducts, mass error tolerance etc.
4. Filter and sort target matches which are color coded based on mass error, isotope score, A-Number etc.

PI	Isotop	Protein	Mass	Intensity	PI Error	A-Number	Score	IsotopeScore	Formula
PI 36:1	M+H	863.5499	863.5499	11700	0.0	1	100	100	C42H81O10P
PI 36:2	M+H	865.5499	865.5499	14715	0.0	2	100	100	C42H81O10P
PI 36:3	M+H	867.5499	867.5499	17730	0.0	3	100	100	C42H81O10P
PI 36:4	M+H	869.5499	869.5499	20745	0.0	4	100	100	C42H81O10P
PI 36:5	M+H	871.5499	871.5499	23760	0.0	5	100	100	C42H81O10P
PI 36:6	M+H	873.5499	873.5499	26775	0.0	6	100	100	C42H81O10P
PI 36:7	M+H	875.5499	875.5499	29790	0.0	7	100	100	C42H81O10P
PI 36:8	M+H	877.5499	877.5499	32805	0.0	8	100	100	C42H81O10P
PI 36:9	M+H	879.5499	879.5499	35820	0.0	9	100	100	C42H81O10P
PI 36:10	M+H	881.5499	881.5499	38835	0.0	10	100	100	C42H81O10P
PI 36:11	M+H	883.5499	883.5499	41850	0.0	11	100	100	C42H81O10P
PI 36:12	M+H	885.5499	885.5499	44865	0.0	12	100	100	C42H81O10P
PI 36:13	M+H	887.5499	887.5499	47880	0.0	13	100	100	C42H81O10P
PI 36:14	M+H	889.5499	889.5499	50895	0.0	14	100	100	C42H81O10P
PI 36:15	M+H	891.5499	891.5499	53910	0.0	15	100	100	C42H81O10P
PI 36:16	M+H	893.5499	893.5499	56925	0.0	16	100	100	C42H81O10P
PI 36:17	M+H	895.5499	895.5499	59940	0.0	17	100	100	C42H81O10P
PI 36:18	M+H	897.5499	897.5499	62955	0.0	18	100	100	C42H81O10P
PI 36:19	M+H	899.5499	899.5499	65970	0.0	19	100	100	C42H81O10P
PI 36:20	M+H	901.5499	901.5499	68985	0.0	20	100	100	C42H81O10P
PI 36:21	M+H	903.5499	903.5499	72000	0.0	21	100	100	C42H81O10P
PI 36:22	M+H	905.5499	905.5499	75015	0.0	22	100	100	C42H81O10P
PI 36:23	M+H	907.5499	907.5499	78030	0.0	23	100	100	C42H81O10P
PI 36:24	M+H	909.5499	909.5499	81045	0.0	24	100	100	C42H81O10P
PI 36:25	M+H	911.5499	911.5499	84060	0.0	25	100	100	C42H81O10P
PI 36:26	M+H	913.5499	913.5499	87075	0.0	26	100	100	C42H81O10P
PI 36:27	M+H	915.5499	915.5499	90090	0.0	27	100	100	C42H81O10P
PI 36:28	M+H	917.5499	917.5499	93105	0.0	28	100	100	C42H81O10P
PI 36:29	M+H	919.5499	919.5499	96120	0.0	29	100	100	C42H81O10P
PI 36:30	M+H	921.5499	921.5499	99135	0.0	30	100	100	C42H81O10P
PI 36:31	M+H	923.5499	923.5499	102150	0.0	31	100	100	C42H81O10P
PI 36:32	M+H	925.5499	925.5499	105165	0.0	32	100	100	C42H81O10P
PI 36:33	M+H	927.5499	927.5499	108180	0.0	33	100	100	C42H81O10P
PI 36:34	M+H	929.5499	929.5499	111195	0.0	34	100	100	C42H81O10P
PI 36:35	M+H	931.5499	931.5499	114210	0.0	35	100	100	C42H81O10P
PI 36:36	M+H	933.5499	933.5499	117225	0.0	36	100	100	C42H81O10P
PI 36:37	M+H	935.5499	935.5499	120240	0.0	37	100	100	C42H81O10P
PI 36:38	M+H	937.5499	937.5499	123255	0.0	38	100	100	C42H81O10P
PI 36:39	M+H	939.5499	939.5499	126270	0.0	39	100	100	C42H81O10P
PI 36:40	M+H	941.5499	941.5499	129285	0.0	40	100	100	C42H81O10P
PI 36:41	M+H	943.5499	943.5499	132300	0.0	41	100	100	C42H81O10P
PI 36:42	M+H	945.5499	945.5499	135315	0.0	42	100	100	C42H81O10P
PI 36:43	M+H	947.5499	947.5499	138330	0.0	43	100	100	C42H81O10P
PI 36:44	M+H	949.5499	949.5499	141345	0.0	44	100	100	C42H81O10P
PI 36:45	M+H	951.5499	951.5499	144360	0.0	45	100	100	C42H81O10P
PI 36:46	M+H	953.5499	953.5499	147375	0.0	46	100	100	C42H81O10P
PI 36:47	M+H	955.5499	955.5499	150390	0.0	47	100	100	C42H81O10P
PI 36:48	M+H	957.5499	957.5499	153405	0.0	48	100	100	C42H81O10P
PI 36:49	M+H	959.5499	959.5499	156420	0.0	49	100	100	C42H81O10P
PI 36:50	M+H	961.5499	961.5499	159435	0.0	50	100	100	C42H81O10P
PI 36:51	M+H	963.5499	963.5499	162450	0.0	51	100	100	C42H81O10P
PI 36:52	M+H	965.5499	965.5499	165465	0.0	52	100	100	C42H81O10P
PI 36:53	M+H	967.5499	967.5499	168480	0.0	53	100	100	C42H81O10P
PI 36:54	M+H	969.5499	969.5499	171495	0.0	54	100	100	C42H81O10P
PI 36:55	M+H	971.5499	971.5499	174510	0.0	55	100	100	C42H81O10P
PI 36:56	M+H	973.5499	973.5499	177525	0.0	56	100	100	C42H81O10P
PI 36:57	M+H	975.5499	975.5499	180540	0.0	57	100	100	C42H81O10P
PI 36:58	M+H	977.5499	977.5499	183555	0.0	58	100	100	C42H81O10P
PI 36:59	M+H	979.5499	979.5499	186570	0.0	59	100	100	C42H81O10P
PI 36:60	M+H	981.5499	981.5499	189585	0.0	60	100	100	C42H81O10P
PI 36:61	M+H	983.5499	983.5499	192600	0.0	61	100	100	C42H81O10P
PI 36:62	M+H	985.5499	985.5499	195615	0.0	62	100	100	C42H81O10P
PI 36:63	M+H	987.5499	987.5499	198630	0.0	63	100	100	C42H81O10P
PI 36:64	M+H	989.5499	989.5499	201645	0.0	64	100	100	C42H81O10P
PI 36:65	M+H	991.5499	991.5499	204660	0.0	65	100	100	C42H81O10P
PI 36:66	M+H	993.5499	993.5499	207675	0.0	66	100	100	C42H81O10P
PI 36:67	M+H	995.5499	995.5499	210690	0.0	67	100	100	C42H81O10P
PI 36:68	M+H	997.5499	997.5499	213705	0.0	68	100	100	C42H81O10P
PI 36:69	M+H	999.5499	999.5499	216720	0.0	69	100	100	C42H81O10P
PI 36:70	M+H	1001.5499	1001.5499	219735	0.0	70	100	100	C42H81O10P
PI 36:71	M+H	1003.5499	1003.5499	222750	0.0	71	100	100	C42H81O10P
PI 36:72	M+H	1005.5499	1005.5499	225765	0.0	72	100	100	C42H81O10P
PI 36:73	M+H	1007.5499	1007.5499	228780	0.0	73	100	100	C42H81O10P
PI 36:74	M+H	1009.5499	1009.5499	231795	0.0	74	100	100	C42H81O10P
PI 36:75	M+H	1011.5499	1011.5499	234810	0.0	75	100	100	C42H81O10P
PI 36:76	M+H	1013.5499	1013.5499	237825	0.0	76	100	100	C42H81O10P
PI 36:77	M+H	1015.5499	1015.5499	240840	0.0	77	100	100	C42H81O10P
PI 36:78	M+H	1017.5499	1017.5499	243855	0.0	78	100	100	C42H81O10P
PI 36:79	M+H	1019.5499	1019.5499	246870	0.0	79	100	100	C42H81O10P
PI 36:80	M+H	1021.5499	1021.5499	249885	0.0	80	100	100	C42H81O10P
PI 36:81	M+H	1023.5499	1023.5499	252900	0.0	81	100	100	C42H81O10P
PI 36:82	M+H	1025.5499	1025.5499	255915	0.0	82	100	100	C42H81O10P
PI 36:83	M+H	1027.5499	1027.5499	258930	0.0	83	100	100	C42H81O10P
PI 36:84	M+H	1029.5499	1029.5499	261945	0.0	84	100	100	C42H81O10P
PI 36:85	M+H	1031.5499	1031.5499	264960	0.0	85	100	100	C42H81O10P
PI 36:86	M+H	1033.5499	1033.5499	267975	0.0	86	100	100	C42H81O10P
PI 36:87	M+H	1035.5499	1035.5499	270990	0.0	87	100	100	C42H81O10P
PI 36:88	M+H	1037.5499	1037.5499	274005	0.0	88	100	100	C42H81O10P
PI 36:89	M+H	1039.5499	1039.5499	277020	0.0	89	100	100	C42H81O10P
PI 36:90	M+H	1041.5499	1041.5499	280035	0.0	90	100	100	C42H81O10P
PI 36:91	M+H	1043.5499	1043.5499	283050	0.0	91	100	100	C42H81O10P
PI 36:92	M+H	1045.5499	1045.5499	286065	0.0	92	100	100	C42H81O10P
PI 36:93	M+H	1047.5499	1047.5499	289080	0.0	93	100	100	C42H81O10P
PI 36:94	M+H	1049.5499	1049.5499	292095	0.0	94	100	100	C42H81O10P
PI 36:95	M+H	1051.5499	1051.5499	295110	0.0	95	100	100	C42H81O10P
PI 36:96	M+H	1053.5499	1053.5499	298125	0.0	96	100	100	C42H81O10P
PI 36:97	M+H	1055.5499	1055.5499	301140	0.0	97	100	100	C42H81O10P
PI 36:98	M+H	1057.5499	1057.5499	304155	0.0	98	100	100	C42H81O10P
PI 36:99	M+H	1059.5499	1059.5499	307170	0.0	99	100	100	C42H81O10P
PI 36:100	M+H	1061.5499	1061.5499	310185	0.0	100	100	100	C42H81O10P

