# **SHIMADZU**

# **Determination of Paraquat in Drinking Water by LCMS-2050**

<u>Marilia Santoro Cardoso<sup>1</sup></u>; Marcos A Pudenzi<sup>1</sup>; Luisa B Ambrosio<sup>1</sup>; Ichiro Hirano<sup>1</sup> <sup>1</sup>Shimadzu do Brasil, Barueri, Brazil

#### 1. Introduction

- Paraquat (Fig. 1) is an herbicide widely used in the world and extremely toxic for humans and animals, capable of controlling the weeds growth in crops.<sup>1</sup>
- According with Portaria GM/MS No. 888 (2021) of the Ministry of Health in Brazil, which establish norms and water potability standards, the maximum value allowed of paraquat in drinking water is 13 μg/L.<sup>2</sup> The aim of this work is to develop a method by direct analyses of paraquat in drinking water using a single quadrupole (LCMS-2050).



**Fig. 1** Molecular structure of paraquat.

#### 2. Methods

Analytical chromatography was performed using an LC-MS (Nexera XSi - LCMS-2050, Shimadzu®, Japan) with the following conditions of HPLC and MS (Table 1 and 2, respectively).

**Table 1**Method parameters of HPLC system.

Column	Zic®-Hilic (100 x 2.1mm x 3.5 µm). Merck - PN:
	1504410001
Mobile Phase	A: Water + 20mM Ammonium formate + 0.5%
	Formic acid
	B: Acetonitrile + 0.3 % formic acid
Flow and Column Temp.	0.4 mL/min - 40 °C
Injection Volume	1 µL



Fig. 2 LCMS 2050 single quadrupole.

Interface Voltage	0.5 kV
Interface	ESI +
Nebulizing Gas Flow	2 L/min
Heating Gas Flow	7 L/min
Drying Gas Flow	4 L/min
Dessolvation Temp.	450 °C
DL Temp.	200 °C

**Table 2** Method parameters of MS system.

## 3. Results

- Method optimization was performed evaluating the analyte intensity with different mobile phases and interface parameters, such as adjust of m/z, interface voltage, nebulizing gas flow, drying gas flow and heating gas flow, in this order (Graphic 1 a-e).
- ◆ Initially, with the aim of increasing the analyte intensity was injected a standard of paraquat at a concentration of 100 ng/mL in different organic mobile phases: acetonitrile versus formic acid (0.3%) in acetonitrile. As a result, a paraquat peak with greater intensity was observed using acetonitrile as mobile phase. Then, a standard of the same concentration was injected in different aqueous mobile phases: ammonium formate (20mM) and formic acid (0.5%), and formic acid (0.5%). As a result, it was observed a better reproducibility using ammonium formate in the phase.
- Following optimizing the chromatography parameters, an optimization of mass spectrometry parameters was performed. According to the results obtained, interface parameters stablished were adjust of m/z was set at 93.1, interface voltage at 0.5 kV, nebulizing gas flow at 2 L/min, drying gas flow at 4 L/min, heating gas flow at 7 L/min, desolvation temperature at 450 °C and desolvation line temperature at 200 °C.



**Graphic.1** Method optimization results using the analyte intensity evaluating interface parameters (adjust of m/z, interface voltage, nebulizing gas flow, drying gas flow and heating gas flow)

- ♦ After mass spectrometry optimization, a triplicate linearity range in water was evaluated by SIM (single ion monitoring) mode with the selected m/z. The method demonstrated to be linear through the calibration curve with a r<sup>2</sup> of 0,99. The results are presented in Fig. 3, which displays the calibration curve and chromatographic peaks at limit inferior of quantification (LIQ) in ultrapure water samples.
- Finally, a triplicate linearity range in real sample (tap water) was evaluated, resulting in the determination of paraquat in the limit inferior of quantification with a r<sup>2</sup> of 0,99 in the calibration curve. In Fig. 4 was presented the chromatographic peaks at LIQ and blank in real samples (tap water).

# **ThP 590**

![](_page_0_Figure_24.jpeg)

![](_page_0_Figure_25.jpeg)

![](_page_0_Figure_26.jpeg)

**Fig. 4** Chromatographic peaks at limit of quantification in and blank in real samples (tap water).

### 4. Conclusion

It was possible to obtain satisfactory results for the determination of Paraquat in tap water using a mass spectrometry single quadrupole, LCMS-2050, which facilitates the work routine and reduces analysis costs. The results obtained reached the limits lower than those recommended in Ordinance GM/MS n° 888 of May 4, 2021, without the need for pre-treatment of the samples.

#### Reference

2) Gabinete do Ministro. Portaria Nº 888, de 04 de maio de 2021.

#### Disclaimer:

The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

The authors are affiliated and funded by Shimadzu do Brasil..

<sup>1)</sup> Martins T. (2013). Semin. Cienc. Biol. Saude. 34(2):175-86.