

Aza-Prilezhaev Aziridination-Enabled Multidimensional Analysis of Isomeric Lipids via High-Resolution U-Shaped Mobility Analyzer-Mass Spectrometry

Yuling Li¹, Yiming Wang¹, Kang Guo¹, Kuofeng Tseng¹, Xiaoqiang Zhang¹, Wenjian Sun¹ (1) Shimadzu Research Laboratory (Shanghai) Co., Ltd.

1. Introduction

- While derivatization-assisted liquid chromatography tandem mass spectrometry remains the gold standard technique in lipidome, it mainly faces challenges in efficiently labeling carbon-carbon double bond (C=C) and differentiating isometric lipids in full dimension.
- We have established an entire workflow by coupling a metal- and additive-free aza-Prilezhaev aziridination (APA)^[1] with U-shaped ion mobility analyzer (UMA)^[2]-MS/MS for comprehensive structural annotation of polar/nonpolar isomeric lipids, achieving high level in multiple aspects simultaneously (sensitivity, specificity and efficiency) (Fig. 1). The user-friendly APA with high conversion, minimal side reaction and good substrate generality is developed to improve the signal intensity by up to 3 orders of magnitude. This method was successfully used to identify C=C location, fatty acyl length and sn-position, headgroup composition and geometric configuration in a single run^[3].





2. Methods

- Lipid solutions were simply mixed with superstoichiometric APA reagent. Without complicated work-up, the final reaction solution was quenched and diluted by acetonitrile prior to UMA-QTOF MS/MS analysis (Fig. 2).
- Two operating modes: Offline & Online

Multiple analyzing modes: MS²: UMA \rightarrow Quad mass selection \rightarrow CID \rightarrow TOF Pseudo MS³: UMA \rightarrow SMP for 1st fragmentation \rightarrow Quadrupole mass selection \rightarrow c.c. for 2nd fragmentation \rightarrow TOF



Workflow of this research. SMP, segmented Fig. 2 multipole; c.c., collision cell; HG, headgroup

3. Results

- ◆ Compared with Prilezhaev epoxidation^[4], APA has higher sensitivity, higher reaction selectivity, less spectrum complexity and better substrate generality (Fig. 3).
- ◆ Nearly full conversion was obtained at 50 °C for about 10 min. Stable mono-aziridination product is dominant, simplifying the spectrum for polyunsaturated lipids.
- Protophilic group was efficiently introduced for unsaturated lipids via APA. The LODs of lipids were decreased by 1 ~ 3 orders of magnitude, achieving nM detection limit.



- products of polyunsaturated lipid. • For shotgun analysis, UMA was successfully applied to roughly distinguish different classes and subclasses of lipids (Fig. 4A-B),
- reduce isotopic type II effect (Fig. 4C-E) and separate geometric isomers (Fig. 4K).
- \blacklozenge Abundant site-specific fragmentation ions indicate the C=C location of C18:1 in MS/MS spectra (Fig. 4F-H). Most classes of lipids were also compatible.
- ◆ sn-Positional isomers were MS² CID of sodiated aziridine-PC 18:1 $(\Delta 9)/16:0$ standard produced sn-1 diagnostic ion at m/z 360.2873 (Fig. 4I), while *sn*-2 diagnostic ion at *m/z* 290.2454 was generated from aziridine-PC 16:0/18:1 (Δ 9) (**Fig. 4J**).





Fig. 4

- level.

Detailed structural analysis of isomeric lipids by APA-UMA-MS/MS. EIM of (A) different classes and (B) subclasses of aziridine-lipids. (C) MS spectrum of a 10:1 mixture of PC 34:2 and PC 34:1. (D) EIMs of m/z 775.5960. (E) MS spectrum isolated from E of 3.60-3.6640 V/mm MS/MS spectra of (F) aziridine-C18:1 (Δ 6), (G) aziridine-C18:1 (Δ 9), (H) aziridine-C18:1 (Δ 11), (I) aziridine-PC 18:1 (9Z)/16:0 and (J) aziridine-PC 16:0/18:1 (9Z), red trace represents C=C positional diagnostic ions, while blue trace represents sn-position diagnostic ions. (K) EIM of sodiated aziridine-C18:1 (9Z/E)

• We established an UMA-based pseudo MS³ strategy to further site-specific enhance fragment regioisomeric of glycerophospholipids. Compared with MS/MS (Fig. 5A-B), the abundance of glycerol backbone structural ions originating from cross-ring cleavage across the dioxolane moiety are improved by around 1 order of magnitude using pseudo MS³ (Fig. 5C-D). In addition, unlike in methods with multiple stage low resolution MS such as ion trap, high resolution TOF provides higher confidence for both precursor and fragment identification of lipids.

For real sample analysis, APA combined with mobility separation is feasible to reduce matrix effect. A total of 70 low polar lipids were detected in olive oil extracts. Without hydrolysis and complicated work-up, considerable number of triacylglycerols, diacylglycerols and sterol esters were identified at deep structure



Fig. 5	MS
_	az
	az
	Pr

4. Conclusion

Reference

- 1) Synthesis (2022) 54, 4513-4520. Jat *et al.*

- Disclaimer

The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures. The authors declare no competing financial interest.

MP 477

IS/MS spectra of (A) aziridine-PC 18:1/16:0 and (B) ziridine-PC 16:0/18:1; Pseudo MS³ spectra of (C) ziridine-PC 18:1/16:0 and (D) aziridine-PC 16:0/18:1; (E) roposed fragmentation scheme of aziridine-PC 16:0_18:1.

APA considerably promotes ionization of lipids with wide polarity and abundance range in positive ESI mode.

Following high-throughput separation of geometric isomers, highly specific cleavage can be performed to generate abundant C=C location and *sn*-position diagnostic ions.

◆ UMA-based pseudo MS³ mode allows enhancement of the backbone structure-specific fragmentation, thus further improving detection sensitivity and increasing lipid coverage.

Development of an integrated platform LC-APA-UMA-MS/MS and corresponding applications are undergoing.

2) Anal. Chem. (2020) 92, 8356-8363. Wang *et al.* 3) Anal. Chem. (2024) 96, 7111-7119. Li *et al.* 4) Anal. Chem. (2019) 91, 11905-11915. Kuo et al.



