



Improving LC-MS Analysis of mRNA: Enzymatic Digestion and LC-MS Analysis for Characterization of mRNA Modifications and Sequence Integrity

Pierre-Alexander Mücke¹; Anke Schnabel¹; Daniel Waldera-Lupa¹; Sascha Rexroth² (1) ProtaGene GmbH, Dortmund, Germany; (2) Shimadzu Europa GmbH, Duisburg, Germany

I. Introduction

- > The rapid advancements in mRNA synthesis and delivery technologies have enabled the development of various mRNAbased vaccines and therapeutics, offering significant benefits in combating human diseases¹⁾. RNA modifications are crucial for ensuring the efficacy of in vitro transcribed therapeutic mRNAs, as enhance translation and stability while reducing they immunogenicity²⁾.
- \succ This study shows the integrated analysis of the mRNA encoding the SARS-CoV-2 (COVID) spike protein based on enzymatic digestion with site-specific endoribonucleases to generate unique oligonucleotides, which can be mapped to the expected mRNA sequence during data analysis (Fig. 1).



2. Methods

Fig. 1 Workflow of this research

Based on an in-silico digest of mRNA encoding the COVID spike protein with various RNases (T1, MC1, Cusativin, 4, A, U2), most appropriate RNases were selected for enzymatic reactions. Inactivation of RNases was performed by magnetic bead removal (RNase T1) or heat (RNase MC1).

The resulting oligonucleotides were separated using ion-pair reversed-phase HPLC with 0.1% TEA and 1% HFIP (Shimadzu Nexera XSi) and introduced into a Q-TOF mass spectrometer (Shimadzu LCMS-9050) running in DDA negative mode. MSconditions were adjusted for effective ion transfer and fragmentation of oligos of about 10-30 nucleotides in length.

Data analysis for confident oligo MS/MS identification and mRNA sequence confirmation was performed with Protein Metrics Byos using the Digested Oligos Workflow. 2'3' cyclic phosphate modifications were allowed as variable modifications at the 3' terminal.



- theoretically generates most oligos with preferred length (Fig. 2)
- \succ As the second enzyme an RNase T1 magnetic bead complex was selected over RNase 4 because of the possibility to easily remove the enzyme from the digestion mix and thus stop the enzymatic reaction



- oligonucleotide backbone (Fig. 4).
- \succ Our data indicate that these optimized conditions resulted in a 30% increase in sequence coverage. Furthermore, the systematic adjustment of CE reduced non-specific fragmentation occurrences, leading to cleaner spectra and fewer ambiguous

TP 599





- > In order to verify the successful capping of the mRNA, XICs were generated for the expected m/z values considering the enzymes used (**Fig. 6**)
- \succ In the T1 digested sample the Cap1 structure could be identified via specific precursor and fragment ion masses

4. Conclusion

- \succ Our findings underscore the significance of tailored enzymatic digestion and dynamic CE optimization in enhancing the quality of LC-MS analysis of mRNA.
- \succ Overall, our approach effectively addresses the technical challenges associated with full-length mRNA characterization and significantly improves both sequence coverage and the detection of modifications, enabling more comprehensive and reliable analysis.

Reference

- 1) The Lancet (2024), Vol. 403, 1192–1204, Parhiz & Atochina-Vasserman, et al.
- 2) Nucleic Acids Research (2016) Vol. 44, 7511–7526, Ramanathan, *et al.*

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

Sascha Rexroth is affiliated and funded by Shimadzu Corporation.

Pierre-Alexander Mücke, Anke Schnabel and Daniel Waldera-Lupa are affiliated and funded by ProtaGene GmbH.