

# Benchtop Matrix Assisted Laser Desorption Ionization Mass Spectrometry of PFAS for preliminary screening of environmental samples

Caroline Jones<sup>1</sup>, Tom Abban<sup>1</sup>, MilaanThirukumaran<sup>2</sup>, Francine Yanchik-Slade<sup>2</sup>, Matthew E. Openshaw<sup>1</sup>  
(1) Shimadzu, Manchester, UK (2) Shimadzu Scientific Instruments, Columbia, USA

## 1. Introduction

- ◆ Per- and polyfluoroalkyl substances (PFAS) are persistent environmental contaminants which negatively impact human health and the environment.
- ◆ High concentrations of PFAS in test samples can cause contamination of sensitive LC-MS/MS instruments requiring extensive decontamination procedures, resulting in significant downtime and costly delays in sample processing.
- ◆ We present the development of a simple, cost-effective method for the rapid screening of PFAS samples using an entry-level linear benchtop MALDI-TOF mass spectrometer to identify samples containing high levels before they enter the analytical workflow.



**Fig. 1** Analysis Conditions for MALDI-8030 analysis of PFAS

System	: MALDI-8030
Polarity	: Negative
Mass Range	: $m/z$ 100-1000
Acquisition	: 980 shots @ 200Hz
Blanking	: 200
Pulsed Extraction	: 400

## 2. Methods


PFOA and PFOS standards were purchased from Merck, UK. Stock solutions (1 ppm prepared in acetonitrile) were diluted in water to the required concentrations.

Samples were tested with a variety of MALDI matrices (CHCA, DHB, DHAP, DAN, NRM and TMGN), with norharmane (NRM) and 1,8-bis(tetramethylguanidino)-naphthalene (TMGN) identified for further investigation. Diluted samples were mixed with NRM (1 mg/mL) or TMGN (5 mg/mL) in 70:30 Acetonitrile:Water. Blank samples of UHP water were prepared to ensure no contamination was present. 1  $\mu$ L of the mixture was spotted onto a FlexiMass-SR48 slide and the samples were analysed using a MALDI-8030 MALDI-TOF mass spectrometer (Shimadzu). Analysis conditions are shown in Fig 1. Peaks at  $m/z$  314.2 and  $m/z$  499.2 were monitored for PFOA and PFOS, respectively.  $^{13}\text{C}_8\text{PFOS}$  ( $m/z$  507) was used as an internal standard for PFOS samples with a final sample concentration of 6.25 ng/mL (see Fig. 2).

MALDI-MSI was used to assess spot homogeneity and limits of detection were established.

Subsequently, PFHxA and PFBS standards were purchased from Merck, UK and analysed using the same method.


## 3. Results




**Fig. 2** Addition of internal standard and matrix is all that is required for analysis

- ◆ With both norharmane and TMGN, PFOS and PFOA were detected at concentrations suitable for preliminary screening to prevent contamination of LC-MS/MS systems<sup>1</sup>. MALDI-MSI spot imaging revealed the migration of PFAS to the edges of the spots during spotting/drying (Fig. 3). Use of an annular acquisition strategy for norharmane spots produced spectra with a higher signal to noise, as expected. Spots prepared with TMGN matrix showed a more homogeneous distribution and work was continued with this matrix (5 mg/mL in 70% acetonitrile).

- ◆ Limits of detection were established for PFOA, PFOS, PFHxA and PFBS standards at 5 ppb, 250 ppt, 200 ppt and 1 ppb respectively (Fig. 4).




**Fig. 3** MALDI-MSI Spot imaging of PFOS and PFOA in Norharmane & TMGN



**Fig. 4** Representative spectra showing LODs for PFOA (5 ppb), PFOS (250 ppt), PFHxA (200 ppt) and PFBS (1 ppb)

- ◆ Analysis of bottled water samples were chosen to be representative of samples received within a PFAS testing environment. A small number of samples showed low levels of possible PFAS present in the sample (Fig. 5).



**Fig. 5** Representative spectra showing analysis of bottled water alongside a mixed PFAS sample containing PFOS (50 ppb), PFOA (100 ppb), PFBS (50 ppb) and PFHxA (100 ppb) and a blank UHP water sample.


- ◆ To demonstrate how MALDI-TOF MS could be applied in a screening workflow, whole slide imaging was performed demonstrating a quick and simple method to visually identify potentially problematic samples.

- ◆ Whole slide imaging (250  $\mu$ m spacing) of a full slide of standard samples (48 spots) with mixed concentration levels was completed in approx. 12 mins.

- ◆ Mixed concentration samples were spotted onto a standard 48 well FlexiMass-SR48 slide and also a 96 well FlexiFocus™ slide for imaging

- ◆ IMAGEREVEAL MS™ (Shimadzu, Japan) software was used to create the slide ion image using a red-green intensity colour bar to display high concentration samples in green and low concentrations in red (see Fig. 7).

- ◆ 100 ng/mL was chosen as a suitable cut off for identifying problematic PFOS samples. Setting the intensity scaling for 100 ng/mL samples as mid-range showing a mixture of red and green, allowed easy identification of high concentration samples as bright green. (See Fig. 7).



**Fig. 7** Whole plate imaging as a rapid screening tool. In the examples shown, PFOS samples at 0, 1, 20, 50, 100 & 200 ng/mL were analysed after spotting onto a) a 96 well target or b) a 48 well target.

## 4. Conclusion

- ◆ We have shown, in principle, that it is possible to use MALDI-TOF mass spectrometry to implement a simple pre-screening method for aqueous samples to prevent widespread contamination in a PFAS monitoring laboratory. Using just a 200  $\mu$ L sample and minimal preparation, a full slide of samples can be acquired in less than 15 minutes
- ◆ Spot imaging during method development highlighted the differences in homogeneity of PFOS and PFOA spots suggesting that multiple internal standards would be necessary to pursue semi-quantitative analysis.
- ◆ Whole slide imaging can provide quick results for multiple samples across a range of PFAS concentrations in a single run with easy visual identification of potentially problematic samples.
- ◆ As a linear MALDI-TOF method, these results do not confirm the presence of PFAS but may be used as an indicative result to modify further testing protocols. Confirmatory LC-MS/MS analysis is required to determine the presence and concentration of any PFAS.

### References

1) Shimadzu Application News LCMS-151

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

All authors are current employees of Shimadzu Group.