


Comparative Analysis of Lipids in a Metabolic Dysfunction-Associated Steatohepatitis (MASH) Mouse Model using Imaging MS and LCMS

UTSouthwestern  
Medical Center  
SHIMADZU  
Excellence in Science

Doug Carlton<sup>1</sup>, Jeffrey Dahl<sup>1</sup>, Toshiya Matsubara<sup>1</sup>, Kyounghee Min<sup>2</sup>, Ei-ichi Matsuo<sup>3</sup>, Philipp E. Scherer<sup>2</sup>, Ruth Gordillo<sup>2</sup>


1. Shimadzu Scientific Instruments, Columbia, MD. 2. University of Texas Southwestern Medical Center, Dallas, TX. 3. Shimadzu Corporation, Kyoto, Japan

Preclinical mouse models of MASH are needed to characterize the disease and develop treatments in humans. Abnormal ratios of lipids such as phosphatidylcholine (PC) to phosphatidylethanolamine (PE) are associated with steatogenesis and inflammation [1]. To characterize the levels of several hundred potential lipid species of interest, we utilized a combined genetic and diet-induced MASH mouse model.



INTRODUCTION


- LCMS analysis was performed on the Shimadzu LCMS-9030 QTOF mass spectrometer (Shimadzu Corporation, Kyoto, Japan) on internal standard (SPLASH LipidoMIX™; Avanti Research, Alabaster, AL) fortified liver homogenate extracts. Samples were analyzed in both positive and negative mode. Binary gradient elution using a Phenomenex Kinetex (Phenomenex, Torrance, CA) column was performed. Targeted lipidomics data analysis was conducted with the Shimadzu Library for phospholipids profiling.
- For lipidomics imaging frozen liver sections (10-μm) were mounted on conductive glass slides. Matrix was deposited using an automated matrix vapor deposition system iMLayer™ (Shimadzu Corporation). 1,5-diaminonaphthalene (DAN) was used for negative mode, and α-cyano-4-hydroxycinnamic acid (CHCA) for positive mode. Atmospheric Pressure IMS analysis was carried out using the Shimadzu iMScope™ QT.




METHODS

MP 429


RESULTS





**A)** Unsupervised PCA components of lipids profiled in positive mode. **B)** Top 20 lipids identified in positive mode with relatively higher VIP score. **C)** Volcano plot TG over control (lipids profiled in negative mode) **D)** Phospholipids ratios associated with the diagnosis and persistence of MASH ( $\Phi \leq 0.001$ ,  $\text{¥} \leq 0.0001$ )[1-4]



Determination of liver lobule functional zones regions of interest (ROI) [5] was performed using IF-staining, overlay of reported lipid molecules localized in the functional zones [6], and computational determination of image clusters **A)** IF- staining: **Glutamine synthetase (pericentral)**, **E-cadherin (periportal)** and **DAPI B)** Overlay of false color molecular images **PI (38:3) (pericentral)**, **PI (36:3) (midzone)**, and **PE (38:6) (periportal)** **C)** Unsupervised image clustering (4 clusters) **pericentral** and **periportal** zones **D)** Optical image obtained with iMScope QT



Unsupervised PCA components of lipids (negative mode) in the different ROIs. Controls PC and PP zones are well differentiated. Control PC and TG PC are also well differentiated



Relative concentration differences **A)** Pericentral vs periportal zones in control livers **B)** Control pericentral vs TG pericentral (high abundance lipids) **C)** control pericentral vs TG pericentral (intermediate abundance lipids)

**References.** [1] [10.1016/j.cmet.2006.03.007](https://doi.org/10.1016/j.cmet.2006.03.007) [2] <https://doi.org/10.1016/j.bbamem.2017.04.006> [3] <https://doi.org/10.3897/folmed.64.e59297> [4] <https://doi.org/10.1016/j.dld.2024.05.015> [5] <https://doi.org/10.1073/pnas.1804203115> [6] <https://doi.org/10.1016/j.jhepr.2023.100725>  
**Acknowledgements.** Shimadzu Corporation (JP) and Shimadzu Scientific Instruments (USA) for the collaborative efforts in mass spectrometry technology resources and useful discussion. Maria G. Di-Lago for her technical assistance in the elaboration of this project (University of Bari, Italy. Italian Ministry of University and Research CUP code: H9I123000410007)