

## 1. Overview

An accurate mass, single quadrupole liquid chromatography mass spectrometry (SQ-LCMS) method was developed to facilitate screening environmentally relevant Per- and polyfluoroalkyl substances (PFAS)

## 2. Introduction

PFAS are ubiquitous and enduring contaminants that have harmful impacts on the environment and public health. Analysis of non-volatile PFAS is commonly performed by triple quadrupole liquid chromatography tandem mass spectrometry (TQ-LCMS); the ability to screen PFAS using SQ-LCMS provides a novel and economical alternative for these laboratories. SQ-LCMS is less selective compared to TQ-LCMS; therefore, accurate mass determination through line-shape calibration is a powerful technique for improving compound confirmation.


## 3. Method

The 40 PFAS targets and 31 internal standards were diluted to match EPA Method 1633A (EPAM1633A) analysis<sup>1</sup>. Neat standards were analyzed from 0.2-1560 ppb in-vial. Cerno Bioscience MassWorks was used to improve mass accuracy. Multiple calibrations were used with differing criteria to evaluate efficacy using the new Target Analysis function.


**Table 1.** HPLC and MS acquisition parameters

Nexera XR		LCMS-2050	
Analytical Column:	Shim-pack Scepter C18-120, 3 μm; 2.1×50 mm	Ionization Type:	ESI
Delay Column:	Shim-pack Scepter C18-120, 3 μm; 2.1×100 mm	Desolvation Temp:	250 °C
Mobile Phase	A: 2 mM Ammonium acetate in water B: Acetonitrile	DL Temp:	200 °C
Flow rate:	0.4 mL/min	Nebulizing Gas:	2.0 L/min
Oven temp (°C):	40 °C	Heating Gas:	7.0 L/min
Injection volume:	5 μL	Probe Position:	+2 mm
		Interface Voltage:	-0.5 kV


## 4. Results and Discussion



**Figure 1.** Chromatogram for all analytes at the mid-range calibration standard (1.25-62.4 ppb).



**Figure 2.** Accuracy (70-130%), S/N > 10, %RSD < 20, and %RSE < 20 were used to determine LOQ. IDL's were calculated from 11 replicates similar to MDLs in CFR Appendix B to Part 136<sup>2</sup>. The screening method LOQ, IDL, and EPA LOQ for the 40 PFAS targets are shown for a) carboxylic acids; b) sulfonates, sulfonamides, sulfonamidoethanols, and sulfonamidoacetic acids; and, c) fluorinated telomer acids and sulfonates.



**Figure 3.** Mass accuracy distribution for the highest calibration standard (62.5-1560 ppb). Compounds with low S/N and mass accuracy ≤ 50 ppm were not identified. The differing criteria and compounds used for analysis were as follows: **Condition 1** – Sigma ESI Tune mix; **Condition 2** – Condition 1 & benzenesulfonamide; **Condition 3** - 7 internal standards with mass accuracy ≤ 20 mDa ; **Condition 4** – Target PFAS with sufficient S/N; **Condition 5** – Target PFAS with mass accuracy ≤ 20 mDa and spectral accuracy ≥ 90%; and, **Condition 6** – Target PFAS with mass accuracy ≤ 10 mDa and spectral accuracy ≥ 95%.

## 5. Conclusion

- Effective screen for regulated PFAS for SQ-LCMS analysis with low LOQs comparable to EPA.
- Utilizing a calibration with more compounds (Condition 4) resulted in the best overall mass accuracies.
- Target Analysis in MassWorks enabled accurate mass determination for most PFAS by SQ-LCMS.

## 6. Reference

1. U.S. Environmental Protection Agency, Method 1633A, Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous, Solid, Biosolids, and Tissue Samples by LC-MS/MS. January 2024.
2. CFR Appendix B to Part 136, Title 40 -- Definition and Procedure for the Determination of the Method Detection Limit—Revision 2