

METHOD 552.3 DETERMINATION OF HALOACETIC ACIDS AND DALAPON IN DRINKING WATER BY LIQUID-LIQUID MICROEXTRACTION, DERIVATIZATION, AND GAS CHROMATOGRAPHY WITH ELECTRON CAPTURE DETECTION

EPA 815-B-03-002

Revision 1.0

July 2003

M. M. Domino and B.V. Pepich (Shaw Environmental and Infrastructure, Inc.)

D.J. Munch and P.S. Fair (US EPA, Office of Ground Water and Drinking Water)

Y. Xie (Penn State University)

D.J. Munch, J.W. Munch (US EPA, Office of Ground Water and Drinking Water) and A. M. Pawlecki-Vonderheide (ICI) Method 552.2, Revision 1.0 (1995)

J.W. Hodgeson (USEPA), D. Becker (Technology Applications, Inc.) Method 552.1, (1992)

J.W. Hodgeson (USEPA), J. Collins and R. E. Barth (Technology Applications, Inc.) Method 552.0, (1990)

**TECHNICAL SUPPORT CENTER
OFFICE OF GROUND WATER AND DRINKING WATER
U. S. ENVIRONMENTAL PROTECTION AGENCY
CINCINNATI, OHIO 45268**

METHOD 552.3

DETERMINATION OF HALOACETIC ACIDS AND DALAPON IN DRINKING WATER BY LIQUID-LIQUID MICROEXTRACTION, DERIVATIZATION, AND GAS CHROMATOGRAPHY WITH ELECTRON CAPTURE DETECTION

1. SCOPE AND APPLICATION

- 1.1 This is a gas chromatography (GC) method for the determination of haloacetic acids and dalapon in drinking waters. Accuracy, precision, and Detection Limit (DL) data are presented for the following compounds in reagent water and finished ground and surface waters.

<u>Analyte</u>	<u>Chemical Abstracts Service (CAS) Registry Number</u>
Bromochloroacetic acid (BCAA)	5589-96-8
Bromodichloroacetic acid (BDCAA)	71133-14-7
Chlorodibromoacetic acid (CDBAA)	5278-95-5
Dalapon	75-99-0
Dibromoacetic acid (DBAA)	631-64-1
Dichloroacetic acid (DCAA)	79-43-6
Monobromoacetic acid (MBAA)	79-08-3
Monochloroacetic acid (MCAA)	79-11-8
Tribromoacetic acid (TBAA)	75-96-7
Trichloroacetic acid (TCAA)	76-03-9

- 1.2 Detection Limits are compound, instrument, and matrix dependent. The DL is defined as the statistically calculated minimum amount that can be measured with 99% confidence that the reported value is greater than zero.⁽¹⁾ Detection Limits for the above listed analytes are provided in Section 17, Table 5. The DL differs from, and is lower than, the Minimum Reporting Level (MRL) (Sect. 3.17). The concentration range for target

analytes in this method was evaluated between 0.5 µg/L and 30 µg/L for a 40-mL sample. Precision and accuracy data are presented in Section 17, Tables 6-11.

- 1.3 This method is restricted to use by or under the supervision of analysts skilled in liquid-liquid extractions, derivatization procedures and the use of GC and interpretation of gas chromatograms.

2. SUMMARY OF METHOD

- 2.1 A 40-mL volume of sample is adjusted to a pH of 0.5 or less and extracted with 4 mL of either methyl tert-butyl ether (MTBE) or tert-amyl methyl ether (TAME) containing an internal standard. The haloacetic acids that have been partitioned into the organic phase are then converted to their methyl esters by the addition of acidic methanol followed by heating for 2 hours. The solvent phase containing the methylated haloacetic acids is separated from the acidic methanol by adding 7 mL of a concentrated aqueous solution of sodium sulfate. The aqueous phase is discarded. The extract is then neutralized with a saturated solution of sodium bicarbonate and the solvent layer is removed for analysis. The target analytes are identified and quantified by capillary column gas chromatography using an electron capture detector (GC/ECD). Analytes are quantified using procedural standard calibration.

3. DEFINITIONS

- 3.1 **EXTRACTION BATCH** – A set of up to 20 Field Samples (not including calibration standards and QC samples) extracted together by the same person(s) during a work day using the same lots of solvents, surrogate solution, and fortifying solutions. Required QC samples include a Laboratory Reagent Blank, Laboratory Fortified Matrix, either a Field Duplicate or Laboratory Fortified Matrix Duplicate and an appropriate number of Continuing Calibration Checks.
- 3.2 **ANALYSIS BATCH** – A set of samples that is analyzed on the same instrument during a 24-hour period that begins and ends with the analysis of the appropriate Continuing Calibration Check standards (CCC). Additional CCCs may be required depending on the length of the analysis batch and/or the number of Field Samples.

- 3.3 INTERNAL STANDARD (IS) – A pure analyte added to an extract or standard solution in a known amount and used to measure the relative responses of other method analytes and surrogates. The internal standard must be an analyte that is not a sample component.
- 3.4 SURROGATE ANALYTE (SUR) – A pure analyte, which chemically resembles target analytes and is extremely unlikely to be found in any sample. This analyte is added to a sample aliquot in a known amount before extraction or other processing, and is measured with the same procedures used to measure other sample components. The purpose of the SUR is to monitor method performance with each sample.
- 3.5 LABORATORY REAGENT BLANK (LRB) – An aliquot of reagent water or other blank matrix that is treated exactly as a sample including exposure to all glassware, equipment, solvents, reagents, sample preservatives, internal standards, and surrogates that are used in the extraction batch. The LRB is used to determine if method analytes or other interferences are present in the laboratory environment, the reagents, or the apparatus.
- 3.6 LABORATORY FORTIFIED BLANK (LFB) – An aliquot of reagent water or other blank matrix to which known quantities of the method analytes and all the preservation compounds are added. The LFB is analyzed exactly like a sample, and its purpose is to determine whether the methodology is in control, and whether the laboratory is capable of making accurate and precise measurements.
- 3.7 LABORATORY FORTIFIED SAMPLE MATRIX (LFSM) – An aliquot of a Field Sample to which known quantities of the method analytes and all the preservation compounds are added. The LFSM is processed and analyzed exactly like a sample, and its purpose is to determine whether the sample matrix contributes bias to the analytical results. The background concentrations of the analytes in the sample matrix must be determined in a separate aliquot and the measured values in the LFSM corrected for background concentrations.
- 3.8 LABORATORY FORTIFIED SAMPLE MATRIX DUPLICATE (LFSMD) – A second aliquot of the Field Sample used to prepare the LFSM which is fortified, extracted and analyzed identically to the LFSM. The LFSMD is used instead of the Field Duplicate to access method precision and accuracy when the occurrence of a target analyte is infrequent.

- 3.9 LABORATORY DUPLICATES (LD1 and LD2) – Two aliquots of the same sample split in the laboratory and analyzed separately with identical procedures. Analyses of LD1 and LD2 indicate precision associated with laboratory procedures, but not with sample collection, preservation, and storage procedures.
- 3.10 FIELD DUPLICATES (FD1 and FD2) – Two separate samples collected at the same time and place under identical circumstances, and treated exactly the same throughout field and laboratory procedures. Analyses of FD1 and FD2 give a measure of the precision associated with sample collection, preservation, and storage, as well as with laboratory procedures.
- 3.11 STOCK STANDARD SOLUTION (SSS) – A concentrated solution containing one or more method analytes prepared in the laboratory using assayed reference materials or purchased from a reputable commercial source.
- 3.12 PRIMARY DILUTION STANDARD (PDS) SOLUTION – A solution containing method analytes prepared in the laboratory from stock standard solutions and diluted as needed to prepare calibration solutions and other needed analyte solutions.
- 3.13 CALIBRATION STANDARD (CAL) – A solution prepared from the primary dilution standard solution and/or stock standard solutions and the internal standards and surrogate analytes. The CAL solutions are used to calibrate the instrument response with respect to analyte concentration.
- 3.14 CONTINUING CALIBRATION CHECK (CCC) – A calibration standard containing the method analytes, internal standard(s) and surrogate(s), which is analyzed periodically to verify the accuracy of the existing calibration for those analytes.
- 3.15 QUALITY CONTROL SAMPLE (QCS) – A sample prepared using a PDS of method analytes that is obtained from a source external to the laboratory and different from the source of calibration standards. The second source PDS and the surrogate PDS are used to fortify the QCS at a known concentration. The QCS is used to check calibration standard integrity.
- 3.16 DETECTION LIMIT – The minimum concentration of an analyte that can be identified, measured and reported with 99% confidence that the analyte concentration is greater than zero (Sect. 9.2.4). This is a statistical determination of precision. Accurate quantitation is not expected at this level.⁽¹⁾

- 3.17 **MINIMUM REPORTING LEVEL (MRL)** – The minimum concentration that can be reported as a quantified value for a target analyte in a sample following analysis. This defined concentration can be no lower than the concentration of the lowest continuing calibration standard for that analyte, and can only be used if acceptable quality control criteria for this standard are met.
- 3.18 **PROCEDURAL STANDARD CALIBRATION** – A calibration method in which aqueous calibration standards are prepared and processed (e.g., extracted, and/or derivatized) in exactly the same manner as the samples. All steps in the process from addition of sampling preservatives through instrumental analyses are included in the calibration. Using procedural standard calibration compensates for any inefficiencies in the processing procedure.
- 3.19 **MATERIAL SAFETY DATA SHEET (MSDS)** – Written information provided by vendors concerning a chemical’s toxicity, health hazards, physical properties, fire, and reactivity data including storage, spill, and handling precautions.

4. INTERFERENCES

- 4.1 All glassware must be meticulously cleaned. Wash glassware with detergent and tap water, rinse with tap water followed by reagent water. A final rinse with solvents may be needed. In place of a solvent rinse, non-volumetric glassware can be muffled at 400 °C for 2 hours. Volumetric glassware should not be heated in an oven above 120 °C. Store glassware inverted or capped with aluminum foil.
- 4.2 Method interferences may be caused by contaminants in solvents, reagents (including reagent water), sample bottles and caps, and other sample processing hardware that lead to discrete artifacts and/or elevated baselines in the chromatograms. All items used in the method must be routinely demonstrated to be free from interferences (less than $\frac{1}{3}$ the MRL for each target) under method conditions by processing and analyzing laboratory reagent blanks as described in Section 9. **Subtracting blank values from sample results is not permitted.**
- 4.2.1 Sodium sulfate can be a source of method interferences. After screening several brands, it was found that a grade suitable for pesticide residue analysis had the lowest amount of method interferants. If the suitability of the available sodium

sulfate is in question, extract and analyze a laboratory reagent blank (Sect. 3.5) to test for interferences prior to sample processing.

- 4.2.2 The ester of bromochloroacetic acid (BCAA) coelutes with a low-level interferant on both the primary and confirmation column. This interferant is present in the LRB and believed to be dimethyl sulfide that originates in the sodium sulfate. Because of difference in peak shape, the BCAA ester tends to “ride on” the interferant allowing accurate quantitation. Concentrations of BCAA at or below 2 μ g/L may require manual integration of the BCAA.
- 4.2.3 TAME, which has higher methylation efficiencies for the trihaloacetic acids than MTBE, is currently only available at 97% purity. Different lots of TAME should be carefully evaluated for potential interferants. Several lots of TAME contained an interferant that coeluted with the ester of monochloroacetic acid (MCAA) on the primary column.
- 4.3 The esters of the surrogate used in Method 552.2 (2,3-dibromopropanoic acid) and dichlorobromoacetic acid are not completely resolved under the chromatographic conditions for the confirmation column listed in Tables 3 and 4. The current surrogate, 2-bromobutanoic acid, is fully resolved from all method compounds on both columns under the conditions listed in Tables 1-4 (Figs. 1-4).
- 4.4 Matrix interferences may be caused by contaminants that are extracted from the sample. The extent of matrix interferences will vary considerably from source to source, depending upon the water sampled.
- 4.5 Interferences by phthalate esters can pose a major problem in analysis when using an electron capture detector (ECD). These compounds generally appear in the chromatogram as large peaks. Common flexible plastics contain varying amounts of phthalates that are easily extracted or leached during laboratory operations. Cross contamination of clean glassware routinely occurs when plastics are handled during extraction steps, especially when solvent-wetted surfaces are handled. Interferences from phthalates can best be minimized by avoiding the use of plastics in the laboratory. Exhaustive purification of reagents and glassware may be required to eliminate background phthalate contamination.^(2,3)

5. SAFETY

- 5.1 The toxicity or carcinogenicity of each reagent used in this method has not been precisely defined. Each chemical compound should be treated as a potential health hazard, and exposure to these chemicals should be minimized. The laboratory is responsible for maintaining a current awareness file of OSHA regulations regarding the safe handling of the chemicals specified in this method. A reference file of MSDSs should also be made available to all personnel involved in the chemical analysis. Additional references to laboratory safety are available.⁽⁴⁻⁶⁾
- 5.2 Pure standard materials and stock standards of these compounds should be handled with suitable protection to skin and eyes, and care should be taken not to breathe the vapors or ingest the materials.
- 5.3 The toxicities of the extraction solvents, MTBE or TAME, have not been well defined. Susceptible individuals may experience adverse effects upon skin contact or inhalation of vapors. Therefore protective clothing and gloves should be used and MTBE or TAME should be used only in a chemical fume hood or glove box. The same precaution applies to pure standard materials.

6. EQUIPMENT AND SUPPLIES

- 6.1 **SAMPLE CONTAINERS** – Amber glass bottles, at least 50-mL, fitted with PTFE (polytetrafluoroethylene) lined screw caps.
- 6.2 **EXTRACTION VIALS** – 60-mL clear glass vials with PTFE-lined screw caps.
- 6.3 **AUTOSAMPLER VIALS** – 2.0-mL amber vials with screw or crimp cap and a PTFE-faced seal.
- 6.4 **STANDARD SOLUTION STORAGE CONTAINERS** – 10- to 20-mL amber glass vials with PTFE-lined screw caps.
- 6.5 **GRADUATED CONICAL CENTRIFUGE TUBES WITH PTFE-LINED SCREW CAPS** – 15-mL with graduated 1-mL markings (Fisher Cat. #: 05-538-30A or equivalent).

- 6.6 PASTEUR PIPETTES – Glass, disposable.
- 6.7 PIPETTES – Class A, 2.0-, 3.0-, 4.0-, and 7.0-mL glass, or adjustable volume dispensers.
- 6.8 VOLUMETRIC FLASKS – Class A, suggested sizes include 5-mL, 10-mL, 100-mL.
- 6.9 MICRO SYRINGES – Various sizes.
- 6.10 HEATER (BLOCK, or SAND BATH or WATER BATH) – A block heater (Thermolyne Model DB16525, Fisher Cat. #: 11-716-50 or equivalent) capable of maintaining regulated, elevated temperature equipped with a sand bath or heating block (Fisher Cat. #: 11-716-27 or equivalent) capable of holding screw cap conical centrifuge tubes in Section 6.5 was used for method development. A thermostated water bath equipped with a surface layer of small plastic spheres (Fisher Cat. #: 14-220-31 or equivalent) to retard evaporation and subsequent heating of tube walls may also be used.
- 6.11 pH PAPER - With a pH range of 0 - 1.5 (Fisher Cat. #: 14-853-55 or equivalent).
- 6.12 BALANCE – Analytical, capable of weighing to the nearest 0.0001 g.
- 6.13 VORTEXER – Used to mix sample extracts (VWR Vortex-Genie, Cat. #:14216-184 or equivalent).
- 6.14 GAS CHROMATOGRAPH – Capillary GC, with a low volume (150 μ L) micro ECD (Agilent Model 6890 or equivalent). The injector system should not allow analytes to contact hot stainless steel or other metal surfaces that promote decomposition. The performance data in Section 17 was collected using hot, splitless injection using a 2-mm i.d. quartz liner. Other injection techniques such as temperature programmed injection, cold on-column injection, and large volume injection may be used if the QC criteria of Sections 9 and 10 are met. Alternate detectors which have equivalent or greater selectivity for the target compounds may be used.
- 6.15 AUTOSAMPLER – Agilent Model 7683 or equivalent.
- 6.16 PRIMARY GC COLUMN – DB-1701, 30-meter length, 0.25-mm i.d., 0.25- μ m film, fused silica capillary with chemically bonded (14% cyanopropylphenyl-methylpolysiloxane), or equivalent bonded, fused silica column.

- 6.17 CONFIRMATION GC COLUMN – DB-5.625, 30-meter length, 0.25-mm i.d., 0.25- μ m film, fused silica capillary with chemically bonded (“equivalent to” 5% phenyl-methylpolysiloxane), or equivalent bonded, fused silica column.

7. REAGENTS AND STANDARDS

- 7.1 REAGENTS AND SOLVENTS – Reagent grade or better chemicals should be used in all analyses. Unless otherwise indicated, it is intended that all reagents shall conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society, where such specifications are available. Other grades may be used, provided it is first determined that the reagent is of sufficiently high purity to permit its use without lessening the quality of the determination.

- 7.1.1 REAGENT WATER – Purified water which does not contain any measurable quantities of any target analytes or interfering compounds greater than 1/3 the MRL for each compound of interest.

- 7.1.2 METHYL tert-BUTYL ETHER (MTBE) – High purity, demonstrated to be free from analytes and interferences (HPLC grade or better).

- 7.1.3 Tert-AMYL METHYL ETHER (TAME) – High purity, demonstrated to be free from analytes and interferences.

7.1.3.1 Currently, TAME is available in > 97% purity. Several manufacturers’ products were evaluated and found to contain a compound that elutes with MCAA on the primary column. Lots evaluated from Fluka (Riedel-de Haen), stored over molecular sieves, were free from this interferant.

- 7.1.4 SODIUM SULFATE, Na_2SO_4 – Pesticide grade, granular, anhydrous. Interferences have been observed when lower quality grades have been used. If interferences are observed, it may be necessary to heat the sodium sulfate in a shallow tray at 400 °C for up to 4 hr to remove phthalates and other interfering organic substances. Store in a capped glass bottle rather than a plastic container.

- 7.1.5 SODIUM BICARBONATE, NaHCO_3 – ACS reagent grade.

- 7.1.6 AMMONIUM CHLORIDE, NH_4Cl – ACS reagent grade.

- 7.1.7 SULFURIC ACID – Concentrated, ACS reagent grade. **Substitution of hydrochloric acid (HCl) is not allowed. Solutions of HCl can contain trace levels of bromide, which can promote the formation of brominated HAAs.**
- 7.1.8 HELIUM (or HYDROGEN)– 99.999% pure or better, GC carrier gas.
- 7.1.9 ARGON 95%/METHANE 5% (or NITROGEN) – 99.999% pure or better, ECD make-up gas.
- 7.1.10 SODIUM SULFATE SOLUTION – Prepare an aqueous solution of sodium sulfate in reagent water at a concentration of 150 g/L. **Substitution of sodium chloride is not allowed. Sodium chloride can contain trace levels of bromide, which can promote the formation of brominated HAAs.**
- 7.1.11 SODIUM BICARBONATE SOLUTION, SATURATED – Add sodium bicarbonate to a volume of water, mixing periodically until the solution has a small amount of undissolved sodium bicarbonate that does not disappear upon further mixing.
- 7.1.12 10% SULFURIC ACID IN METHANOL SOLUTION – Add 10 mL of sulfuric acid dropwise (due to heat evolution) to 50-60 mL of methanol contained in a 100-mL volumetric flask that has been placed in a cooling bath. Mix, let cool, and dilute to the 100-mL mark with methanol.
- 7.2 STANDARD SOLUTIONS – Standard Solutions may be prepared from certified, commercially available solutions or from neat compounds. When preparing from neat material, compound purity needs to be 96% or greater. When a compound purity is assayed to be 96% or greater, the weight can be used without correction to calculate the concentration of the stock standard. Solution concentrations listed in this section were used to develop this method and are included as an example. Standards for sample fortification generally should be prepared in the smallest volume that can be accurately measured to minimize the addition of organic solvent to aqueous samples. Laboratories should use standard QC procedures to determine when Standard Solutions described in this section need to be replaced.
- 7.2.1 STANDARD PREPARATION AND STORAGE TECHNIQUES – All Standard Solutions are prepared from either neat or solid materials following the general procedure outlined below. This procedure assumes that the standard stock

solutions are prepared in 10-mL volumetric flasks. The procedure should be scaled accordingly if larger volumetric flasks are used. Most standards in this method are made up in MTBE, whether the extracting solvent is MTBE or TAME. This is because MTBE is more water-soluble than TAME.

7.2.1.1 For analytes which are solids in their pure form, prepare stock standard solutions by accurately weighing the suggested amount (refer to the sections below for suggested weights) of pure material into a clean, tared 10-mL volumetric flask. Dilute the flask to volume with MTBE and mix thoroughly.

7.2.1.2 For analytes which are liquid in their pure form at room temperature, place about 9.8 mL of MTBE into a 10-mL volumetric flask. Allow the flask to stand, unstoppered, for about 10 minutes to allow solvent film to evaporate from the inner walls of the volumetric, and weigh to the nearest 0.1 mg. Using a 10- μ L (or gastight) syringe add the desired volume of the neat standard material to the flask by keeping the syringe needle just above the surface of the MTBE. Be sure that the standard material falls dropwise directly into the MTBE without contacting the inner wall of the volumetric. Record the weight and calculate the concentration of the stock standard from the net gain in weight. Dilute to volume, stopper, then mix by inverting the flask several times.

7.2.1.3 Transfer the stock solutions to amber glass vials or amber bottles with PTFE-lined caps and store at ≤ 0 °C.

7.2.2 INTERNAL STANDARD (IS) SOLUTIONS – 1,2,3-trichloropropane (99+%) is used as an internal standard for the method. This compound has been shown to be an effective internal standard for the method analytes, but other compounds may be used if the QC requirements in Section 9 are met.

7.2.2.1 INTERNAL STANDARD STOCK SOLUTION (2.0 mg/mL) – Prepare an internal standard stock solution by accurately transferring approximately 0.0200 g of neat 1,2,3-trichloropropane (weighed to the nearest 0.1mg) into a 10-mL volumetric flask containing methanol as described in Section 7.2.1. The resulting concentration of the stock internal standard solution will be approximately 2.0 mg/mL.

7.2.2.2 INTERNAL STANDARD PRIMARY DILUTION STANDARD

(100µg/mL) – Prepare an internal standard fortification solution at 100µg/mL (or other suitable concentration) by the addition of 500 µL (or other appropriate volume) of the stock standard to a 10-mL volumetric flask containing MTBE. Dilute to volume, mix thoroughly and transfer to an amber glass vial with a PTFE-lined screw cap and store at $\leq 0^{\circ}\text{C}$.

7.2.2.3 MTBE (or TAME) EXTRACTION SOLVENT WITH INTERNAL

STANDARD (1.00 µg/mL) – The internal standard 1,2,3-trichloropropane is added to the extraction solvent prior to analyte extraction to compensate for any volumetric differences encountered during sample processing. Have enough working solvent available to extract all calibration standards, QC samples and field samples in each extraction batch. The volume of fortified solvent should be determined by the sample workload. Never use two different lots of working solvent for one extraction batch. The following example illustrates preparation of 100 mL of fortified solvent. Add 1 mL of the primary dilution standard (100 µg/mL) to a 100-mL volumetric flask containing MTBE (or TAME), dilute to volume and mix thoroughly. Transfer the standard to an amber bottle for storage. This results in a final internal standard concentration of 1.00 µg/mL. This solution is used to extract the samples (Sect. 11.1).

7.2.3 SURROGATE (SUR) ANALYTE STANDARD SOLUTION – 2-bromobutanoic acid (99+%) is used as a surrogate compound in this method to evaluate the extraction and derivatization procedures. This compound has been shown to be an effective surrogate for the method analytes and is well resolved from all target analytes and common interferants on both chromatographic columns. The Method 552.2 surrogate, 2,3-dibromopropanoic acid, can be used with this method, but requires a longer chromatographic run to separate it from DCBAA (Sect. 4.3). Other SUR analytes may be used as long as they are halogenated carboxylic acids and the QC requirements of Section 9.7 are met. All surrogates must be added as free acids. Alternate candidates should be confirmed to have adequate esterification efficiencies (e.g. $\geq 80\%$) by comparison to a commercially available premethylated standard. Lower esterification efficiencies may lead to poor SUR precision.

Note: 2-bromo-2-methylpropanoic acid was ruled out as a potential SUR because it was found to degrade in the Field and QC Samples if they were not processed

immediately. Nearly 50% of the fortified concentration was lost when samples were fortified and then allowed to set for 3 hours prior to acidification. Alternate surrogate compounds MUST be carefully evaluated.

7.2.3.1 SURROGATE STOCK SOLUTION (10 mg/mL) – Prepare surrogate stock standard solutions of 2-bromobutanoic acid [80-58-0] by accurately transferring approximately 0.100 g of the neat material (weighed to the nearest 0.1 mg.) into a clean, tared 10-mL volumetric flask as described in Section 7.2.1. The resulting concentration of the Surrogate Stock Solution will be approximately 10 mg/mL.

7.2.3.2 SURROGATE PRIMARY DILUTION STANDARD (20 µg/mL) – Prepare a primary dilution standard (PDS) at a concentration of 20 µg/mL (or other suitable concentration) by adding 50 µL (or suitable volume) of the SUR stock standard to a 25-mL volumetric flask containing MTBE. Dilute the flask to volume, and mix thoroughly. Transfer the SUR PDS to an amber glass vial with a PTFE-lined screw cap and store at ≤ 0 °C. Addition of 20 µL of the primary dilution standard to the 40-mL aqueous sample results in a surrogate concentration of 10 ng/mL.

7.2.4 ANALYTE STANDARD SOLUTIONS – Obtain the analytes listed in the table in Section 1.1 as neat or solid free acid standards or as commercially prepared ampulized solutions from a reputable standard manufacturer. The use of premethylated standards is not allowed for the preparation of analyte standards. Prepare the Analyte Stock and Primary Dilution Standards as described below.

7.2.4.1 ANALYTE STOCK SOLUTION – Prepare separate stock standard solutions for each analyte of interest at a concentration of 1-5 mg/mL in MTBE. Method analytes are available as neat materials or ampulized solutions (> 99% purity) from a number of commercial suppliers.

7.2.4.1.1 For analytes which are solids in their pure form, prepare stock standard solutions by accurately weighing approximately 0.01 to 0.05 grams of pure material in a 10-mL volumetric flask using the technique described in Section 7.2.1.1.

7.2.4.1.2 For analytes which are liquid in their pure form at room temperature, add 10.0 μL of standard material to the flask using the technique described in Section 7.2.1.2.

7.2.4.2 ANALYTE PRIMARY DILUTION STANDARD (PDS)– Prepare the Analyte PDS solution by combining and diluting the Analyte Stock Solutions prepared above (Sect. 7.2.4.1) with MTBE. This solution will be used to prepare the Analyte Secondary Dilution Standard. As a guideline to the analyst, the concentration used in the primary dilution standard solution during method development was 100 $\mu\text{g}/\text{mL}$ for all ten analytes.

7.2.4.3 ANALYTE SECONDARY DILUTION STANDARD (SDS)– The Analyte SDS is used to fortify reagent water for calibration standards. Prepare at least one Analyte SDS by diluting the Analyte PDS with methanol. Two Analyte SDSs were used during method development. The first was prepared at a concentration of 5.00 $\mu\text{g}/\text{mL}$ (1/20 dilution of the Analyte PDS) and the second at a concentration of 16.7 $\mu\text{g}/\text{L}$ (1/6 dilution of the Analyte PDS). The SDS should be made daily. It should not be stored, since methanol will derivatize the acids that it contains. The lowest calibration standard concentration must be at or below the MRL of each analyte. The concentrations of the other standards should span the range of concentration expected in the Field Samples or the working range of the instrument.

7.2.5 CALIBRATION STANDARDS – At least 5 calibration standards are required to prepare the initial calibration curve (Sect. 10.2). Fortify an appropriate number of reagent water solutions with varying amounts of the Analyte PDS over the concentration range of interest. During method development, the Calibration Standards typically ranged from 1.0 $\mu\text{g}/\text{L}$ to 20 $\mu\text{g}/\text{L}$. The lowest standard must be at or below the MRL, which will depend upon instrument sensitivity. Because this method employs a procedural calibration technique, these standards must be treated like samples and require the addition of all preservation and other reagents. They are extracted by the procedure described in Section 11.

8. SAMPLE COLLECTION, PRESERVATION, AND STORAGE

8.1 SAMPLE BOTTLE PREPARATION

- 8.1.1 Grab samples must be collected in accordance with conventional sampling practices⁽⁷⁾ using amber glass containers with PTFE-lined screw caps and capacities of at least 50 mL.
- 8.1.2 Prior to shipment to the field, add crystalline or granular ammonium chloride (NH_4Cl) to the sample containers to produce a concentration of 100 mg/L in the Field Sample. For a typical 50-mL sample, this requires 5 mg of ammonium chloride.

Note: Enough ammonium chloride must be added to the sample to convert the free chlorine residual in the sample matrix to combined chlorine. Chloramines, formed by the reaction of hypochlorite with the ammonium ion, do not react further to produce additional haloacetic acids at significant concentrations and protect against microbiological degradation.⁽⁸⁾ This concentration of ammonium chloride was determined to convert 8 mg/L of free chlorine residual to combined chlorine.

8.2. SAMPLE COLLECTION

- 8.2.1 Fill sample bottles but take care not to flush out the ammonium chloride. Because the target analytes of this method are not volatile, it is not necessary to ensure that the sample bottles are completely headspace free.
- 8.2.2 When sampling from a water tap, remove the aerator. Open the tap and allow the system to flush until the water temperature has stabilized (usually about 3-5 minutes). Collect samples from the flowing system.
- 8.2.3 When sampling from an open body of water, fill a 1-quart wide-mouth bottle or 1 liter beaker with sample from a representative area, and carefully fill sample vials from the container.
- 8.2.4 After collecting the sample, seal the bottle and agitate by hand for 15 seconds.

- 8.3 **SAMPLE SHIPMENT AND STORAGE** – All samples must be chilled during shipment and must not exceed 10 °C during the first 48 hours after collection. Samples must be confirmed to be at or below 10 °C when they are received at the laboratory. Samples stored in the lab must be held at or below 6 °C and protected from light until extraction. Samples should not be frozen.
- 8.4 **SAMPLE AND EXTRACT HOLDING TIMES** – Chlorinated Field Samples that are preserved according to the method should not exhibit biological degradation of analytes during the allotted 14-day storage time. Samples must be extracted within 14 days. Unchlorinated Field Samples should be extracted as soon as possible to prevent biological degradation of analytes. Extracts must be stored at ≤ -10 °C or less and protected from light in glass vials with PTFE-lined caps. MTBE extracts must be analyzed within 21 days of extraction. TAME extracts must be analyzed within 28 days of extraction.
- 8.4.1 The brominated trihaloacetic acids can exhibit degradation during storage. Tribromoacetic acid (TBAA), the least stable HAA ester, is degraded to bromoform. This is thought to occur as a result of peroxide contamination in the solvent. Low concentrations of peroxides that were barely detectable with the commercially available peroxide test strips were found to generate significant levels of bromoform in standard solutions. Bromoform can be chromatographed under method conditions (see Tables 1-4), and can be monitored as an indication of high peroxide levels in the solvent. If large bromoform peaks are observed in HAA-fortified reagent water samples (or calibration standards), the ether solvent should be purified or replaced. Storing these solvents under nitrogen minimizes peroxide formation.
- 8.4.2 Extracts that were prepared using the TAME procedure showed much lower tendency to form bromoform during the method development studies. This is thought to be associated with lower levels of peroxides. This may be due to TAME being stored over molecular sieves (Sect. 7.1.3).

9. QUALITY CONTROL

- 9.1 Each laboratory using this method is required to operate a formal Quality Control (QC) program. The requirements of this program consist of an Initial Demonstration of Capability (IDC) , and subsequent analysis in each analysis batch of a Laboratory

Reagent Blank (LRB), Continuing Calibration Check Standards (CCCs), a Laboratory Fortified Sample Matrix (LFSM), and either a Laboratory Fortified Sample Matrix Duplicate (LFSMD) or a Field Duplicate Sample. This section details the specific requirements for each QC parameter. The QC criteria discussed in the following sections are summarized in Section 17, Tables 15 and 16. These criteria are considered the minimum acceptable QC criteria, and laboratories are encouraged to institute additional QC practices to meet their specific needs.

9.1.1 Process all quality control samples through all steps of Section 11, including methylation. Sample preservatives as described in Section 8.1 must be added prior to extracting and analyzing the quality control samples.

9.2 INITIAL DEMONSTRATION OF CAPABILITY (IDC) – Requirements for the Initial Demonstration of Capability are described in the following sections and summarized in Section 17, Table 15.

9.2.1 INITIAL DEMONSTRATION OF LOW SYSTEM BACKGROUND – Before any samples are analyzed, it must be demonstrated that a laboratory reagent blank (LRB) is reasonably free of contamination.

9.2.2 INITIAL DEMONSTRATION OF ACCURACY – Prior to the analysis of the IDC samples, verify calibration accuracy with the preparation and analysis of a mid-level QCS as defined in Section 9.10. If the analyte recovery is not $\pm 30\%$ of the expected value, the accuracy of the method is unacceptable. The source of the problem must be identified and corrected. After the accuracy of the calibration has been verified, prepare and analyze 4-7 replicate LFBs (or CCCs in this method) fortified at 10 $\mu\text{g/L}$, or near the mid-range of the initial calibration curve, according to the procedure described in Section 11. Sample preservatives as described in Section 8.1.2 must also be added to these samples. The average recovery of the replicate values must be within $\pm 20\%$ of the expected value.

9.2.3 INITIAL DEMONSTRATION OF PRECISION – Using the same set of replicate data generated for Section 9.2.2, calculate the standard deviation and relative standard deviation of the replicate recoveries. The relative standard deviation (RSD) of the results of the replicate analyses must be less than 20%.

9.2.4 DETECTION LIMIT (DL) DETERMINATION – Prepare, extract and analyze at least seven replicate LFBs at a concentration estimated to be near the Detection

Limit over a period of at least three days (both extraction and analysis should be conducted over at least three days) using the procedure described in Section 11. The fortification level may be estimated by selecting a concentration with a signal of 2 to 5 times the noise level. The appropriate concentration will be dependent upon the sensitivity of the GC/ECD system being used. Sample preservatives as described in Section 8.1.2 must be added to these samples. Calculate the Detection Limit using the equation

$$DL = St_{(n-1, 1-\alpha=0.99)}$$

where

$t_{(n-1, 1-\alpha=0.99)}$ = Student's t value for the 99% confidence level with n-1 degrees of freedom (3.143 for 7 replicates)

n = number of replicates, and

S = standard deviation of the replicate analyses.

NOTE: Do not subtract blank values when performing DL calculations. The DL is a statistical determination of precision only.⁽¹⁾ If the DL replicates are fortified at a low enough concentration, it is likely that they will not meet precision and accuracy criteria, and may result in a calculated DL that is higher than the fortified concentration. Therefore no precision and accuracy criteria are specified for the DL.

- 9.2.5 Minimum Reporting Level (MRL) – The MRL is the threshold concentration of an analyte that a laboratory can expect to accurately quantitate in an unknown sample. The MRL should be established at an analyte concentration that is at least 3 times the DL or at a concentration which would yield a signal-to-noise (S/N) ratio of greater than or equal to five. Depending upon the study's data quality objectives it may be set at a higher concentration. **The concentration of the lowest calibration standard must be at or below the MRL.**
- 9.2.6 METHOD MODIFICATIONS – The analyst is permitted to modify GC columns, GC conditions (see Tables 1-4), and internal standards (see Sect. 7.2.2) or surrogate standards (see Sect. 7.2.3). Other detectors (see Sect. 6.14) may be used if they have equivalent or better selectivity and have sufficient sensitivity. Each time such method modifications are made, the analyst must repeat the procedures of the IDC (Sect. 9.2).

- 9.3 LABORATORY REAGENT BLANK (LRB) – An LRB is required with each extraction batch (Sect. 3.1) of samples to determine the background system contamination. If the LRB produces a peak within the retention time window of any analyte that would prevent the determination of that analyte, determine the source of contamination and eliminate the interference before processing samples. Background contamination must be reduced to an acceptable level before proceeding. Background from method analytes or contaminants that interfere with the measurement of method analytes must be below 1/3 of the MRL. If the target analytes are detected in the LRB at concentrations equal to or greater than this level, then all data for the problem analyte(s) must be considered invalid for all samples in the extraction batch.
- 9.4 CONTINUING CALIBRATION CHECK (CCC) – CCC Standards are prepared in the same extraction batch as the samples of interest. They must contain all compounds of interest, and they are extracted in the same manner as the Field Samples and calibration solutions used to prepare the initial calibration curve. Calibration checks, prepared with the samples being analyzed, are required at the beginning of each analysis batch, after every ten samples, and at the end of the analysis batch. See Section 10.3 for concentration requirements, frequency requirements, and acceptance criteria.
- 9.5 LABORATORY FORTIFIED BLANK (LFB) – Since this method utilizes procedural calibration standards, which are fortified reagent waters, there is no difference between the LFB and the continuing calibration check standard. Consequently, the analysis of an LFB is not required; however, the acronym LFB is used for clarity in the IDC.
- 9.6 INTERNAL STANDARDS (IS) – The analyst must monitor the peak area of each internal standard in all injections during each analysis day. The IS response (as indicated by peak area) for any chromatographic run must not deviate by more than $\pm 50\%$ from the average area measured during the initial calibration for that IS. A poor injection could cause the IS area to exceed these criteria. Inject a second aliquot of the suspect extract to determine whether the failure is due to poor injection or instrument response drift.
- 9.6.1 If the reinjected aliquot produces an acceptable internal standard response, report results for that aliquot.
- 9.6.2 If the internal standard area for the reinjected extract deviates greater than 50% from the initial calibration average, the analyst should check the continuing calibration check standards that ran before and after the sample. If the continuing calibration check fails the criteria of Section 10.3, recalibration is in order per

Section 10. If the calibration standard is acceptable, extraction of the sample should be repeated provided the sample is still within holding time. Otherwise, report results obtained from the reinjected extract, but annotate as suspect.

- 9.7 SURROGATE RECOVERY – The surrogate standard is fortified into the aqueous portion of all samples, LRBs, CCCs, LFSMs, and LFSMDs prior to extraction. It is also added to the calibration standards. The surrogate is a means of assessing method performance from extraction to final chromatographic measurement. Calculate the recovery (R) for the surrogate using the equation

$$R = \left(\frac{A}{B} \right) \times 100\%$$

where

A = calculated surrogate concentration for the QC or Field Sample, and

B = fortified concentration of the surrogate.

- 9.7.1 When surrogate recovery from a sample, blank, or CCC is less than 70% or greater than 130%, check (1) calculations to locate possible errors, (2) standard solutions for degradation, (3) contamination, and (4) instrument performance. Correct the problem and reanalyze the extract.
- 9.7.2 If the extract reanalysis meets the surrogate recovery criterion, report only data for the reanalyzed extract.
- 9.7.3 If the extract reanalysis fails the 70-130% recovery criterion, the analyst should check the calibration by injecting the last calibration standard that passed. If the calibration standard fails the criteria of Section 9.7.1, recalibration is in order per Section 10.2. If the calibration standard is acceptable, extraction of the sample should be repeated provided the sample is still within the holding time. If the re-extracted sample also fails the recovery criterion, report all data for that sample as suspect/surrogate recovery to inform the data user that the results are suspect due to surrogate recovery.
- 9.8 LABORATORY FORTIFIED SAMPLE MATRIX (LFSM) – Analysis of an LFSM (Sect. 3.7) is required in each extraction batch and is used to determine that the sample matrix does not adversely affect method accuracy. If the occurrence of target analytes in the samples is infrequent, or if historical trends are unavailable, a second LFSM or

LFSMD must be prepared, extracted, and analyzed from a duplicate field sample used to prepare the LFSM to assess method precision. Extraction batches that contain LFSMDs do not require the analysis of a Field Duplicate (Sect. 9.9). If a variety of different sample matrices are analyzed regularly, for example, drinking water from groundwater and surface water sources, method performance should be established for each. Over time, LFSM data should be documented for all routine sample sources for the laboratory.

9.8.1 Within each extraction batch, a minimum of one Field Sample is fortified as an LFSM for every 20 samples extracted. The LFSM is prepared by spiking a sample with an appropriate amount of Analyte PDS (Sect. 7.2.4.2). Select a spiking concentration that is greater than or equal to the matrix background concentration, if known. Use historical data and rotate through the designated concentrations when selecting a fortifying concentration.

9.8.2 Calculate the recovery (R) for each analyte using the equation

$$R = \frac{(A - B)}{C} \times 100\%$$

where

A = measured concentration in the fortified sample,
B = measured concentration in the unfortified sample, and
C = fortification concentration.

9.8.3 Analyte recoveries may exhibit matrix bias. For samples fortified at or above their native concentration, recoveries should range between 70 and 130%, except for low-level fortification near or at the MRL (within a factor of 2-times the MRL concentration) where 50 to 150% recoveries are acceptable. If the accuracy of any analyte falls outside the designated range, and the laboratory performance for that analyte is shown to be in control in the CCCs, the recovery is judged to be matrix biased. The result for that analyte in the unfortified sample is labeled suspect/matrix to inform the data user that the results are suspect due to matrix effects.

9.8.3.1 Because HAAs are disinfection by-products, many Field Samples will contain a number of HAAs at varying concentrations. Field Samples that have native HAA concentrations above the DL but below the MRL and are fortified at concentrations at or near the MRL should be corrected for

the native levels in order to obtain meaningful R values. This is the only permitted use of analyte results below the MRL.

9.9 FIELD DUPLICATE OR LABORATORY FORTIFIED SAMPLE MATRIX

DUPLICATE (FD or LFSMD) – Within each extraction batch, a minimum of one Field Duplicate (FD) or Laboratory Fortified Sample Matrix Duplicate (LFSMD) must be analyzed. Duplicates check the precision associated with sample collection, preservation, storage, and laboratory procedures. HAAs are typically found in waters disinfected with chlorine; however, if target analytes are not routinely observed in Field Samples, an LFSMD should be periodically analyzed rather than an FD.

9.9.1 Calculate the relative percent difference (RPD) for duplicate measurements (FD1 and FD2) using the equation

$$RPD = \frac{|FD1 - FD2|}{(FD1 + FD2)/2} \times 100 \%$$

9.9.1.1 RPDs for Field Duplicates should fall in the range of $\pm 30\%$. Greater variability may be observed when Field Duplicates have analyte concentrations that are within a factor of 2 of the MRL. At these concentrations Field Duplicates should have RPDs that fall in the range of $\pm 50\%$. If the RPD of any analyte falls outside the designated range, and the laboratory performance for that analyte is shown to be in control in the CCC, the recovery is judged to be matrix biased. The result for that analyte in the unfortified sample is labeled suspect/matrix to inform the data user that the results are suspect due to matrix effects.

9.9.2 If an LFSMD is analyzed instead of a Field Duplicate, calculate the relative percent difference (RPD) for duplicate LFSMs (LFSM and LFSMD) using the equation

$$RPD = \frac{|LFSM - LFSMD|}{(LFSM + LFSMD)/2} \times 100 \%$$

9.9.2.1 RPDs for duplicate LFSMs should fall in the range of $\pm 30\%$ for samples fortified at or above their native concentration. Greater variability may be observed when LFSMs are fortified at analyte concentrations that are within a factor of 2 of the MRL. LFSMs fortified at these concentrations

should have RPDs that fall in the range of $\pm 50\%$ for samples fortified at or above their native concentration. If the RPD of any analyte falls outside the designated range, and the laboratory performance for that analyte is shown to be in control in the CCC, the recovery is judged to be matrix biased. The result for that analyte in the unfortified sample is labeled suspect/matrix to inform the data user that the results are suspect due to matrix effects.

- 9.10 **QUALITY CONTROL SAMPLES (QCS)** – As part of the IDC (Sect. 9.2), each time a new Analyte PDS (Sect. 7.2.4.2) is prepared, and at least quarterly, analyze a QCS sample from a source different from the source of the calibration standards. If a second vendor is not available then a different lot of the standard should be used. The QCS should be prepared and analyzed just like a CCC. Acceptance criteria for the QCS is identical to the CCCs; the calculated amount for each analyte must be $\pm 30\%$ of the expected value. If measured analyte concentrations are not of acceptable accuracy, check the entire analytical procedure to locate and correct the problem.

10. CALIBRATION AND STANDARDIZATION

- 10.1 An acceptable initial calibration must be established during the IDC and prior to analyzing Field or QC Samples. After initial calibration is successful, a Continuing Calibration Check (CCC) is required at the beginning and the end of each analysis batch, and after every tenth sample (Sect. 10.3). Because this is a procedural standard method, the analyst will need to make a decision to include either an appropriate number of CCCs or an entire initial calibration curve with each extraction batch. Initial calibration must be repeated each time a major instrument modification is made or maintenance is performed. Failure to meet CCC criteria may also require recalibration.
- 10.2 **INITIAL CALIBRATION** – This method uses the procedural calibration technique to compensate for incomplete methylation of some of the target compounds. This is most pronounced for the brominated trihaloacetic acids. Many of the QC criteria throughout this method are expressed in terms of percent recovery. It should be noted that these recoveries are relative to the initial procedural curve rather than absolute recoveries.
- 10.2.1 Establish GC operating parameters equivalent to the suggested specifications in Section 17, Table 1 (or Table 2). The GC system must be calibrated using the

internal standard (IS) technique. Other columns or conditions may be used if equivalent or better performance can be demonstrated.

10.2.2 Prepare a set of at least 5 calibration standards as described in Section 7.2.5. The lowest concentration calibration standard must be at or below the MRL, which will depend on system sensitivity. It is recommended that at least four of the calibration standards be at concentrations greater than the MRL.

10.2.3 CALIBRATION – Use the GC data system software to generate a linear regression or quadratic calibration curve using the internal standard. The analyst may choose whether or not to force zero to obtain a curve that best fits the data. Examples of common GC system calibration curve options are: 1) A_x/A_{is} versus Q_x/Q_{is} ; and 2) RRF versus A_x/A_{is} . Establish a relative response factor by using the equation

$$RRF = \frac{(A_x)(Q_{is})}{(A_{is})(Q_x)}$$

where

A_x = integrated peak area of the analyte,

A_{is} = integrated peak area of the internal standard,

Q_x = quantity of analyte injected in ng or concentration units,

Q_{is} = quantity of internal standard injected in ng or concentration units, and

RRF = relative response factor.

10.2.4 Acceptance criteria for the calibration of each analyte is determined by calculating the concentration of each analyte and surrogate in each of the analyses used to generate the calibration curve. Each calibration point, except the lowest point, for each analyte should calculate to be 70-130% of its expected value. The lowest point should calculate to be 50-150% of its expected value. Laboratories that have difficulty achieving these criteria will have trouble meeting the QC requirements summarized in Section 10.3.2. These laboratories should reanalyze the calibration standards, restrict the range of calibration, or select an alternate method of calibration.

10.3 CONTINUING CALIBRATION CHECK (CCC) – An appropriate number of CCCs must be prepared with each extraction batch. The CCC verifies the initial calibration at the beginning and the end of each analysis batch, and after every 10th sample during analyses. In this context, a “sample” is considered to be a Field Sample. LFSMs, LFSMDs, FSDs, and CCCs are not counted as samples. The beginning CCC for each analysis batch must be at or below the MRL in order to verify instrument sensitivity prior to any analyses. If standards have been prepared such that all analytes are not in the same Calibration Standard, it may be necessary to analyze two Calibration Standard solutions to meet this requirement. Subsequent CCCs should alternate between a medium and a high concentration standard.

10.3.1 Inject an aliquot of the appropriate concentration calibration check standard solution prepared with the extraction batch and analyze using the same conditions used during the initial calibration.

10.3.2 Calculate the concentration of each analyte and surrogate in the check standard. The calculated amount for each analyte for medium and high level CCCs must be $\pm 30\%$ of the expected value. The calculated amount for the lowest calibration level for each analyte must be within $\pm 50\%$ of the expected value. If these conditions do not exist, then all data for the problem analyte must be considered invalid, and remedial action should be taken which may require recalibration. Any Field Sample extracts that have been analyzed since the last acceptable calibration verification should be reanalyzed after adequate calibration has been restored. For Extraction Batches with more than 10 Field Samples the analyst may either extract and analyze a third CCC or reanalyze the mid-level or high-level CCC as the final CCC.

11. PROCEDURE

11.1 SAMPLE EXTRACTION

11.1.1 Remove the samples from storage (Sect. 8.3) and allow them to equilibrate to room temperature.

11.1.2 Place 40 mL of the water sample into a precleaned 60-mL glass vial with a PTFE-lined screw cap using a clean, graduated cylinder for each sample.

11.1.3 Add 20 μL of surrogate standard (20 $\mu\text{g}/\text{mL}$ of 2-bromobutanoic acid in MTBE per Section 7.2.3.2) to the aqueous sample.

NOTE: When fortifying an aqueous sample with components that are contained in MTBE solution, be sure that the needle of the syringe is well below the surface of the water. After injection, cap the sample and invert once. This will insure that the standard solution is mixed well.

11.1.4 Adjust the pH to less than or equal to 0.5 by adding up to 2 mL of concentrated sulfuric acid. Cap, mix and then check the pH with narrow range pH paper (Sect. 6.11). **Substitution of other acids is not allowed.**

11.1.5 Add approximately 18 g of muffled sodium sulfate (Sect. 7.1.4) and immediately shake until almost all is dissolved. **Substitution of other salts for sodium sulfate is not allowed.** Sodium sulfate is added to increase the ionic strength of the aqueous phase and thus further drive the haloacetic acids into the organic phase. The addition of salt also decreases the solubility of MTBE or TAME in the aqueous phase and allows greater volumetric recovery of the extraction solvent. The addition of this salt should be done immediately after the addition of the sulfuric acid so that the heat generated from the addition of the acid (Sect. 11.1.4) will help dissolve the salt.

11.1.6 Add exactly 4.0 mL MTBE with IS or TAME with IS (Sect. 7.2.2.3) and shake vigorously for three minutes by hand. This may be accomplished for a entire Extraction Batch using a test tube rack.

11.1.7 Allow the phases to separate for approximately 5 minutes.

11.2 SAMPLE METHYLATION WITH ACIDIC METHANOL

11.2.1 Using a pasteur pipet, transfer 3 mL of the upper MTBE or TAME layer to a 15-mL graduated conical centrifuge tube.

11.2.2 Add 3 mL of 10% sulfuric acid in methanol (Sect. 7.1.12) to each centrifuge tube. Cap the tube.

11.2.3 Methylation of the method analytes is accomplished during this step. Careful control of both reaction time and reaction temperature are critical to method

precision and accuracy.⁽⁹⁾ Place the tubes in a heating block (or sand bath) at 50 ± 2 °C (for MTBE) or 60 ± 2 °C (for TAME) and heat for 2 hours (± 10 min). The tubes should fit snugly into the heating block to ensure proper heat transfer. Verify the reaction temperature by placing a thermometer into a tube containing water rather than inserting it into the block to ensure an accurate reading. Placing the thermometer directly in the heating block well or sand bath will give a higher temperature reading than the actual sample temperature. The MTBE reaction temperature is set at its highest practical limit (50 °C), since MTBE boils at 55 °C. Similarly, the TAME reaction temperature has been set at 60 °C, since methanol boils at 65 °C. This helps prevent solvent loss from the reaction tubes.

11.2.3.1 Even at 2 hours, methylation for some of the more sterically hindered compounds like TBAA, CDBAA and BDCAA is not complete. Shortening reaction time will decrease methylation efficiencies for compounds and result in lower precision.

11.2.3.2 Methylation efficiencies are increased with increasing temperature (and time). The TAME procedure has higher methylation efficiencies than MTBE procedure for the sterically hindered HAAs like the brominated trihaloacetic acids (Sect. 11.2.3.1) and is more precise for these targets. **Methylation temperature, however, should not be increased above the recommended reaction temperatures (Sect. 11.2.3).**

11.2.3.3 **Lower reaction temperatures or times are not allowed.** Care should be taken to ensure that the calibration standards are heated identically to Field Samples in the extraction batch.

11.2.3.4 Methylation may be accomplished by heating the reaction tubes with a water bath, provided the water bath is not covered in a way as to cause the entire tube to be heated. If tube walls are heated, the tubes can lose some of their contents leading to higher variability in analytical results. A water bath that is covered with a layer of small, floating plastic spheres may be used (Sect. 6.10).

11.2.4 Remove the centrifuge tubes from the heating source and allow them to cool before removing their caps.

- 11.2.5 Add 7 mL of a 150 g/L sodium sulfate solution (Sect. 7.1.10) to each centrifuge tube. Vortex each tube to ensure full equilibration between the phases. Allow the two phases to settle fully, but do not allow the tubes to sit more than a few minutes. Adding the sodium sulfate solution may cause some loss of the formed HAA-esters through acid-catalyzed hydrolysis over prolonged periods.
- 11.2.6 Remove and discard the lower (acidic aqueous methanol) phase from each tube with a long pasteur pipet. Leave no more than 0.3 mL of aqueous phase to ensure complete neutralization in the following step.
- 11.2.7 Add 1 mL of saturated sodium bicarbonate solution (Sect. 7.1.11). Vortex each centrifuge tube for several seconds at least four times to complete the neutralization reaction. Loosen the tube caps after the first agitation to release the evolved CO₂.
- 11.2.8 Transfer 1 mL of the upper ether layer to an auto-sampler vial. A duplicate vial should be filled using the excess extract.
- 11.2.9 Analyze the samples as soon as possible. Store the extracts at ≤ -10 °C (Sect. 8.4).

11.3 GAS CHROMATOGRAPHY

- 11.3.1 The instrument used in the development of this method was equipped with a low volume (150 μ L) micro electron capture detector (ECD). Other configurations are allowed as described in Section 6.14.
- 11.3.2 COLUMN SELECTION AND INSTALLATION – Strict attention should be paid to established column installation guidelines with regard to the proper cutting and placement of the capillary columns within the instrument. The brominated trihaloacetic acids are particularly sensitive to the condition of the injection port and GC column. If the response for these or other method analytes diminish, trimming approximately 0.5 m from the head of the column and replacing the GC inlet liner often restores the response for these analytes. Quartz inlet liners and inlet liners with standard deactivation were found to require less frequent maintenance than those with Siltek™ deactivation for this method. If conditions in the laboratory necessitate frequent column trimming, a guard column is recommended.

11.3.3 Method development was conducted with the GC set in constant pressure mode. If constant velocity mode is used by the analyst, linear velocities should be reduced from those listed in Tables 1-4 to obtain similar resolution.

11.4 ANALYSIS OF EXTRACTS

11.4.1 Establish operating conditions as described in Section 17, Table 1 or Table 2 (Table 3 or 4 if performing confirmation). Confirm that retention times, compound separation, and resolution are similar to those summarized in Tables 1-4 and Figures 1-4.

11.4.2 Establish an appropriate relative retention time window for each target and surrogate to identify them in the QC and Field Samples. This should be based on measurements of actual relative retention time variation for each compound in standard solutions analyzed on the GC over the course of time. Plus or minus three times the standard deviation of the relative retention time observed for each compound while establishing the initial calibration can be used to calculate a suggested window size; however, the experience of the analyst should weigh heavily on the determination of the appropriate retention window.

11.4.3 Calibrate the instrument as described in Section 10.2 or confirm the calibration is still valid by analyzing CCCs as described in Section 10.3. Begin analyzing Field and QC Samples at their appropriate frequency by injecting aliquots under the same conditions used to establish the initial calibration.

11.4.4 The analyst must not extrapolate beyond the established calibration range. If an analyte result exceeds the range of the initial calibration curve, the extract may be diluted with MTBE or TAME containing the internal standard (Sect. 7.2.1.3), and the diluted extract injected. Acceptable surrogate performance (Sect. 9.7) should be determined from the undiluted sample extract. Incorporate the dilution factor into final concentration calculations. The dilution will also affect analyte MRLs.

12. DATA ANALYSIS AND CALCULATION

12.1 Identify the method analytes in the sample chromatogram by comparing the retention time of the suspect peak to retention time of an analyte peak in a calibration standard.

Surrogate retention times should be confirmed to be within acceptance limits (Sect. 11.4.2) even if no target compounds are detected.

- 12.2 Calculate the analyte concentrations using the initial calibration curve generated as described in Section 10.2. Quantitate only those values that fall between the MRL and the highest calibration standard. Samples with target analyte responses that exceed the highest standard require dilution and reanalysis (Sect. 11.4.4).
- 12.3 Analyte identifications should be confirmed using the confirmation column specified in Table 3 (or Table 4) or another column that is dissimilar to the primary column. GC/MS confirmation is acceptable if the analyte concentrations are sufficient.
- 12.4 Adjust the calculated concentrations of the detected analytes to reflect the initial sample volume and any dilutions performed.
- 12.5 Analyte concentrations are reported in $\mu\text{g/L}$ as the total free acid (usually to 2 significant figures); however, calculations should use all available digits of precision.

13. METHOD PERFORMANCE

- 13.1 **PRECISION, ACCURACY, AND DETECTION LIMITS** – Tables for these data are presented in Section 17. Detection Limits are presented in Table 5 and were calculated using the formula presented in Section 9.2.4. Single laboratory precision and accuracy data are presented in Tables 6-11.
- 13.2 **SAMPLE STORAGE STABILITY STUDIES** – An analyte storage stability study was conducted by fortifying the analytes (10 $\mu\text{g/L}$ of each analyte) into a chlorinated surface water that was collected, preserved, and stored as described in Section 8. The average of triplicate analyses, conducted on days 0, 3, 7, and 14 are presented in Table 12.
- 13.3 **EXTRACT STORAGE STABILITY STUDIES** – Extract storage stability studies were conducted on TAME and MTBE extracts obtained from a chlorinated surface water fortified at 10 $\mu\text{g/L}$. The average of triplicate injections are reported in Tables 13 and 14. Extract storage stability can decline in the presence of peroxides. This is discussed in Section 8.4.

14. POLLUTION PREVENTION

- 14.1 This method utilizes liquid-liquid microextraction to extract analytes from water. It requires the use of very small volumes of organic solvent and very small quantities of pure analytes, thereby minimizing the potential hazards to both the analyst and the environment as compared to the use of large volumes of organic solvents in conventional liquid-liquid extractions.
- 14.2 For information about pollution prevention that may be applicable to laboratory operations, consult “Less is Better: Laboratory Chemical Management for Waste Reduction” available from the American Chemical Society’s Department of Government Relations and Science Policy, 1155 16th Street N.W., Washington, D.C., 20036.

15. WASTE MANAGEMENT

- 15.1 The analytical procedures described in this method generate relatively small amounts of waste since only small amounts of reagents and solvents are used. The matrices of concern are finished drinking water or source water. However, the Agency requires that laboratory waste management practices be conducted consistent with all applicable rules and regulations, and that laboratories protect the air, water, and land by minimizing and controlling all releases from fume hoods and bench operations. Also, compliance is required with any sewage discharge permits and regulations, particularly the hazardous waste identification rules and land disposal restrictions. For further information on waste management, consult “The Waste Management Manual for Laboratory Personnel” also available from the American Chemical Society at the address in Section 14.2.

16. REFERENCES

1. Glaser, J.A., Foerst, D.L., McKee, G.D., Quave, S.A., and Budde, W.L., “Trace Analyses for Wastewaters,” Environ. Sci. Technol., 15 (1981) 1426-1435.
2. Giam, C.S., Chan, H.S., and Nef, G.S., "Sensitive Method for Determination of Phthalate Ester Plasticizers in Open-Ocean Biota Samples," Analytical Chemistry, 47, 2225 (1975).

3. Giam, C.S., and Chan, H.S., "Control of Blanks in the Analysis of Phthalates in Air and Ocean Biota Samples," U.S. National Bureau of Standards, Special Publication 442, 701-708, (1976).
4. "OSHA Safety and Health Standards, General Industry," (29CFR1910). Occupational Safety and Health Administration, OSHA 2206, (Revised, Jan. 1976).
5. "Carcinogens-Working with Carcinogens," Publication No. 77-206, Department of Health, Education, and Welfare, Public Health Service, Center for Disease Control, National Institute of Occupational Safety and Health, Atlanta, Georgia, August 1977.
6. "Safety In Academic Chemistry Laboratories," 3rd Edition, American Chemical Society Publication, Committee on Chemical Safety, Washington, D.C., (1979).
7. ASTM Annual Book of Standards, Part II, Volume 11.01, D3370-82, "Standard Practice for Sampling Water," American Society for Testing and Materials, Philadelphia, PA, 1986.
8. Pepich, B.V., Domino, M.M., Munch, D.J. and Fair, P.S. "Validating Sample Preservation Techniques and Holding Times for the Approved Compliance Monitoring Methods for Haloacetic Acids Under the U.S. EPA's Stage 1 D/DBP Rule," Submitted to Wat. Res., Sept. 2002.
9. Xie, Y., Rashid, I., Zhou, H., and Gammie, L. "Acidic Methanol Methylation for HAA Analysis: Limitations and Possible Solutions" JAWWA, 94:11, 115-122 (Nov. 2002).

17. TABLES, DIAGRAMS, FLOWCHARTS, AND VALIDATION DATA

TABLE 1: CHROMATOGRAPHIC CONDITIONS AND AVERAGE RETENTION TIME DATA FOR THE PRIMARY COLUMN (DB-1701) USING MTBE

Compound	Average T _r (min)*	RSD (%)
Monochloroacetic acid (MCAA)	9.29	0.03
Monobromoacetic acid (MBAA)	14.19	0.02
Dalapon	14.69	0.02
Dichloroacetic acid (DCAA)	15.07	0.02
Bromoform**	17.12	0.01
Trichloroacetic acid (TCAA)	18.89	0.01
1,2,3 Trichloropropane (IS)	20.69	0.01
Bromochloroacetic acid (BCAA)	21.22	0.01
2-bromobutanoic acid (SUR)	21.63	0.01
Bromodichloroacetic acid (BDCAA)	23.60	0.01
Dibromoacetic acid (DBAA)	24.01	0.005
Chlorodibromoacetic acid (CDBAA)	25.55	0.004
Tribromoacetic acid (TBAA)	26.88	0.003

*The average retention time represents the average of 8 injections of fortified reagent water extracts that had been fortified at 10 µg/L of each analyte.

**Bromoform is not a target analyte. RT value is provided for information only.

Primary Column: DB-1701, 30 m x 0.25 mm i.d., 0.25 µm film thickness.

Injector: Injector temperature, 210 °C; 2-mm straight quartz liner; injection volume 1 µL; splitless injection hold for 45 sec then purge @ 30 mL/min.

GC Program: (for MTBE) 40 °C initial held for 10 minutes, program at 2.5 °C/min to 65 °C, then 10 °C/min to 85 °C, then 20 °C/min to 205 °C, hold for 7 min.

Detector: Agilent Micro ECD (150 µL volume); detector temperature, 290 °C; detector make up gas 95% Ar/5% CH₄ at 20 mL/min.

Carrier Gas: Helium (UHP), set at constant pressure. Initial carrier gas velocity (at 40 °C) 33 cm/sec.

Data Collection: Agilent GC Chemstation with a digitization rate of 20 Hz.

TABLE 2: CHROMATOGRAPHIC CONDITIONS AND AVERAGE RETENTION TIME DATA FOR THE PRIMARY COLUMN (DB-1701) USING TAME

Compound	Average T _r (min)*	RSD (%)
Monochloroacetic acid (MCAA)	5.94	0.08
Monobromoacetic acid (MBAA)	8.76	0.04
Dalapon	9.17	0.04
Dichloroacetic acid (DCAA)	9.34	0.04
Bromoform**	11.06	0.02
Trichloroacetic acid (TCAA)	12.28	0.02
1,2,3 Trichloropropane (IS)	13.56	0.01
Bromochloroacetic acid (BCAA)	13.87	0.01
2-bromobutanoic acid (SUR)	14.20	0.01
Bromodichloroacetic acid (BDCAA)	15.78	0.004
Dibromoacetic acid (DBAA)	16.14	0.003
Chlorodibromoacetic acid (CDBAA)	17.58	0.002
Tribromoacetic acid (TBAA)	18.89	0.02

*The average retention time represents the average of 8 injections of fortified reagent water extracts that had been fortified at 10 µg/L of each analyte.

**Bromoform is not a target analyte. RT value is provided for information only.

Primary Column: DB-1701, 30 m x 0.25 mm i.d., 0.25 µm film thickness.

Injector: Injector temperature, 210 °C; 2-mm straight quartz liner; injection volume 1 µL; splitless injection hold for 45 sec then purge @ 30 mL/min.

GC Program: (for TAME) 55 °C initial hold 8 minutes, program at 2.5 °C/min to 65 °C, then 10 °C/min to 85 °C, then 20 °C/min to 205 °C, hold for 7 min.

Detector: Agilent Micro ECD (150 µL volume); detector temperature, 290 °C; detector make up gas 95% Ar/5% CH₄ at 20 mL/min.

Carrier Gas: Helium (UHP), set at constant pressure. Initial carrier gas velocity (at 55 °C) 32 cm/sec.

Data Collection: Agilent GC Chemstation with a digitization rate of 20 Hz.

TABLE 3: CHROMATOGRAPHIC CONDITIONS AND AVERAGE RETENTION TIME DATA FOR THE CONFIRMATION COLUMN (DB-5.625) USING MTBE

Compound	Average T _r (min)*	RSD (%)
Monochloroacetic acid (MCAA)	5.28	0.04
Monobromoacetic acid (MBAA)	8.31	0.01
Dichloroacetic acid (DCAA)	9.05	0.01
Dalapon	10.63	0.02
Bromoform**	12.74	0.02
Trichloroacetic acid (TCAA)	14.73	0.02
Bromochloroacetic acid (BCAA)	15.25	0.02
1,2,3 Trichloropropane (IS)	15.65	0.01
2-bromobutanoic acid (SUR)	17.41	0.02
Dibromoacetic acid (DBAA)	21.08	0.01
Bromodichloroacetic acid (BDCAA)	21.39	0.01
Chlorodibromoacetic acid (CDBAA)	24.35	0.002
Tribromoacetic acid (TBAA)	25.94	0.000

*The average retention time represents the average of 8 injections of fortified reagent water extracts that had been fortified at 10 µg/L of each analyte.

**Bromoform is not a target analyte. RT value is provided for information only.

Confirmation

Column: DB-5.625, 30 m x 0.25 mm i.d., 0.25 µm film thickness.

Injector: Injector temperature, 210 °C; 2-mm straight quartz liner; injection volume 1 µL ; splitless injection hold for 45 sec then purge @ 30 mL/min.

GC Program: (for MTBE) 40 °C initial held for 10 minutes, program at 2.5 °C/min to 65 °C, then 10 °C/min to 85 °C, then 20 °C/min to 205 °C,. Post run 210 °C hold for 7 min.

Detector: Agilent Micro ECD (150 µL volume); detector temperature, 290 °C; detector make up gas 95% Ar/5% CH₄ at 20 mL/min.

Carrier Gas: Helium (UHP), set at constant pressure. Initial carrier gas velocity (at 40 °C) 32 cm/sec.

Data Collection: Agilent GC Chemstation with a digitization rate of 20 Hz.

TABLE 4: CHROMATOGRAPHIC CONDITIONS AND AVERAGE RETENTION TIME DATA FOR THE CONFIRMATION COLUMN (DB-5.625) USING TAME

Compound	Average T _r (min)*	RSD (%)
Monochloroacetic acid (MCAA)	4.07	0.26
Monobromoacetic acid (MBAA)	5.47	0.14
Dichloroacetic acid (DCAA)	5.81	0.12
Dalapon	6.61	0.10
Bromoform**	7.97	0.07
Trichloroacetic acid (TCAA)	9.11	0.06
Bromochloroacetic acid (BCAA)	9.44	0.05
1,2,3 Trichloropropane (IS)	9.79	0.05
2-bromobutanoic acid (SUR)	11.05	0.03
Dibromoacetic acid (DBAA)	13.77	0.01
Bromodichloroacetic acid (BDCAA)	14.01	0.01
Chlorodibromoacetic acid (CDBAA)	16.45	0.008
Tribromoacetic acid (TBAA)	17.96	0.004

*The average retention time represents the average of 8 injections of fortified reagent water extracts that had been fortified at 10 µg/L of each analyte.

**Bromoform is not a target analyte. RT value is provided for information only.

Confirmation

Column: DB-5.625, 30 m x 0.25 mm i.d., 0.25 µm film thickness.

Injector: Injector temperature, 210 °C; 2-mm straight quartz liner; injection volume 1 µL; splitless injection hold for 45 sec then purge @ 30 mL/min.

GC Program: (for TAME) 55 °C initial held for 8 minutes, program at 2.5 °C/min to 65 °C, then 10 °C/min to 85 °C, then 20 °C/min to 205 °C,. Post run 210 °C hold for 7 min.

Detector: Agilent Micro ECD (150 µL volume); detector temperature, 290 °C; detector make up gas 95% Ar/5% CH₄ at 20 mL/min.

Carrier Gas: Helium (UHP), set at constant pressure. Initial carrier gas velocity (at 55 °C) 32 cm/sec.

Data Collection: Agilent GC Chemstation with a digitization rate of 20 Hz.

TABLE 5: DETECTION LIMITS IN REAGENT WATER USING THE MTBE AND TAME PROCEDURES

Compound	MTBE		TAME	
	Fortification Level (µg/L)	Detection Limit ^a (µg/L)	Fortification Level (µg/L)	Detection Limit ^a (µg/L)
Monochloroacetic acid (MCAA)	1.00	0.17	1.00	0.20
Monobromoacetic acid (MBAA)	0.50	0.027	0.50	0.13
Dalapon	0.50	0.024	0.50	0.14
Dichloroacetic acid (DCAA)	0.50	0.020	0.50	0.084
Trichloroacetic acid (TCAA)	0.50	0.019	0.50	0.024
Bromochloroacetic acid (BCAA)	0.50	0.016	0.50	0.029
Bromodichloroacetic acid (BDCAA)	0.50	0.034	0.50	0.031
Dibromoacetic acid (DBAA)	0.50	0.012	0.50	0.021
Chlorodibromoacetic acid (CDBAA)	0.50	0.054	0.50	0.035
Tribromoacetic acid (TBAA)	0.50	0.11	0.50	0.097

^aFortified reagent waters were extracted and analyzed over 3 days for 7-9 replicates following the procedure outlined in Section 9. MTBE detection limits were determined with 9 replicates. TAME detection limits were determined with 7 replicates.

TABLE 6: PRECISION AND ACCURACY IN REAGENT WATER FORTIFIED AT 1.0 µg/L

Compound	MTBE			TAME		
	Mean Recovery (%)	RSD (%) (n=8)	S/N Ratio ^a	Mean Recovery (%)	RSD (%) (n=8)	S/N Ratio ^a
Monochloroacetic acid (MCAA)	95.8	4.0	12	81.4	5.1	17
Monobromoacetic acid (MBAA)	92.2	1.8	81	90.7	3.7	92
Dalapon	97.5	1.6	110	92.8	2.1	150
Dichloroacetic acid (DCAA)	93.8	1.6	190	97.8	2.2	190
Trichloroacetic acid (TCAA)	105	0.52	580	107	0.90	520
Bromochloroacetic acid (BCAA)	102	0.36	580	103	0.94	740
Bromodichloroacetic acid (BDCAA)	117	1.2	1400	113	1.1	1800
Dibromoacetic acid (DBAA)	105	0.63	1300	105	0.86	1800
Chlorodibromoacetic Acid (CDBAA)	125	2.3	650	112	1.5	1600
Tribromoacetic acid (TBAA)	128	4.1	970	109	1.8	1100

^aSignal-to-noise ratios were calculated for each target compound peak by dividing the peak height for each compound by the peak-to-peak noise, which was determined for each component from the method blank over a period of time equal to the full peak width in the target analyte's retention time window.

TABLE 7: PRECISION AND ACCURACY IN REAGENT WATER FORTIFIED AT 10 µg/L

Compound	MTBE		TAME	
	Mean Recovery ^a (%)	RSD (%) (n=8)	Mean Recovery ^a (%)	RSD (%) (n=8)
Monochloroacetic acid (MCAA)	101	3.5	102	2.5
Monobromoacetic acid (MBAA)	101	2.2	101	1.6
Dalapon	99.8	0.88	99.9	0.50
Dichloroacetic acid (DCAA)	98.2	0.73	100	0.33
Trichloroacetic acid (TCAA)	102	1.2	101	1.1
Bromochloroacetic acid (BCAA)	101	0.72	101	0.65
Bromodichloroacetic acid (BDCAA)	107	2.4	103	1.2
Dibromoacetic acid (DBAA)	102	0.52	101	0.59
Chlorodibromoacetic acid (CDBAA)	111	3.8	104	1.6
Tribromoacetic acid (TBAA)	113	4.7	104	2.0

TABLE 8: PRECISION AND ACCURACY IN CHLORINATED SURFACE WATER FORTIFIED AT 1.0 µg/L

Compound	MTBE		TAME	
	Mean Recovery ^a (%)	RSD (%) (n=8)	Mean Recovery ^b (%)	RSD (%) (n=8)
Monochloroacetic acid (MCAA)	112	6.2	131	6.2
Monobromoacetic acid (MBAA)	99.5	2.3	98.6	4.2
Dalapon	107	2.2	103	3.5
Dichloroacetic acid (DCAA)	103	1.5	106	3.8
Trichloroacetic acid (TCAA)	89.2	1.4	89.0	1.1
Bromochloroacetic acid (BCAA)	99.5	1.6	102	2.2
Bromodichloroacetic acid (BDCAA)	91.7	4.0	87.5	1.8
Dibromoacetic acid (DBAA)	103	2.8	103	4.5
Chlorodibromoacetic acid (CDBAA)	95.5	7.1	94.4	3.3
Tribromoacetic acid (TBAA)	110	8.1	111	4.3

^aRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.14 µg/L), DCAA (0.56 µg/L), TCAA (0.32 µg/L), BCAA (1.1 µg/L), BDCAA (0.68 µg/L), DBAA (2.2 µg/L), CDBAA (0.91 µg/L), and TBAA (1.0 µg/L).

^bRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.22 µg/L), DCAA (0.78 µg/L), TCAA (0.36 µg/L), BCAA (2.0 µg/L), BDCAA (0.62 µg/L), DBAA (4.0 µg/L), CDBAA (0.82 µg/L), and TBAA (0.79 µg/L).

TABLE 9: PRECISION AND ACCURACY IN CHLORINATED SURFACE WATER FORTIFIED AT 10 µg/L

Compound	MTBE		TAME	
	Mean Recovery ^a (%)	RSD (%) (n=8)	Mean Recovery ^b (%)	RSD (%) (n=8)
Monochloroacetic acid (MCAA)	97.8	3.1	97.1	3.8
Monobromoacetic acid (MBAA)	98.8	2.0	98.8	2.0
Dalapon	99.6	0.90	97.7	0.42
Dichloroacetic acid (DCAA)	98.0	1.1	97.6	0.42
Trichloroacetic acid (TCAA)	99.8	0.80	100	1.4
Bromochloroacetic acid (BCAA)	101	1.4	101	0.93
Bromodichloroacetic acid (BDCAA)	103	4.3	104	2.3
Dibromoacetic acid (DBAA)	102	1.4	101	0.93
Chlorodibromoacetic acid (CDBAA)	105	5.5	105	3.2
Tribromoacetic acid (TBAA)	108	6.3	106	3.4

^aRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.14 µg/L), DCAA (0.56 µg/L), TCAA (0.32 µg/L), BCAA (1.1 µg/L), BDCAA (0.68 µg/L), DBAA (2.2 µg/L), CDBAA (0.91 µg/L), and TBAA (1.0 µg/L).

^bRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.22 µg/L), DCAA (0.78 µg/L), TCAA (0.36 µg/L), BCAA (2.0 µg/L), BDCAA (0.62 µg/L), DBAA (4.0 µg/L), CDBAA (0.82 µg/L), and TBAA (0.79 µg/L).

TABLE 10: PRECISION AND ACCURACY IN CHLORINATED GROUND WATER FORTIFIED AT 1.0 µg/L

Compound	MTBE		TAME	
	Mean Recovery ^a (%)	RSD (%) (n=8)	Mean Recovery ^b (%)	RSD (%) (n=8)
Monochloroacetic acid (MCAA)	126	4.5	124	9.5
Monobromoacetic acid (MBAA)	113	3.8	98.9	2.6
Dalapon	112	1.1	87.3	2.7
Dichloroacetic acid (DCAA)	102	3.7	107	3.2
Trichloroacetic acid (TCAA)	92.2	1.2	103	2.1
Bromochloroacetic acid (BCAA)	106	3.8	104	2.4
Bromodichloroacetic acid (BDCAA)	105	6.1	106	4.0
Dibromoacetic acid (DBAA)	111	5.3	105	1.7
Chlorodibromoacetic acid (CDBAA)	103	8.8	103	4.5
Tribromoacetic acid (TBAA)	99.2	7.3	103	3.2

^aRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.07 µg/L), Dalapon (0.19 µg/L), DCAA (2.6 µg/L), TCAA (1.1 µg/L), BCAA (3.0 µg/L), BDCAA (1.8 µg/L), DBAA (3.3 µg/L), CDBAA (1.8 µg/L), and TBAA (1.1 µg/L).

^bRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.28 µg/L), Dalapon (0.35 µg/L) DCAA (2.7 µg/L), TCAA (1.6 µg/L), BCAA (3.2 µg/L), BDCAA (2.4 µg/L), DBAA (2.5 µg/L), CDBAA (1.7 µg/L), and TBAA (0.54 µg/L).

TABLE 11: PRECISION AND ACCURACY IN CHLORINATED GROUND WATER FORTIFIED AT 10 µg/L

Compound	MTBE		TAME	
	Mean Recovery ^a (%)	RSD (%) (n=8)	Mean Recovery ^b (%)	RSD (%) (n=8)
Monochloroacetic acid (MCAA)	104	1.8	99.7	1.7
Monobromoacetic acid (MBAA)	100	1.4	98.7	1.1
Dalapon	97.8	1.1	97.9	0.63
Dichloroacetic acid (DCAA)	95.9	1.0	97.1	0.83
Trichloroacetic acid (TCAA)	98.2	0.81	101	1.4
Bromochloroacetic acid (BCAA)	100	0.53	99.8	1.4
Bromodichloroacetic acid (BDCAA)	106	1.8	106	2.3
Dibromoacetic acid (DBAA)	101	0.89	101	1.6
Chlorodibromoacetic acid (CDBAA)	110	2.8	106	2.8
Tribromoacetic acid (TBAA)	116	3.4	105	2.6

^aRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.07 µg/L), Dalapon (0.19 µg/L), DCAA (2.6 µg/L), TCAA (1.1 µg/L), BCAA (3.0 µg/L), BDCAA (1.8 µg/L), DBAA (3.3 µg/L), CDBAA (1.8 µg/L), and TBAA (1.1 µg/L).

^bRecoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.28 µg/L), Dalapon (0.35 µg/L), DCAA (2.7 µg/L), TCAA (1.6 µg/L), BCAA (3.2 µg/L), BDCAA (2.4 µg/L), DBAA (2.5 µg/L), CDBAA (1.7 µg/L), and TBAA (0.54 µg/L).

TABLE 12: SAMPLE HOLDING TIME DATA FOR SAMPLES FROM A CHLORINATED SURFACE WATER FORTIFIED WITH METHOD ANALYTES^a

Compound	Day 0 R ^b (%)	Day 3 R ^b (%)	Day 7 R ^b (%)	Day 14 R ^b (%)
Monochloroacetic Acid (MCAA)	97.9	96.8	98.5	104
Monobromoacetic Acid (MBAA)	107	99.6	105	109
Dalapon	104	104	103	99.7
Dichloroacetic Acid (DCAA)	101	99.0	102	98.8
Trichloroacetic Acid (TCAA)	101	102	98.8	96.2
Bromochloroacetic Acid (BCAA)	101	99.3	100	97.8
Bromodichloroacetic Acid (BDCAA)	112	111	113	110
Dibromoacetic Acid (DBAA)	102	99.1	100	96.9
Chlorodibromoacetic Acid (CDBAA)	117	115	119	116
Tribromoacetic Acid (TBAA)	123	118	124	121

^a All samples were fortified at 10 µg/L. All samples were stored at 10 °C for 48 hours, and at 6 °C thereafter.

^b Recoveries were corrected for haloacetic acid (HAA) concentrations in the unfortified matrix, based on an average value of three measurements. HAAs detected in the unfortified matrix included MBAA (0.30 µg/L), DCAA (0.73 µg/L), TCAA (0.33 µg/L), BCAA (2.0 µg/L), BDCAA (0.67 µg/L), DBAA (4.3 µg/L), CDBAA (0.91 µg/L), and TBAA (1.1 µg/L).

TABLE 13: MTBE EXTRACT HOLDING TIME DATA FOR SAMPLES FROM A CHLORINATED SURFACE WATER FORTIFIED WITH METHOD ANALYTES^a

Compound	Initial R (%)	Day 7 R (%)	Day 14 R (%)	Day 21 R (%)
Monochloroacetic Acid (MCAA)	99.0	98.6	100	99.3
Monobromoacetic Acid (MBAA)	99.7	98.2	98.1	99.0
Dalapon	99.0	96.7	98.1	97.8
Dichloroacetic Acid (DCAA)	97.9	96.3	97.5	97.3
Trichloroacetic Acid (TCAA)	99.3	99.6	98.5	99.1
Bromochloroacetic Acid (BCAA)	100	101	100	101
Bromodichloroacetic Acid (BDCAA)	99.3	106	102	103
Dibromoacetic Acid (DBAA)	101	102	101	102
Chlorodibromoacetic Acid (CDBAA)	99.7	107	102	103
Tribromoacetic Acid (TBAA)	102	110	102	106

^aExtract storage stability was conducted on three of the eight extracts from the precision and accuracy study conducted in chlorinated surface water (fortified at 10 µg/L) reported in Table 9. Sample storage stability is expressed as a recovery value (%) calculated as described in Section 9.8.2. All values have been corrected for background levels in the unfortified sample based on triplicate measurements as described in Table 9. Extracts were stored at < -10 °C.

TABLE 14: TAME EXTRACT HOLDING TIME DATA FOR SAMPLES FROM A CHLORINATED SURFACE WATER FORTIFIED WITH METHOD ANALYTES^a

Compound	Initial R (%)	Day 7 R (%)	Day 14 R (%)	Day 21 R (%)	Day 28 R (%)
Monochloroacetic Acid (MCAA)	98.2	95.8	99.2	98.3	98.3
Monobromoacetic Acid (MBAA)	98.9	98.7	99.3	98.8	98.8
Dalapon	100	101	101	100	101
Dichloroacetic Acid (DCAA)	97.5	96.4	97.3	97.4	97.0
Trichloroacetic Acid (TCAA)	99.6	101	100	102	102
Bromochloroacetic Acid (BCAA)	101	102	101	102	102
Bromodichloroacetic Acid (BDCAA)	103	104	102	108	108
Dibromoacetic Acid (DBAA)	100	101	100	102	102
Chlorodibromoacetic Acid (CDBAA)	103	106	102	110	110
Tribromoacetic Acid (TBAA)	104	108	103	110	111

^aExtract storage stability was conducted on three of the eight extracts from the precision and accuracy study conducted in chlorinated surface water (fortified at 10 µg/L) reported in Table 9. Sample storage stability is expressed as a recovery value (%) calculated as described in Section 9.8.2. All values have been corrected for background levels in the unfortified sample based on triplicate measurements as described in Table 9. Extracts were stored at < -10 °C.

TABLE 15: INITIAL DEMONSTRATION OF CAPABILITY (IDC) REQUIREMENTS

Method Reference	Requirement	Specification and Frequency	Acceptance Criteria
Section 9.2.1	Initial Demonstration of Low Method Background	Analyze LRB prior to any other IDC steps	Demonstrate that LRB is reasonably free of contamination.
Section 9.2.2	Initial Demonstration of Accuracy	Analyze QCS Analyze 4-7 replicate LFBs fortified at midrange concentration.	Recovery for analytes must be $\pm 30\%$ of expected value. Mean recovery must be within $\pm 20\%$ of expected value.
Section 9.2.3	Initial Demonstration of Precision	Calculate standard deviation and RSD for replicates used in the Initial Demonstration of Accuracy (Sect. 9.2.2)	RSD must be $\leq 20\%$.
Section 9.2.4	Detection Limit Determination	Over a period of three days, prepare a minimum of 7 replicate LFBs fortified at a concentration estimated to be near the Detection Limit. Analyze the replicates through all steps of the analysis. Calculate the Detection Limit using the equation in Section 9.2.4.	Note: Data from DL replicates are not required to meet method precision and accuracy criteria. If the DL replicates are fortified at a low enough concentration, it is likely that they will not meet precision and accuracy criteria.
Section 9.2.5	Minimum Reporting Limit	Estimate after determining DL and prior to analysis of Field Samples	At least 3 times the DL or at a concentration which yields a signal-to-noise ratio of at least 5.

TABLE 16: QUALITY CONTROL REQUIREMENTS (SUMMARY)

Method Reference	Requirement	Specification and Frequency	Acceptance Criteria
Section 8.3	Sample Shipment and Storage	Samples must not exceed 10 °C within 48 hours after collection. Samples stored in laboratory must not exceed 6 °C.	Sample results are valid only if samples are properly collected, preserved, and stored.
Section 8.4	Sample Holding Time	14 days with appropriate preservation and storage	Sample results are valid only if samples are extracted within sample hold time.
Section 8.4	Extract Holding Time	21 days at ≤ -10 °C protected from light (MTBE) 28 days at ≤ -10 °C protected from light (TAME)	Sample results are valid only if extracts are properly stored and analyzed within extract hold time.
Section 9.3	Laboratory Reagent Blank (LRB)	With each extraction batch of up to 20 samples	Demonstrate that all target analytes in LRB are below $\frac{1}{3}$ the MRL and that interferences do not prevent the identification and quantification of method analytes. If targets exceed $\frac{1}{3}$ the MRL, results for all problem analytes in extraction batch are invalid.

TABLE 16: QUALITY CONTROL REQUIREMENTS (SUMMARY)

Method Reference	Requirement	Specification and Frequency	Acceptance Criteria
Sections 9.4 and 10.3	Continuing Calibration Check (CCC)	Verify initial calibration by analyzing a calibration standard at the beginning of each analysis batch, after every 10 samples, and after the last sample. Low CCC - \leq MRL Mid CCC - near midpoint in initial calibration curve High CCC - near highest calibration standard	1) The result for each analyte must be 70-130% of the expected value for all but the lowest standard. The lowest standard must be 50-150% of the expected value. 2) The peak area of internal standards must be 50-150% of the average peak area calculated during the initial calibration. Results for analytes that do not meet IS criteria or are not bracketed by acceptable CCCs are invalid.
Section 9.6	Internal Standard (IS)	1,2,3 Trichloropropane is added to the extraction solvent.	Peak area counts for the IS in all injections must be within 50-150% of the average peak area calculated during the initial calibration. If the IS does not meet criteria, corresponding target results are invalid.
Section 9.7	Surrogate Standard (SUR)	2-bromobutanoic acid [80-58-0] is added to all calibration standards and samples, including QC samples.	Surrogate recovery must be 70-130% of the expected value. If surrogate fails this criterion, report all results for sample as suspect/surrogate recovery.
Section 9.8	Laboratory Fortified Sample Matrix (LFSM) and Laboratory Fortified Matrix Duplicate (LFSMD)	Analyze one LFSM per analysis batch (20 samples or less) fortified with method analytes at a concentration greater than the native concentration. LFSMD should be used in place of Field Duplicate if frequency of detects for target analytes is low.	Recoveries at mid and high levels not within 70-130% or low-level recoveries not within 50-150% of the fortified amount may indicate a matrix effect. Target analyte RPDs for LFSMD not within $\pm 30\%$ at mid and high levels of fortification and within $\pm 50\%$ near MRL may indicate a matrix effect.

Method Reference	Requirement	Specification and Frequency	Acceptance Criteria
Section 9.9	Field Duplicates (FD)	Extract and analyze at least one FD with each extraction batch (20 samples or less). A LFSMD may be substituted periodically for a FD when the frequency of detects for target analytes is low.	Target analyte RPDs for FD should be within $\pm 30\%$ at mid and high level concentrations and within $\pm 50\%$ near MRL.
Section 9.10	Quality Control Sample (QCS)	Analyzed when new Primary Dilution Standards (PDS) are prepared, during the IDC, or quarterly.	Results must be 70-130% of the expected values.
Section 10.2	Initial Calibration	Use internal standard calibration technique to generate a calibration curve. Use at least 5 standard concentrations.	When each calibration standard is calculated as an unknown using the calibration curve, the result must be 70-130% of the expected value for all except the lowest standard, which must be 50-150% of the expected value.

FIGURE 1: CHROMATOGRAM OF THE HALOACETIC ACIDS ON A DB-1701 COLUMN USING THE CHROMATOGRAPHIC CONDITIONS GIVEN IN TABLE 1 USING THE MTBE PROCEDURE




FIGURE 2: CHROMATOGRAM OF THE HALOACETIC ACIDS ON A DB-1701 COLUMN USING THE CHROMATOGRAPHIC CONDITIONS GIVEN IN TABLE 2 USING THE TAME PROCEDURE




FIGURE 3: CHROMATOGRAM OF THE HALOACETIC ACIDS ON A DB-5.625 COLUMN USING THE CHROMATOGRAPHIC CONDITIONS GIVEN IN TABLE 3 USING THE MTBE PROCEDURE




FIGURE 4: CHROMATOGRAM OF THE HALOACETIC ACIDS ON A DB-5.625 COLUMN USING THE CHROMATOGRAPHIC CONDITIONS GIVEN IN TABLE 4 USING THE TAME PROCEDURE

