

Screening System for Phthalate Esters

Py-Screener




Making the Difficult Simple

The Py-Screener™ system is designed to screen for phthalate esters in polymers. The use of phthalate esters in toys and food packaging is currently restricted. Moving forward, they are expected to be regulated as restricted substances under the Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS (II) Directive) in Europe. The pyrolyzer GC/MS (Py-GC/MS) is used to selectively detect and quantify phthalate esters thermally extracted from samples. This screening system consists of a sampling toolkit, special standards, and special software and can be easily operated even by novices.


Easy to Operate Even for Novices

Organics Solvents Are Not Required for Sample Preparation.


Analytical standards and test samples can be prepared without using organic solvents. To prepare a sample, just use the cutter to remove a portion from the test material, place it in the sample cup, and weigh it. Sample preparation videos provide support so that even novices can easily prepare samples.



Preparation of a Phthalate Ester Standard



Weighing Using an Electronic Balance



Preparation of a Test Sample



Analytical Balances AP Series



Shimadzu balances product lineup

All Required Items Are Available.

Special Standards Developed in Cooperation with SGS Japan Toolkit Required for Sample Preparation

The analytical standards for this system were developed in cooperation with SGS Japan, the market leader for RoHS tests. Standards for sensitivity confirmation, quantitation, and blank tests can be prepared simply by punching out a portion of a standard material using the micro puncher. A toolkit used for preparing samples has been created with Frontier Laboratories Ltd.



Standards Containing Phthalate Esters for Py-GC/MS




Sampling Toolkit

SHIMADZU Py-Screener



Sample Preparation



Analysis

Easy to Operate Using Special Software


Using customized software, operations are easy, even for novices. To automatically start continuous analyses, just place the prepared standards and test samples in the autosampler, and enter the number of samples, the sample names, and their weights. Phthalate ester inspections compliant with the IEC62321-8 international analysis standard, the industry standard for RoHS inspections, can be performed easily by anyone. Shimadzu offers both a simultaneous inspection method for phthalate esters and brominated flame retardants, which supports inspections for a wide range of regulated compounds, and a special high-speed inspection method for phthalate esters, which significantly reduces the inspection time.

	Vial#	Sample Name	Sample Amt.
1	1	Blank_Cup	0.5
2	2	Phthalate_STD_Blank	0.51
3	3	Phthalate_STD_100	0.51
4	4	Phthalate_STD_1000	0.51
5	5	ERM-EC591	0.51
6	6	Test_Sample	0.51
7	7	Test_Sample	0.51
8	8	Test_Sample	0.51
9	9	Test_Sample	0.51
10	10	Test_Sample	0.51

Easy to Operate Even for Novices

Tabular Display of Concentrations and Criteria Clarifies the Results.

The concentrations of target components detected in continuous measurements are displayed in a table and color-coded using criteria based on concentration ranges. The results for continuously measured test samples can be checked at a glance. Also, the system is equipped with accuracy control functions in order to ensure the reliability of blank concentrations, instrument sensitivity, and other data, so even novices can feel confident that they are reporting reliable measurement results.




Py-Screener Software

The special software displayed on the monitor helps you navigate the required procedures. Even novices can operate the system using the software.

Ample Maintenance Support

Maintenance Navigation Supports Long-Term Operation with Periodic Replacement Kits

Using the Maintenance Navigator, the procedures appropriate for pyrolyzer and GC-MS maintenance can be performed easily and confidently. Procedures for locating and resolving leaks are included, and kits with common replacement parts simplify maintenance and troubleshooting procedures.




Maintenance Navigator Windows

A Total Solution Proposed by Shimadzu for the RoHS (II) Directive

This system can be applied to polybrominated biphenyls (PBBs) and polybrominated diphenyl ethers (PBDEs), two types of brominated flame retardants already regulated under the RoHS (II) Directive.

In addition to this system, Shimadzu provides an Energy Dispersive X-ray (EDX) fluorescence spectrometer for inorganic compound screening and various other analytical systems for accurate quantitation. These systems provide a total solution for everything from screening to the accurate quantitation of substances already regulated under the RoHS (II) Directive, and substances for which regulation is anticipated. It is possible to create a report by integrating the inspection data obtained by EDX and GCMS, and since Shimadzu offers RoHS-compatible instruments, a comprehensive management system can be constructed.



Reports can be created by integrating the inspection data obtained by EDX and GCMS.

A comprehensive management system can be constructed by using Shimadzu's RoHS-compatible instruments.

Applicable Systems and Software

GC-MS : GCMS-QP2020 NX, GCMS-QP2020, GCMS-QP2010 Ultra

Pyrolyzer : EGA-PY-3030D multi-shot pyrolyzer

Autosampler : AS-1020E auto-shot sampler

GC/MS Workstation : GCMSsolution™ (Ver.4.30 or later, ver.4.50 or later for GCMS-QP2020 NX) + LabSolutions Insight™

Py Workstation : EGA-PY3030 program (Ver. 1.54 or later)

Caution

- Note that there are no guarantees regarding the accuracy of the information contained in the method files, or the usefulness of the information obtained from the results of their use.
- In order to accurately identify the registered substances, perform the measurements using the system conditions in the method files contained in the product.

Py-Screener, GCMS-QP, GCMSsolution and LabSolutions Insight are trademarks of Shimadzu Corporation.



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

Shimadzu Corporation

www.shimadzu.com/an/