SHIMADZU

Simplified molecular imaging analysis of excreted microbial metabolites using a benchtop MALDI-TOF system

Sidrah Rahman¹, Michael Nairn², <u>Tom Abban²</u>, Simona Salivo², Matt Openshaw², Rian L. Griffiths¹ ¹School of Pharmacy, University of Nottingham, Nottingham, UK. ²Shimadzu, Manchester, UK.

1. Overview

There has always been a high barrier to the application of MALDI imaging as this type of research has historically been the realm of high-end, high-cost MALDI-TOF instruments. There have also been challenges with the sample topography in imaging due to variable colony thickness. This work aims to demonstrate the utility of an affordable, entry-level linear benchtop MALDI-TOF mass spectrometer and its application to locate target biomolecules relevant to polymicrobial biofilms.



2. Introduction

Biofilms (Figure 1) are common throughout healthcare, industrial and environmental settings, with an estimated 65% of bacterial infections biofilm related [1]. Bacterial metabolites such as quorum sensing molecules produced by Pseudomonas aeruginosa (PA) have been validated as lung infection biomarkers. This is particularly relevant in Cystic Fibrosis (CF) [2]. Our aim was to use mass spectrometry (MS) approaches to better understand biomolecules produced by mono- and polymicrobial biofilms, changes that are relevant macroscopically and spatial distributions at the site of the biofilms microscopically. Here we show that accurate mass LESA MS metabolite and lipid profiling can be combined with spatial information from MALDI imaging to offer new insights into the behaviours of bacterial metabolites and lipids produced by biofilms.

3. Methods and Materials

Fleximass-DS MALDI target slides were lightly sanded to enhance agar adhesion. The MALDI target slide was loaded into a target mask and 500 µL of RPMI media poured on to each slide to create a thin agar layer over the MALDI target. Overnight cultures were used to inoculate discs on top of the agar-coated MALDI targets, or 6-well plate for LESA. After 24 hr incubation, the discs were removed (Figure 2) [3]. The MALDI targets were dried prior to applying 9-aminoacridine (9-AA) MALDI matrix using the iMLayer[™] automated matrix sublimation device (Shimadzu, Figure 3). Imaging data was acquired in negative ion mode using a MALDI-8030 dual polarity MALDI-TOF mass spectrometer (Shimadzu, Figure 4). 1:1 methanol: water (quorum sensing metabolites) or 2:1 chloroform: methanol (rhamnolipids) was used for LESA (Advion) MS and MSMS on an Orbitrap Exactive or Q-Exactive (Thermo) in both polarities.







4. Results

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

The University of Nottingham

Assignment	PA	PA-SA	PA-CA	PA- SA-CA
НQ	*	*	*	*
QNO	*	*	X	*
9:1 quinolone	*	*	Х	*
9 quinolone	*	*	Х	*
9 PQS	*	*	Х	Х
11 quinolone	*	*	*	X
yocyanin	*	*	*	Х
ha-C8-C10	*			
ha-C10-C10	*			
ha-C12-C10/C10-C12	*			
ha-C16:1-C10	*			

MP-525