


LCMSsolution Ver.3

LCMSsolution Ver.3 software Helps Provide MS to HPLC Users

The LabSolutions Workstation Series provides a standard operating environment for all laboratory analyzers. LCMSsolution is fully integrated with the LCsolution HPLC workstation software.



Seamless Addition of MS (Just Add the MS Software Unit)

- LCsolution can be updated to LCMSsolution by adding the MS software module. It is the cheapest way to add MS to an existing HPLC instrument and requires no new operations to be memorized.
- All LCsolution and LCMSsolution data files are fully compatible.
- Upgrade your HPLC workstation software with confidence; the data conversion functions handle CLASS-LC10 and CLASS-VP data files.

Using LCMSsolution is easy!

