




- | | | |
|---|----------------------------|--------------|
| 1 | Suction Filter Assembly | 228-17421-91 |
| 2 | Suction Filter (no tubing) | 228-20031-00 |
| 3 | Inlet Check Valve | 228-20256-91 |
| 4 | GFP Plunger Seal | 228-18092-00 |
| 5 | Sapphire Plunger Assembly | 228-16837-91 |
| 6 | Rinse Seal | 228-18306-00 |
| 7 | Outlet Check Valve | 228-09054-93 |
| 8 | Line Filter Assembly | 228-18405-91 |
| 9 | Line Filter Packing | 228-12564-00 |