


LC-20AD LC-20ADXR



- 1 Suction Filter Assembly
- 2 Suction Filter (no tubing)
- 3 Inlet Check Valve
- 4 Graphite-filled PTFE Plunger seal (LC-20AD)
Polyethylene Plunger Seal (LC-20ADXR)
- 5 Sapphire Plunger
- 6 Diaphragm (2 pack)
- 7 Outlet Check Valve
- 8 Line Filter

228-45708-91
 228-45707-91
 228-48249-91
 228-35146-00
 228-32628-91
 228-35601-93
 228-32784-91
 228-45705-91
 228-35871-96