


SIL-10ADvp



- | | | |
|----|---------------------------------|--------------|
| 1 | Suction Filter Assembly | 228-45708-91 |
| 2 | Suction Filter (no tubing) | 228-45707-91 |
| 3 | 6PV Stator Assembly* | 228-36730-91 |
| 4 | 6PV Rotor* | 228-36457-00 |
| 5 | 6PV PEEK Stator | 228-41477-01 |
| 6 | 6PV PEEK Rotor | 228-41476-00 |
| 7 | 6PV Flat Packing | 228- |
| 8 | 6V Rotor | 228-21217-97 |
| 9 | 6V Stator | 228-21220-94 |
| 10 | Needle--for PEEK or Vespel Seal | 228-33518-94 |
| 11 | Vespel Needle Seal | 228-33518-94 |
| | PEEK Needle Seal | 228-33355-01 |
| 12 | SUS Outlet Tubing 0.31 x 600 mm | 228-39830-91 |
| 13 | Rinse Port Septa (10/pack) | 228-36215-91 |
| 14 | Metering Pump Seal | 228-35145-00 |
| 15 | Metering Pump Plunger Assembly | 228-35010-91 |
| 16 | 50 uL Sample Loop | 228-39436-91 |

6-Port Valve (6V)

6-Position Valve (6PV)

*(For instruments with Serial Numbers prior to C2105109409)

6-Position Valve (6PV)
(For instruments with Serial Number C2105109409 and later)

